CANINE PANCREATIC LIPASE TEST KIT

The VetScan cPL Rapid Test is a rapid test that must be interpreted by the VetScan VUE for the semi-quantitative determination of pancreas-specific lipase levels in **canine serum or plasma samples**.

INTRODUCTION

Pancreatitis in dogs is a painful and potentially fatal disease that can be difficult to diagnose due to non-specific signs of disease that are also present in many other ailments.

Signs of pancreatitis may include:

- Fever
- Vomiting
- Diarrhea
- Anorexia
- Dehydration
- Lethargy
- Restlessness
- Depression
- Abdominal pain
- Elevated heart rate
- Respiratory distress

Inflammation of the pancreas can disrupt the flow of enzymes into the digestive tract and these enzymes can leak out of the pancreas into the abdominal cavity causing severe discomfort. The enzymes may begin to break down fat and proteins in the pancreas and may also affect nearby organs, such as the liver and kidneys. Pancreatitis can progress rapidly in dogs if not recognized and treated quickly. With timely treatment and monitoring, the illness may be cured without lasting damage to the organs.

Inflammation of the pancreas leads to leakage of a substance, “pancreatic lipase,” into the peripheral blood stream. Published studies have shown abnormally elevated concentrations of pancreatic lipase in plasma and serum may be indicative of pancreatitis or another whole body inflammatory condition. Pancreatic lipase concentration in the peripheral blood is well accepted as a more reliable method of diagnosis than serum amylase or lipase activities (1,2). There are various sources for both lipase and amylase in the body which makes measuring these analytes an unreliable method for detecting pancreatitis (3). Thus, a pancreatic lipase assay is useful in determining the health of the pancreas; an elevated cPL is suggestive of pancreatitis.

The VetScan cPL Rapid Test employs affinity purified antibodies that bind to canine pancreatic specific lipase to assess its concentration. These purified polyclonal antibodies are bound to colloidal gold particles and use a cross-species double antibody sandwich assay to visualize the presence of pancreatic lipase. Once bound to antibody-coated gold particles, the lipase flows through the strip and is captured by immobilized antibodies on the test strip. The accumulation of the captured gold particle/enzyme complex causes color to become visible on the Test line (T). Signal is further amplified by use of a competitive antibody calibration scheme employed on the Control line (C). The darkness of the lines is quantified by densitometric analysis in the VetScan VUE. The darkness of the test and control line correlate predictably with the concentration of pancreatic lipase in the sample. A Control Line (C) will always appear whether the sample is high or low. A Test line (T) may appear even in normal or healthy dogs – the presence of a test line is not necessarily indicative of abnormal levels of cPL. This test cannot be read visually – it must be read with a VetScan VUE.

INSTRUCTIONS FOR USE

- This test is for the immunological detection of canine pancreatic lipase in canine plasma or serum. **This test is not calibrated for use with whole blood.**
- Samples must be at room temperature (15-30 °C/59-86 °F) before running the test.
- Serum or plasma, fresh, or stored refrigerated at 2-8 °C (35-46 °F) for up to one week may be used.
- For longer storage, serum or plasma, may be frozen (-20 °C/-4 °F or colder). Frozen samples must be thawed and centrifuged at > 1600g to remove any particulate material before use.
- Excessive hemolysis may obscure the results and will result in an error reported by the VetScan VUE when attempting to read the test. Icteric or lipemic samples will not affect test results.
- EDTA, sodium and lithium heparin or ACD in plasma will not affect the results

PRECAUTIONS AND WARNINGS

- Do not remove the Test Device from the pouch until ready for use.
- The Test Device must be used as soon as possible after removing from the pouch and within a maximum of 15 minutes.
- For veterinary use only.
- Do not use the components after expiration date.
- The Test Device should be in a horizontal position, avoid tilting when inserting into the VetScan VUE.
- Use a separate transfer pipette for each test.
- The chase buffer is not interchangeable from lot to lot.
- Do not use a Test Device from a pouch that appears tampered with, torn or damaged.
- Do not use a Test Device if it appears cracked, broken, or otherwise damaged.
- The kit components must not be frozen.
- Important: Results must be interpreted by the VetScan VUE. Refer to VetScan VUE user instructions. DO NOT ATTEMPT TO READ THE TEST VISUALLY.**

STORAGE

- The test devices and chase buffer can be stored at room temperature or refrigerated (2-30°C/35-86°F) but should never be frozen.
- Test devices and chase buffers are stable until expiration when stored at recommended temperatures.

KIT COMPONENTS

- Test Devices
- Chase Buffer Bottle
- Transfer Pipettes
- Instructions for Use

TEST PROCEDURE

- Remove the Test Device from the protective pouch and place on a flat surface. Label the Test Device with the patient I.D. or control identification. Only use the provided white space for this purpose. Do not cover up barcodes or the test label.
- Using the Transfer Pipette provided, hold the pipette vertically and dispense one drop of sample (serum or plasma) into the sample well.
- Hold the Chase Buffer Bottle vertically, add 3 drops of the chase buffer into the sample well.
- Insert cassette into the VetScan VUE and initiate a Timed Scan; Results will be sent electronically to the VetScan VUE app. **DO NOT** read results visually.

INTERPRETATION OF TEST RESULTS

Test results must be read by the VetScan VUE. The test cannot be interpreted visually.

The VUE will display the pancreatic lipase serum concentration in micrograms per liter. Results < 200µg/L are considered normal and within the reference range. Results from 200µg/L to 400µg/L are abnormal and may be indicative of pancreatitis. Patients should be re-evaluated in 2-3 weeks if possible. Results higher than 400µg/L are consistent with pancreatitis.


The range indicated on the visual display encompasses one standard deviation above and below the calculated value. Each individual rapid test is subject to variability and the results should be considered with this in mind.

Invalid Results

Results scanned more than 15 minutes after addition of the chase buffer are not diagnostic and should be ignored.

Damage to the cassette or debris on the reading window will interfere with reading of the device. This will likely result in an error and no result given. If any of these things are observed the sample should be re-run on another test device.

This test can only indicate the concentration of pancreatic lipase. It cannot diagnose disease. All results should be considered in the context of all available clinical data. Any treatment administered to a patient should be subject to conclusions drawn from all clinical findings. It is possible a patient may have concurrent disease.

REFERENCES

- McCord, K. (2012). “A multi-institutional study evaluating the diagnostic utility of the spec cPLTM and SNAP® cPLTM in clinical acute pancreatitis in 84 dogs.” J Vet Intern Med 26(4); pp. 888-96.
- Steiner, J. and Newman, S. (2008). “Sensitivity of Serum Markers for Pancreatitis in Dogs with Macroscopic Evidence of Pancreatitis.” Veterinary Therapeutics Winter 9.4; pp. 263-73
- Steiner, J. (2017). “Chapter 290: Pancreatitis: Diagnosis and Treatment” in Ettinger, S., Feldman, E., Cote, E. ed. Textbook of Veterinary Internal Medicine Expert Consult, 8th ed. 2 vols. Saunders, pp. 1684

■ Abaxis, Inc. 3240 Whipple Road, Union City, CA 94587 USA | +1 800 822 2947 | www.abaxis.com
ABAXIS Europe GmbH, Bunsenstr. 9-11, 64347 Griesheim, Germany | +49 6155 780 210
For patent information, see www.abaxis.com/about_us/patents

Learn more about all of our products and services at www.abaxis.com
Abaxis and VetScan are registered trademarks of Abaxis, Inc. © Abaxis 2017.

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
REFERENCIAS

- McCord, K. (2012). “A multi-institutional study evaluating the diagnostic utility of the spec cPLTM and SNAP® cPLTM in clinical acute pancreatitis in 84 dogs.” J Vet Intern Med 26(4); págs. 888-96.
- Steiner, J. and Newman, S. (2008). “Sensitivity of Serum Markers for Pancreatitis in Dogs with Macroscopic Evidence of Pancreatitis.” Veterinary Therapeutics Winter 9.4; págs. 263-73
- Steiner, J. (2017). “Chapter 290: Pancreatitis: Diagnosis and Treatment” in Ettinger, S., Feldman, E., Cote, E. ed. Textbook of Veterinary Internal Medicine Expert Consult, 8.ª ed. 2 vols. Saunders, pág. 1684

■ Abaxis, Inc. 3240 Whipple Road, Union City, CA 94587 (Estados Unidos) | +1 800 822 2947 | www.abaxis.com
ABAXIS Europe GmbH, Bunsenstr. 9-11, 64347 Griesheim (Alemania) | +49 6155 780 210
Para obtener información sobre patentes, consulte www.abaxis.com/about_us/patents.

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
REFERÊNCIAS

- McCord, K. (2012). “A multi-institutional study evaluating the diagnostic utility of the spec cPLTM and SNAP® cPLTM in clinical acute pancreatitis in 84 dogs.” J Vet Intern Med 26(4); págs. 888-96.
- Steiner, J. and Newman, S. (2008). “Sensitivity of Serum Markers for Pancreatitis in Dogs with Macroscopic Evidence of Pancreatitis.” Veterinary Therapeutics Winter 9.4; págs. 263-73
- Steiner, J. (2017). “Chapter 290: Pancreatitis: Diagnosis and Treatment” in Ettinger, S., Feldman, E., Cote, E. ed. Textbook of Veterinary Internal Medicine Expert Consult, 8th ed. 2 vols. Saunders, pp. 1684

■ Abaxis, Inc. 3240 Whipple Road, Union City, CA 94587 EUA | +1 800 822 2947 | www.abaxis.com
ABAXIS Europe GmbH, Bunsenstr. 9-11, 64347 Griesheim, Alemanha | +49 6155 780 210
Para informaçoês relativas à patente, visite www.abaxis.com/about_us/patents

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Abaxis e VetScan são marcas comerciais registadas da Abaxis, Inc. © Abaxis 2017.