# SWINE INFLUENZA VIRUS TYPE A ANTIGEN TEST KIT

# FIUDFTFCT\* Swine

## **ENGLISH**

results indicate

### I. GENERAL INFORMATION FluDETECT™ Swine is an in vitro, lateral flow immunoassay designed to aid in the qualitative detection of Swine Influenza (SIV) Type A virus in

swine nasal samples. Positive results may be submitted to a reference lab for confirmation and subtype determination. Negative that no detectable SIV is present. II. TEST PRINCIPLES

FlubETECT" Swine is based on Rapid Immuno Migration (RIM") technology. The test strip uses two antibodies that are specific to the p56 nucle-oprotein of Swine Influenza Virus. An anti-influenza A antibody bound to Influenza A antigen present in the sample forms an antigen-antibody complex which migrates along a strip and is captured on a sensitized reaction line by the second antibody. The accumulation of the complex causes the formation of a dearly visible pink/purple band. The presence of a control band, located above the reaction line, ensures that the test was performed correctly. III. SAMPLE COLLECTION

- Use the provided swabs to collect nasal samples from swine. See Section VI. PRECAUTIONS. Samples obtained while swine exhibit clinical signs will contain the highest detectable amount of virus. If additional testing is to be performed on the sample, it is possible to collect the sample in a viral transport medium. These alternate media
- re not provided in this test kit. Approved media that may be used, in order of preference are: Brain Heart Infusion Broth porcine orig BRID, ris Buffered Tryptone Broth (TBTB), Nutrient Broth (NB) or Peptone Broth (PB). If further testing is planned, follow Section VILB. ample Extraction Method.
- IV. SAMPLE STORAGE Samples should preferably be tested immediately after collection. If testing is delayed, samples should be kept refrigerated (up to 48 hours at

For long term storage, samples should be kept frozen (-70 °C or colder). Do not store samples at -20 °C. Do not store samples in a self-defrosting freezer. Avoid multiple freeze-thaw cycles. V KIT CONTENTS

## 1 Vial containing 20 test strips and desiccant.


- 1 Extraction Ruffer hottle (6.0 ml) 20 Swabs
- 20 Test tubes
- 20 Test tube caps
- 1 Test tube rack
- Instructions for use.
- VI. PRECAUTIONS

- Do not use this kit or any of its components after the expiration date.

- Kit should be stored at 2 °C -30 °C. Kit should not be frozen
- Do not mix materials from different test kits. The vial holding the test strips contains a desiccant and should be kept tightly closed when not in use.
- Use the test strips within 10 minutes after removal from the desiccant vial. Test strips should only be handled in the upper, labelled region. Avoid contact with the surface of the test strip.
- The test-strip should be placed in the test-tube-vertically.
- Use a separate swab for each sample. Swabs with wooden handles or containing calcium alginate may interfere with the test and must not be used.
- Do not centrifuge samples prior to use

Swabs containing visible blood may partly obscure a weak positive band due to hemoglobin background.

- Handle all reagents and samples as biohazardous material. Extraction buffer is preserved with sodium azide
- For veterinary use only. SAMPLE EXTRACTION
- Zoetis recommends two methods of sample ext be performed on the sample, follow METHOD B. xtraction. If the sample will be tested using this test kit only, follow METHOD A. If further testing is to
- Extraction Buffer Provided in kit box
- Place 8 drops (approximately 0.25 mL) of Extraction Buffer in the test tube provided (Figure 1, Step 1). Place the swab containing the sample in the tube and rotate the swab 5 - 10 times in the buffer (Figure 1, Step 2).
- When removing the swab from the tube, press the swab against the side of the tube repeatedly until no more liquid comes from the swab (Figure 1, Step 3). Discard the swab in an appropriate biohazard container.
- If the extracted sa SAMPLE STORAGE.. acted samples will not be tested immediately, cap the tube with the provided cap and store the sample according to Section IV.
- FIGURE 1 Sample Extraction and Test Procedure Method "A



- Discard the swab in an appropriate biohazard container or break the handle of the swab below the top of the tube such that the tube containing the swab tip can be sealed with the provided cap.
- If the extracted sample will not be tested immediately, cap the tube with the provided cap and store the sample according to Section IV.
- SAMPLE STORAGE.
- NOTE: 0.2 mL of extracted sample is required for each test; remaining volume can be used in alternate test methods.
- Allow samples and kit to come to 15 °C 30 °C before testing. Testing samples - Use either Method A or B as appropriate.

Remove a test strip from the desiccant vial for each sample to be tested. Handle the test strip on the labeled portion of the strip, (Figure I, Step 4) **Note:** Prior to use, test and control bands appear yellow. The bands are dyed yellow for quality control purposes. The dye does no interfere with the test results and will wash away while the test is developing.

Place the test strip directly into the test tube containing the sample. Place test strip so that the pink pad is submerged in the extracted sample. Incubate the test strip in the sample for 15 minutes. (Figure 1, Step 5) Remove the test strip from the test tube to read.

B. Sample Extracted in Viral Transport Media Place 0.2 mL of the viral transport media into the test tube provided.

A. Sample Extracted in Buffer:

- Add 3 drops of Extraction Buffer to tube; tap side of tube to mix. Remove a test strip from the desiccant vial for each sample to be tested. Handle the test strip on the labeled portion of the strip. (Figure 1, Step 4) Note: Prior to use, test and control bands appear yellow. The bands are dyed yellow for quality control purposes. The dye does not
- - interfere with the test results and will wash away while the test is developing. Place the test strip directly into the test tube containing the sample. Place test strip so that the pink pad is submerged in the extracted
- sample. Incubate the test strip in the sample for 15 minutes. (Figure 1, Step 5) Remove the test strip from the test tube to read.
- IX. READING TEST After 15 minutes, (Figure 1, Step 6). , observe the presence or absence of pink/purple bands in the center of the test strip between the two absorption pads

Discard the test strip in an appropriate biohazard container

- The control band appears in the upper end of the test strip (closest to the handle), while the sample test results are read in the lower part of the test strip
- Test is **VALID** if a pink/purple band (Control Line) is present in the upper part of the test strip. The absence of the Control Line indicates that the test is invalid and must be repeated (Figure 1, Step 6).

X. RESULTS Valid Results

- Interpretation of Results POSITIVE for Swine Influenza Virus: Two pink/purple bands (Control Line and Test Line) are clearly visible on the test strip (C & T). A POSITIVE result indicates that a detectable level of SIV is present in the sample. Positive samples can be submitted to an approved reference laboratory for confirmation and subtype determination.
- **NEGATIVE** for Swine Influenza Virus: One pink/purple band (Control Line) is present in the upper part of the test strip (C). A **NEGATIVE** result indicates that no detectable SIV is present in the sample.

Very faint lines may be due to non-specific binding and should be further investigated.

SYMBOL DESCRIPTIONS

Use by (expiration date)

Temperature limitations

(storage temperature range)

Note: The Control Line on the upper part of the test strip may appear prior to the end of the incubation period. This does not mean that the test is complete, as a test band may appear more slowly than the control band. The test strip must incubate for a full 15 minutes before a sample is interpreted as Negative. The test can be considered complete if the Test Line on the low part of the strip appears before the 15 minute incubation period is over. This sample is interpreted as Positive. If the test strip remains in the test tube for more than 20 minutes a false positive ghost band could appear in place of the reaction band (f).

EC REP

 $\prod_{i}$ 

IVD

Authorized Representative in the European Community

Consult Instructions for use

Manufacture

In vitro diagnostic medical device

## zoetis Zoetis Inc. Kalamazoo, MI 49007, USA

**LOT** Batch code

Serial number

SN

VLN/PCN 190/5181.51 1-888-963-8471

EC REP ZOETIS FRANCE 23 Rue Pierre Gilles de Gennes

69007 Lyon, FRANCE