

i-STAT[®]Alinity v
Hospital Resource Guide





The i-STAT Alinity v delivers blood gas, acid-base, electrolyte, chemistry, and hematology measurements in a completely portable, handheld package. Accuracy is ensured by extensive quality checks and calibrations that occur automatically with each cartridge run. Results are obtained in as little as three minutes - making it the ideal solution for critical care situations, anesthetic monitoring, and fieldwork.

Cartridge Storage:

Refrigerate at 2 °C to 8 °C (35 °F to 46 °F).

Cartridge Stability:

Cartridges may be stored at room temperature 18-30 °C (64-86 °F), but this will decrease the shelf life. Refer to the cartridge box for room storage shelf life information. Once a cartridge has been warmed to room temperature, do not return it to the refrigerator.

Allow the cartridge to warm for 5 minutes at room temperature before removing from the pouch for analysis.

Use cartridges immediately after opening pouch.

Sample Preparation and Considerations:

- Whole blood samples without anticoagulant or whole blood collected into a lithium heparin tube may be used.
- Blood may be either venous or arterial, depending on the analytes to be measured.
- Venous samples are typically performed for acid-base, electrolyte, and hematologic studies.
- Samples for iCa should be collected in balanced heparin.
- For most accurate results, run samples immediately after collection.
 - Samples for pH, pCO₂, pO₂, TCO₂, and iCa should be tested within 10 minutes if stored anaerobically.
 - All other analytes should be tested within 30 minutes.

For additional information regarding individual cartridges and tests sample collection and handling, see Cartridge & Test Information sheets: www.pointofcare.abbott



Acid-Base Utilization

Acid-base analysis is vital to your diagnostic protocols¹

Chemical reactions, especially those occurring *in vivo*, are dependent on many factors, none more important than optimal pH. Illness, whether acute or chronic, often results in pH abnormalities. Failure to recognize and address these abnormalities may result in:

- Missed diagnoses
- Inappropriate treatment
- Delayed or poor patient response to therapy
- Increased time in hospital
- Frequent relapse
- Inability to thrive
- Patient death

| Acid-base definitions | |
|------------------------------------|---|
| pH | Measurement of the H ⁺ ion concentration |
| pCO₂ | Partial pressure of the carbon dioxide; reflects the amount of carbonic acid present |
| HCO₃⁻ | Bicarbonate, the body's major buffer |
| Anion Gap | Represents the concentration of all unmeasured anions in the plasma; the difference between measured cations and measured anions (Na ⁺ + K ⁺) - (Cl ⁻ + HCO ₃ ⁻); helpful in determining the cause of acid-base abnormalities. |
| Base Excess | mEq/L of strong base or acid needed to return the pH to 7.40. |
| Electrolytes | Na ⁺ , K ⁺ , Cl ⁻ |
| TCO₂ | Total carbon dioxide, which is primarily HCO ₃ ⁻ (95%) |
| pO₂ | Partial pressure of oxygen; measurement of the tension or pressure of oxygen dissolved in blood |

Note: A venous sample is acceptable for interpretation of acid-base parameters. For detailed information on pO₂, an arterial sample is recommended.

Acid-Base Diagnostic Chart³

Step 1


Determine whether the pH is high or low

Step 2

Determine the primary disturbance

Step 3

Determine whether there is metabolic or respiratory compensation



Common Disease States Where Acid-Base Analysis Is Beneficial

| EXPECTED ACID-BASE ABNORMALITIES (depending on species) | | | | | |
|---|---|--|--|---|--|
| ACIDEMIA ⁴ | | | ALKALEMIA ⁵ | | |
| pH < 7.35 (canine and equine) pH < 7.25 (feline) | | | pH > 7.45 (canine and equine) pH > 7.40 (feline) | | |
| Metabolic acidosis ↑ H ⁺ >>> ↓ pH (Most common presentation) ↓ HCO ₃ ⁻ >>> ↓ pH (rare in small animals) | | | Metabolic alkalosis ↑ HCO ₃ ⁻ or ↓ H ⁺ >>> ↑ pH | | Respiratory alkalosis ↓ O ₂ >>> hyperventilation >>> ↓ pCO ₂ >>> ↑ pH Reduced ability to uptake or exchange O ₂ |
| LACTIC ACIDOSIS | VOMITING/DIARRHEA | RENAL FAILURE | DIABETIC KETOACIDOSIS | UPPER GI OBSTRUCTION | RESPIRATORY |
| <ul style="list-style-type: none"> • An increase in lactic acid production as a result of decreased tissue perfusion and/or decreased oxygenation • Occurs in many disease states, most commonly: <ul style="list-style-type: none"> • Hypovolemia/shock • Vomiting/diarrhea • Colic • Gastric torsion (GDV) | <ul style="list-style-type: none"> • Lactic acidosis secondary to hypovolemia • +/- loss of sodium bicarbonate (NaHCO₃) • Electrolyte abnormalities • Anion gap often normal | <ul style="list-style-type: none"> • Uremic toxins increase acid levels • Loss of sodium bicarbonate (NaHCO₃) OR hydrogen ion retention (H⁺) • Electrolyte abnormalities • Lactic acidosis with anemia and/or severe dehydration | <ul style="list-style-type: none"> • Ketoacids • Lactic acidosis • Electrolyte abnormalities • High/normal anion gap, depending on severity | <ul style="list-style-type: none"> • Loss of Cl⁻ in the form of HCl (hydrochloric acid) • Hypochloremia is common • Potential loss of free body water | <ul style="list-style-type: none"> • Hyperventilation • Pain • Iatrogenic (mechanical ventilation) • Decreased tissue perfusion (due to anemia, dehydration, other) • Compensation for metabolic acidosis (hyperventilation) • Head trauma |
| CARTRIDGE CHOICES | | | | | |
| CG4+: Acid-base, lactate, pO₂, TCO₂ Helpful with GDV and other severe GI cases Diagnosis and monitoring for emergencies and/or severe cases | | | CG8+: Acid-base, pO₂, HCT, glucose, Na, K, iCa Helpful for monitoring diabetic and kidney disease patients Diagnosis and monitoring for emergencies and/or severe cases Neoplasia screening | | |
| <p>Disclaimer: Cartridge examples are for suggestive purposes only. Diagnostic testing choices should be based on medical history, physical examination and the patient's response to treatment.</p> | | | | | |

i-STAT Alinity v Cartridge Test Menu

The i-STAT Alinity v uses a wide range of disposable, single-use cartridges that contain the necessary reagents to provide reference lab quality results, while improving efficiency throughout the animal health continuum of care.

| | | CHEM8+ | CG4+ | CG8+ | G | Crea |
|--------------|--|--------|------|------|---|------|
| Hematology | Hematocrit (Hct) | ● | | ● | | |
| | Hemoglobin (Hb)* | ● | | ● | | |
| Chemistry | Blood Urea Nitrogen (BUN) | ● | | | | |
| | Creatinine (Crea) | ● | | | | ● |
| | Ionized Calcium (iCa) | ● | | ● | | |
| | Glucose (Glu) | ● | | ● | ● | |
| Electrolytes | Chloride (Cl) | ● | | | | |
| | Sodium (Na) | ● | | ● | | |
| | Potassium (K) | ● | | ● | | |
| Acid Base | pH | | ● | ● | | |
| | Partial Pressure of Carbon Dioxide (<i>PCO</i> ₂) | | ● | ● | | |
| | Bicarbonate (HCO ₃)* | | ● | ● | | |
| | Total Carbon Dioxide (TCO ₂)* | ● | ● | ● | | |
| | Anion Gap (AnGap)* | ● | | | | |
| | Base Excess (BE)* | | ● | ● | | |
| Blood Gas | Partial Pressure of Oxygen (<i>PO</i> ₂) | | ● | ● | | |
| | Oxygen Saturation (sO ₂)* | | ● | ● | | |
| Specialty | Lactate (Lac) | | ● | | | |

*Calculated Value. Note TCO₂ is a measured values on the CHEM8+ cartridge, but is a calculated value on the CG4+and CG8+, cartidges.

i-STAT Alinity v System and Reference Interval⁶

| | | Units | System interval | Reference interval*** | | |
|------------------------|--|--------|-----------------|-----------------------|--------------|-------------|
| | | | | Canine | Feline | Equine |
| Hematology | Hematocrit (Hct) | % PCV | 15 - 75 | 35 - 57 | 26 - 50 | 25 - 44 |
| | Hemoglobin (Hb)* | g/dL | 5.1 - 25.5 | 12 - 19 | 9 - 17 | 8 - 15 |
| Chemistry | Blood Urea Nitrogen (BUN) | mg/dL | 3 - 140 | 7 - 26 | 17 - 35 | 4 - 27 |
| | Creatinine (Crea) | mg/dL | 0.2 - 20.0 | 0.5 - 1.4 | 0.8 - 2 | 0.7 - 2 |
| | Ionized Calcium (iCa) | mmol/L | 0.25 - 2.50 | 1.21 - 1.45 | 1.04 - 1.44 | 1.31 - 1.83 |
| | Glucose (Glu) | mg/dL | 20 - 700 | 81 - 118 | 70 - 161 | 71 - 111 |
| Electrolytes | Chloride (Cl) | mEq/L | 65 - 140 | 109 - 121 | 116 - 127 | 95 - 105 |
| | Sodium (Na) | mEq/L | 100 - 180 | 141 - 150 | 145 - 157 | 132 - 139 |
| | Potassium (K) | mEq/L | 2.0 - 9.0 | 3.3 - 4.9 | 3.4 - 4.9 | 2.6 - 5.8 |
| Acid-Base | pH | | 6.5 - 8.2 | 7.32 - 7.44 | 7.28 - 7.46 | 7.37 - 7.46 |
| | Partial Pressure of Carbon Dioxide (<i>PCO</i> ₂) | mmHg | 5 - 130 | 26 - 45 | 25 - 42 | 39 - 52 |
| | Bicarbonate (HCO ₃)* | mEq/L | 1.0 - 85.0 | 16 - 26 | 15 - 24 | 25 - 33 |
| | Total Carbon Dioxide (TCO ₂) | mEq/L | 5 - 50 | 16 - 26 | 16 - 24 | 25 - 33 |
| | Anion Gap (AnGap)* | mEq/L | (-10) - (+99) | 8 - 21 | 8 - 20 | 5 - 17 |
| | Base Excess (BE)* | mEq/L | (-30) - (+30) | (-9) - (+1) | (-11) - (-1) | 0 - 9 |
| Blood Gas (arterial)** | Partial Pressure of Oxygen (<i>PO</i> ₂) | mmHg | 5 - 800 | 85 - 100 | 90 - 110 | 62 - 170 |
| | Oxygen Saturation (sO ₂)* | % | 0 - 100 | 95 - 100 | 95 - 100 | 96 - 100 |
| Blood Gas** (venous) | Partial Pressure of Oxygen (<i>PO</i> ₂) | mmHg | 5 - 800 | 25 - 70 | 27 - 51 | 22 - 80 |
| | Oxygen Saturation (sO ₂)* | % | 0 - 100 | 49 - 100 | 52 - 91 | 49 - 100 |
| Specialty | Lactate (Lac) | mmol/L | 0.30 - 20.00 | 0.4 - 2.8 | 0.4 - 2.6 | 0.3 - 1.1 |

*Calculated Value

**Arterial blood gas ranges are built into software. Venous blood gas ranges are not available in the software at this time.

***Reference interval are for venous samples unless specified

Highlighted cells reflect interval for arterial samples. Equine arterial ranges developed for i-STAT Alinity v. Canine and feline arterial interval developed for i-STAT 1.

Reference intervals are provided only as a guideline. The most definitive reference intervals are those established for your patient population and using individualized patient trends. Test results should be interpreted in conjunction with the patient's clinical signs.

10 Sylvan Way
Parsippany, NJ 07054 USA
Tel +1 888 963 8471

3240 Whipple Road
Union City, CA 94587 USA
Tel +1 800 822 2947

www.zoetisUS.com/istatalinityv

*Arterial blood gas ranges are built into software. Venous blood gas ranges are not available in the software at this time.

¹Monnig AA. Practical Acid-Base in Veterinary Patients. *Veterinary Clinics of North America: Small Animal Practice*. 2013; 43: 1273-1286. doi:10.1016/j.cvsm.2013.07.009.

²George JW, Zabolotzky SM. Water, Electrolytes, and Acid Base [Chapter 5]. *Duncan & Prasse's Veterinary Laboratory Medicine*. 2011: 147-150.

³Kerl ME. Acid-Base, Oximetry, and Blood Gas Analysis [Chapter 128]. *Textbook of Veterinary Internal Medicine Expert Consult*. Eighth Edition. 2016: 531-535.

⁴Flaherty D, Blackwood L. Blood gas analysis and acid-base disorders [Chapter 9]. *BSAVA Manual of Canine and Feline Clinical Pathology*, Third Edition. 2016: 169-171.

⁵Flaherty D, Blackwood L. Blood gas analysis and acid-base disorders [Chapter 9]. *BSAVA Manual of Canine and Feline Clinical Pathology*, Third Edition. 2016: 172-173.

⁶Data on File, Study report DH65R-US-19-084, Zoetis Inc., 2021. Arterial ranges for canine and feline reflect historical data and were not updated in this study.

For Cartridge & Test Information sheets, please visit: www.pointofcare.abbott

i-STAT and Alinity are trademarks of the Abbott Group of Companies in various jurisdictions. All trademarks are the property of Zoetis Services LLC or a related company or a licensor unless otherwise noted.

© 2022 Zoetis Services LLC. All rights reserved. ABX-00075 R3

The Zoetis logo, featuring the word "zoetis" in a lowercase, orange, sans-serif font. The letter 'z' is stylized with a curved underline that extends under the 'o'.