# FELINE LEUKEMIA VIRUS ANTIGEN TEST KIT

ViraCHEK™ | FeLV

ENGLISH

For the detection of Feline Leukemia Virus (FeLV) antigen in feline whole blood, serum, and plasma.

## GENERAL INFORMATION AND INTENDED USES

The plastic wells contain antibodies directed against the FeLV group specific antigen, p27. A monoclonal antibody directed against p27 is labeled with the enzyme horseradish peroxidase (HRP). The specimen (whole blood, plasma or serum) is incubated simultaneously with the solid phase and enzyme-labeled antibodies. If present, p27 group specific antigen is bound to the antibody coated well and enzyme-linked antibody at the same time. Any free enzyme linked antibody is washed away and a chromogenic substrate is added. The development of a distinctly blue color specifically indicates the presence of FeLV. In the absence of FeLV, no color change will be observed.

ViraCHEK™/FeLV is highly specific, sensitive and simple to use. Test results for circulated antigen can be obtained in 10 minutes or less. Most cats initially become antigen positive within 2-3 weeks after virus exposure. Progressive infections will remain positive while regressive infections may revert to antigen negative 2-8 weeks, or more, later.

## KIT COMPOSITION AND CONSERVATION

Contains materials sufficient to test 48 to 140 samples.

ITEM	REAGENT NATURE	VOLUME	RECONSTITUTION AND CONSERVATION
<b>A</b>	Anti-FeLV Antibody Coated Wells	3 sets of 4x12 wells	Ready to use.
<b>CONTROL +</b>	Positive Control; preserved with Gentamicin Sulfate and Amphotericin B	3.0 mL	Ready to use. Red Cap.
<b>CONTROL -</b>	Negative Control; preserved with Gentamicin Sulfate and Amphotericin B	3.0 mL	Ready to use. Gray Cap.
<b>1</b>	HRP Monoclonal Antibody Conjugate; preserved with Phenol and Microcide III	7.5 mL	Ready to use. Blue Cap.
<b>2</b>	Substrate Buffer	2 x 7.0 mL	Ready to use. Purple Cap.
	Well holder		

Store all reagents provided in the kit at 2 – 7 °C. Reagents should not be frozen.

## REAGENTS REQUIRED TO PERFORM 48 TESTS

- 48 Anti-FeLV Antibody Coated Wells
- 1.0 mL Positive Control
- 1.0 mL Negative Control
- 2.5 mL Conjugate
- 1 Bottle Substrate Buffer

## EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- Marking Pen
- Deionized or distilled water or normal saline
- Wash Bottle
- Precision Pipettes
- Timer

## WARNINGS TO THE USERS OF REAGENTS AND ANTIBODY COATED PLATES

- Handle all reagents and samples as biohazardous material. It is recommended to dispose reagents and contaminated material according to the applicable regulations.
- Wear suitable protective clothing.
- Irritating to skin and eyes. Keep all reagents away from skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Take care not to contaminate any test reagents with samples or bacterial agents.
- The best results are achieved by following the protocols described below, using good, safe laboratory techniques.
- Do not use this kit or any of its contents after the expiration date and do not intermix components from different serial numbers.
- Use a separate pipette tip for each sample.
- Follow instructions exactly. Improper washing or contamination of reagents may produce nonspecific color development.
- Do not expose kit to direct sunlight.
- For veterinary use only.
- NEVER PIPETTE BY MOUTH. Harmful if swallowed.

**NOTE: Allow all components to come to 21 – 25 °C before starting.**

## SAMPLE COLLECTION AND STORAGE

- Follow proper sample collection procedures.
- 50 µL (0.05 mL) of feline anticoagulated whole blood, plasma or serum is required.
- Harvest sample and store properly (up to seven days at 2 - 7 °C for serum and plasma, up to 24 hours for whole blood).
- For prolonged storage, serum and plasma samples should be kept frozen (-20 °C or colder).
- Test only good quality samples (i.e. avoid bacterial contamination, heavy hemolysis or lipemia). When in doubt, obtain a better quality sample.

**FOR THE TEST PROCEDURE PLEASE SEE DETAILED INSTRUCTIONS IN THE BOTTOM PANEL (BELOW).**

ViraCHEK™

ViraCHEK™

DANGER / DANGER / PELIGRO / PERICOLO / PERIGO / GEFAHR / NIEBEZPIECZEŃSTWO / อันตราย

**Substrate Buffer:** Danger. Causes serious eye irritation. May damage the unborn child. Harmful to aquatic life with long-lasting effects. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves / protective clothing / eye protection / face protection. Wash thoroughly after handling. Avoid release to the environment. If exposed or concerned: get medical advice / attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical advice / attention. Store locked up. Dispose of contents / container in accordance with local / regional / national / international regulations.

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## INTERPRETATION OF RESULTS

### Controls:


- POSITIVE control should be distinctly blue.
- NEGATIVE control should be completely clear.

### Samples:

- POSITIVE samples will be blue. Color intensity will vary with level of FeLV antigen present.
- NEGATIVE samples will be clear. Compare directly with the negative control against a white background.

### NOTES

- When mixing well contents in Steps 4 and 11, tap microplate rapidly for a full 15 seconds.
- Washing is the most important step. Microwells cannot be overwashed.** Underwashing will result in nonspecific blue color development in the negative control and sample wells.
- Remember! Use saline as the wash solution with whole blood samples.
- Read results at 5 minutes. If no color is seen at 5 minutes, the sample is negative.

## SYMBOL DESCRIPTIONS

	Use by Date (expiration date)		Authorized Representative in the European Community
	Batch Code		Consult Instructions for Use
	Serial Number		In Vitro diagnostic medical device
	Temperature limitations (storage temperature range)		Manufacturer

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Zoetis Inc.