













vetscan VS2

Vetscan VS2 testing options

Test 12 chemistry profiles with results in 12 minutes on one analyzer



Baseline Profiles	 Comprehensive Diagnostic Profile (CDP) ALB • ALP • ALT • AMY • BUN • Ca • CRE • GLOB* • GLU • K+ • Na+ • PHOS • TBIL • TP	Canine, feline, equine, ferret, rabbit
	 Preventive Care Profile Plus (PCPP) ALB • ALP • ALT • AST • BUN • Ca • Cl- • CRE • GLOB* • GLU • K+ • Na+ • TBIL • tCO ₂ • TP	Canine, feline, equine
	 Prep Profile II (PPII) ALP • ALT • BUN • CRE • GLU • TP	Canine, feline, equine
Preanesthetic + Critical Care Profiles	 Electrolyte Plus (EP) Cl- • K+ • Na+ • tCO ₂	Canine, feline, equine
	 Critical Care Plus (CCP) ALT • BUN • Cl- • CRE • GLU • K+ • Na+ • tCO ₂	Canine, feline, equine
Organ-specific Profiles	 Kidney Profile Plus (KPP) ALB • BUN • Ca • Cl- • CRE • GLU • K+ • Na+ • PHOS • tCO ₂	Canine, feline, equine
	 Mammalian Liver Profile (MLP) ALB • ALP • ALT • BA • BUN • CHOL • GGT • TBIL	Canine, feline, equine
Specialty Testing Chemistry Profiles	 T_i /Cholesterol Profile (T_i /CHOL) T _i • CHOL • GLU • K+ • Na+ • tCO ₂	Canine, feline
	 Phenobarbital Profile (PHB) ALB • ALP • ALT • AST • BUN • GGT • PHB • TBIL	Canine, feline
Species-specific chemistry profiles	 Avian/Reptilian Profile Plus (A/RPP) ALB • AST • BA • Ca • CK • GLOB* • GLU • K+ • Na+ • PHOS • TP • UA	Avian, reptile
	 Equine Profile Plus (EPP) ALB • AST • BUN • Ca • CK • CRE • GGT • GLOB* • GLU • K+ • Na+ • TBIL • tCO ₂ • TP	Equine
	 Large Animal Profile (LAP) ALB • ALP • AST • BUN • Ca • CK • GGT • GLOB* • Mg • PHOS • TP	Bovine

* Calculated value

LOOK DEEPER

Learn more at zoetisdiagnostics.com/VS2

