

ORDER FORM

Canine/Feline Most Common Tests

REQUIRED

Clinic Account #: _____ Collection Date: _____

Hospital/Clinic: _____ Zip Code: _____

Veterinarian: _____ Owner Last Name: _____

Animal Name: _____ Species: _____ Breed: _____

Age: _____ Sex: M M(N) F F(S) Patient ID/MRN: _____

Visit www.zoetisdx.com
to see our full directory
of tests or to order online.

Staff Pet: Please check if this is for an
employee's personal pet.*

ClinChem Panels

- 516 Pre-anesthetic ClinChem** 1 mL serum in SST or RTT
- CHEM18 ClinChem 18** 0.5 mL serum in SST or RTT
- 1000 ClinChem 25** 1 mL serum in SST or RTT
- 513 Mini Renal ClinChem** 1 mL serum in SST or RTT
- 481 Renal ClinChem** 1 mL serum in SST or RTT
- 479 Renal ClinChem + SDMA** 2 mL serum in SST or RTT
- 510 Liver ClinChem** 1 mL serum in SST or RTT
- 514 NSAID ClinChem** 1 mL serum in SST or RTT
- 488 Critical Care ClinChem** 1 mL serum in SST or RTT
- 512 Metabolic ClinChem** 2 mL serum in SST or RTT
- 515 Phenobarbital ClinChem** 2 mL serum in RTT
- 482 Senior Cat ClinChem** 1 mL serum in SST or RTT
- 507 Feline Hyperthyroid ClinChem** 2 mL serum in SST or RTT

Individual ClinChem Tests

- 7037 SDMA** 1 mL serum in SST or RTT
- 1023 Triglycerides** 0.5 mL serum in SST or RTT

Hematology

- 2000 Basic CBC** for wellness screens in non-geriatric, healthy patients
2 mL EDTA whole blood, air-dried slides
- CBCC Comprehensive CBC** recommended for sick or geriatric patients
2 mL EDTA whole blood, air-dried slides
- PATH Blood Smear Review** by a pathologist May be added on to
Comprehensive CBC only

GI/Pancreatic Testing

- 1030 DGGR Lipase** 0.5 mL serum in SST or RTT
- 6180 Bile Acids Pre/Resting** 0.5 mL serum in SST or RTT
- 6181 Bile Acids Post** 0.5 mL serum in SST or RTT
- 2073 Canine Pancreatic Lipase (cPL)** 0.5 mL serum in SST or RTT. Fast for
12-18hrs (6hrs for diabetic patients)
- 6206 TLI, Canine** 0.5 mL serum in SST or RTT. Fast for 12-18hrs (6hrs for
diabetic patients)
- 7253 PLI and TLI, Feline (TX GI)** 1 mL separated serum in SST or RTT. Fast for
12-18hrs (6hrs for diabetic patients)
- 6207 Cobalamin (Vitamin B12)** 1 mL separated serum in SST or RTT. Freeze
and protect sample from light
- 6208 Folate** 1 mL serum in SST or RTT. Freeze and protect sample from light
- KFCFOF Cobalamin and Folate** 2 mL serum in SST or RTT. Freeze and protect
sample from light
- 7250 cPL, TLI, Cobalamin, and Folate, Canine** 2 mL serum in SST or RTT. Fast
for 12-18hrs (6hrs for diabetic patients). Freeze and protect sample from light
- 7251 PLI, TLI, Cobalamin, and Folate, Feline (TX GI)** 1 mL separated serum in
SST or RTT. Fast for 12-18hrs (6hrs for diabetic patients). Freeze and protect
sample from light

Fecal/Parasitology

- 2320A Ova and Parasite (Centrifugal Fecal Flotation) Test** 3 gm fresh feces
- 2321 Giardia Antigen Test** 3 gm fresh feces
- 2320 Ova and Parasite (Centrifugal Fecal Flotation) Test and Giardia
Antigen Test** 5 gm of fresh feces, submitted in a fecal container
- 6001 Heartworm Antigen Test** 0.25 mL serum in SST or RTT or EDTA plasma
- 6225 Modified Knott's Test for Microfilariae (Heartworm)** 2 mL EDTA whole
blood

Bundles

If submitting a Bundle with urine, please indicate the source:

Sample Collection Type: Cysto Cath Free/Void

- 496 Pre-anesthetic Bundle** 1 mL serum in SST or RTT; 2 mL whole blood
in LTT, air-dried slides
- JWP Junior Wellness Bundle** 1 mL serum in SST or RTT; 2 mL EDTA whole
blood, air-dried slides; 5 mL urine in sterile container
- 501 Adult Basic Bundle** 2 mL serum, separated, in SST or RTT; 2 mL EDTA
whole blood, air-dried slides
- 502 Adult Basic + HW Bundle** 2 mL serum in SST or RTT; 2 mL EDTA whole
blood, air-dried slides
- AWP Adult Wellness Bundle** 1 mL serum in SST or RTT; 2 mL EDTA whole
blood, air-dried slides; 5 mL urine in sterile container
- SBASIC Senior Basic Bundle** 1 mL serum in SST or RTT; 2 mL EDTA whole
blood, air-dried slides
- SWP Senior Wellness Bundle** 1 mL serum in SST or RTT; 2 mL EDTA whole
blood, air-dried slides; 5 mL urine in sterile container
- 508 Critical Care Bundle** 1 mL serum in SST or RTT; 2 mL EDTA whole blood,
air-dried slides; 5 mL urine in sterile container
- 478 Renal Bundle** 1 mL serum in SST or RTT; 2 mL EDTA whole blood,
air-dried slides; 5 mL urine in sterile container
- 493 Metabolic Bundle** 1 mL serum in SST or RTT; 2 mL EDTA whole blood,
air-dried slides; 5 mL urine in sterile container
- 490 Feline Hyperthyroid Bundle** 2 mL serum in SST or RTT; 2 mL EDTA whole
blood, air-dried slides; 5 mL urine in sterile container
- 500 Senior Cat Bundle** 2 mL serum in SST or RTT; 2 mL EDTA whole blood,
air-dried slides; 5 mL urine in sterile container
- 492 Feline GI Bundle** 2 mL serum in SST or RTT, frozen and protected from
light, patient fasted 12-18 hours (6hrs for diabetic patients)
- 491 Canine GI Bundle** 2 mL serum in RTT, frozen and protected from light,
patient fasted 12-18 hours (6hrs for diabetic patients)

Endocrinology

- 5102 Total T4** 0.5 mL serum in SST or RTT
- 484 Thyroid Panel** 1 mL serum in SST or RTT
- 5101 Cortisol (Baseline/Resting)** 0.5 mL serum in SST or RTT
- 1024 Fructosamine** 0.5 mL serum in SST or RTT
- 5105 Progesterone** 0.5 mL serum in RTT

Serology/Infectious Diseases

- 6002 Tick Panel Screening Test** 0.25 mL serum in SST or RTT or EDTA plasma
or whole blood
- 6005 FeLV/FIV Rapid Test** 0.25 mL serum in SST or RTT or EDTA plasma or
whole blood

Urinalysis

Sample Collection Type: Cysto Cath Free/Void

- 3000 Urinalysis** 5 mL urine in sterile container
- 1555 Urine Protein: Creatinine Ratio** 1 mL urine in sterile container
- 521 Urinalysis with Reflex UPC** 5 mL urine in sterile container
- 538 Urinalysis with Aerobic Culture and Sensitivity** 5 mL urine in sterile
container and 1 mL urine in sterile container

Microbiology

- 7140 Urine Culture and Sensitivity** 1 mL urine in WTT or sterile tube

Visit www.zoetisdx.com to see our full directory of tests or to order online. For questions, or to order by phone, please call 1-888-965-9652.

*Subject to program guidelines

WHITE COPY: LAB | YELLOW COPY: CUSTOMER

All trademarks are the property of Zoetis Services LLC or a related company or a licensor unless otherwise noted.

© 2023 Zoetis Services LLC. All rights reserved. ZRL-00299 R1