

ZOETIS DIAGNOSTICS

vetscan Rapid™ Flex4

Quickly and easily screen for vector-borne diseases in dogs

Make in-house diagnostics simple with this convenient, easy-to-use visual rapid test.



Detects antibodies to *Anaplasma phagocytophilum*, *A. platys*, *Borrelia burgdorferi*, *Ehrlichia canis*, *E. chaffeensis*, *E. ewingii*, and *Dirofilaria immitis* antigen



Reliable, point-of-care rapid test **help make the diagnoses**

- ✓ Highly accurate¹⁻⁴ results in as little as **8 minutes**
- ✓ Delivered through **lateral flow technology**
- ✓ Distinguishes between **Lyme vaccination and Lyme exposure⁵**



Easy operation **requires minimal hands-on time**

- ✓ **Simple, intuitive 3-step process** helps minimize errors
- ✓ **Automatic reading and reporting** with the Vetscan Vue®
- ✓ **Frees staff** to spend time where they're needed most



Efficient operation and storage **minimizes waste**

- ✓ 12-month shelf life helps **ensure every test counts**
- ✓ Room temperature storage **saves valuable refrigerator space**

¹ Additional Study: Abaxis Anaplasma Antibody Test Kit Sensitivity and Specificity Report

² Additional Study: Canine Lyme Disease S&S Clinical Trial Submission Report

³ Additional Study: Abaxis Ehrlichia Sensitivity and Specificity Report to USDA

⁴ Henry, LG, et al., Vet Parasitol, 2018

⁵ Study Report No. DHX6Z-US-18-051

LOOK DEEPER

Learn more at zoetisdiagnostics.com/rapids

zoetis

ZOETIS DIAGNOSTICS

vetscan Rapid™ Flex4

Results you can trust

The Vetscan Flex4 Rapid Test offers highly accurate¹⁻⁴ results through lateral flow technology, a reliable tool in veterinary and human medicine for efficient diagnosis and patient evaluation for more than 30 years.



	Sensitivity	Specificity
Ehrlichia	97.4% ¹	97.5% ¹
Lyme	100% ²	100% ²
Anaplasma	93.3% ³	96.4% ³
Heartworm	98.5% ⁴	94% ⁴

Deeper diagnostic insights that elevate care

Vetscan rapid tests are part of an ecosystem of innovative diagnostic products and services, all seamlessly connected by the Zoetis Virtual Laboratory.

Comprehensive

Combine rapid tests with other Virtual Laboratory products and services to provide a complete picture that helps you individualize care.

Connected

Connect with ZoetisDx.com, and view results online when you use the Vetscan Vue® to scan, interpret and share results.

Supported

Schedule free consultations with a specialist at ZoetisDx.com, or request confirmatory testing on certain tests online or contact Dx tech support.

* Heartworm testing is always available, completely complimentary. Confirmation testing for Lyme, Anaplasma, and Ehrlichia results is on a case-by-case basis, and additional charges may apply.

¹ Additional Study: Abaxis Anaplasma Antibody Test Kit Sensitivity and Specificity Report

² Additional Study: Canine Lyme Disease S&S Clinical Trial Submission Report

³ Additional Study: Abaxis Ehrlichia Sensitivity and Specificity Report to USDA

⁴ Henry, LG, et al., Vet Parasitol, 2018

LOOK DEEPER

Learn more at zoetisdiagnostics.com/rapids

zoetis