

ZOETIS DIAGNOSTICS

reference
laboratories
by *zoetis*

Directory of Services

2025





Dear Valued Veterinary Partner,

At Zoetis Reference Laboratories, we understand the unique challenges you face in today's veterinary landscape — from the increasing demand for care coupled with workforce shortages to the daily stresses of treating ill and injured animals — and we are here to support you.

Our mission is to help you elevate your practice by providing a laboratory and diagnostic partnership that ensures you have the tests you need, with consistent, accurate and timely results. Every decision we make is centered around your clinic and your patients, and our medical roster is deep with expertise and committed to assisting you at every step.

We continually strive for excellence in all that we do, even incorporating gold standards from human laboratory diagnostic medicine into our Quality program. Our dedication to quality and continuous improvement supports you in your goal of providing the highest level of care to your patients.

Together, we can navigate the challenges of modern veterinary medicine and deliver exceptional care to the animals that depend on us.

Sincerely,

Kavita Bhat

Kavita Bhat
Vice President of US Zoetis Reference Laboratories, Global Diagnostics

Quality is Key

Expect excellence from our rigorous Quality Management System (QMS)

The Zoetis Reference Laboratories (ZRL) Quality Management System incorporates American Society for Veterinary Clinical Pathology (ASVCP) quality standards¹ as best practice but doesn't stop there. ZRL goes the extra mile by incorporating gold standards from human laboratory diagnostic medicine into our Quality program, including Clinical and Laboratory Standards Institute (CLSI), ISO 15189 and College of American Pathologists (CAP). From timely test results to technical oversight to comprehensive tracking of metrics, we continually strive for excellence in all that we do.



Rigorous quality standards upheld in all locations and departments and through all phases of testing

Our commitment to quality and continuous improvement supports veterinarians in providing a high level of care to their patients. The rigorous ZRL quality program includes:

- Dedicated, experienced Quality team with National and Regional leadership
- Medical leadership, ensuring both service excellence and focus on patient safety and well-being
- Rigorous analytical performance oversight, including robust validation, calibration, Quality Control (QC) and Proficiency Testing (PT)
- Comprehensive staff training and competency program
- Semi-annual internal auditing program
- Comprehensive and proactive management of non-conforming events
- Monitoring and incorporation of customer feedback for improvement projects
- Continuous improvement program, incorporating Lean and Six Sigma² principles and tools

The ZRL laboratory team is assembled to ensure quality

Our diligent, well-rounded and diverse laboratory team includes:

- Board-Certified (DACVP) Clinical and Anatomic Pathologists
- Board-Certified (DACVP) Parasitologist
- Clinical Microbiologist (PhD)
- Board-Certified Specialty Consultants
- Extensively credentialed Quality leadership
- Qualified Technicians (MLS/MT/CLS, RVT/CVT and BS/MS)

References: 1. American Society of Veterinary Clinical Pathology. Quality Assurance and Laboratory Standards (QALS) Guidelines. https://www.asvcp.org/page/QALS_Guidelines. Accessed September 30, 2024.
2. American Society for Quality. What is Lean Six Sigma? <https://asq.org/quality-resources/six-sigma>. Accessed September 30, 2024.

A Guide to This Directory

Welcome to the 2025 Directory of Services for Zoetis Reference Laboratories!

To ensure our Directory is as helpful as possible, we've rearranged it to make the information you need more easily accessible, including:

- A simplified look at [turnaround times](#) that includes sample collection and submission guidelines
- Information to help you [determine which of our test options best fits your patient's needs](#)
- How you can easily [order tests and view the results digitally](#) anytime, anywhere
- How to [order supplies](#)

Refer to the [Table of Contents](#) for a comprehensive list of the resources available to you in this guide.

Using the interactive features

- Use the Table of Contents or the navigation bar across the bottom of each page to jump directly to sections.
- Click test codes throughout the Directory to navigate to that test.
- **Adobe Acrobat Reader® users:** Click KEY in the upper right corner throughout the Test Menu section to access the table key.
All other users: View the [key](#) here.
- Not sure where to find a test? Locate it within the alphabetical test listing at the back of this Directory, then click the test code to navigate to it.

Need extra support?

We're happy to help. Please feel free to contact us using the numbers listed in our [Dedicated Customer Service](#) section.

Contents

Select a section to jump directly to that page.

Everything you need to know, from our next-day turnaround times to how to choose a test, prepare and submit your sample, and access your results.

A listing of our comprehensive test menu, organized by species, then laboratory section, then by test name.

.....
An overview of our policies pertaining to species, add-on tests, sample retention, payment and processing fees, and more.
.....

An alphabetical listing of our tests by species.

About Zoetis Reference Laboratories

Your Laboratory Testing Partner

Zoetis Reference Laboratories gives you deeply informed diagnostic insights with a comprehensive test menu and reliable turnaround times, backed by expert support, all easily accessible through innovative connectivity solutions.



Diagnose with **expert medical insights**

- ✓ Dedicated anatomic and clinical pathologists
- ✓ Comprehensive test menu with next-day results
- ✓ Accurate results supported by specialist insights for enhanced patient outcomes



Tailor your diagnostic approach with **flexible solutions**

- ✓ Fully customizable diagnostic workflows, from test to treatment selection
- ✓ Seamlessly integrated with Vetscan Imagyst® and other in-clinic analyzers
- ✓ End-to-end ordering and support managed via your PIMS and the ZoetisDx portal



Optimize your point-of-care with **high-value diagnostics**

- ✓ Maximize efficiency by integrating Zoetis Reference Laboratories into your practice workflow
- ✓ Standardize monthly costs with tiered pricing options
- ✓ Expand your diagnostics and grow your practice with My Zoetis Rewards

Test on Your Terms

From the clinic to the reference lab and everything in between, we've got you covered.

Our reference laboratory service is always there as foundational support for your diagnostics practice, with quality, next-day results you can rely on. It's just one part of an innovative Zoetis Diagnostics offering that enables a uniquely integrated, fully customizable workflow for your practice and your patients.

This powerful portfolio of products and services gives you the tools to confidently diagnose and treat any case in the way that best suits your patients' needs. And it's all supported by the Zoetis Virtual Laboratory — an integrated support network of board-certified specialists paired with diagnostic insights that enhances every element of your diagnostic practice.

Wherever and however you test, you'll get the deep insights you need to elevate patient care.



Comprehensive Diagnostics

Full-service reference laboratories and breakthrough point-of-care capabilities for cytology, hematology, fecal testing and more.



Expert Support

Convenient expert review[†] and complimentary specialist consultations on any case via Zoom or email.



Connected Insights

Anytime, anywhere access to digital ordering, test results, expert support and more on the ZoetisDx portal — conveniently integrated with your PIMS*.

* Compatible with select practice information management software (PIMS).

[†] Additional costs may apply.

Turnaround Times

Fast, high-quality results are more important than ever in veterinary diagnostics. Our turnaround times reflect this, helping to ensure you have the results you need when you need them — easily accessed through your ZoetisDx portal — so you can make informed treatment decisions as quickly as possible.

Most common tests performed in our laboratories provide next-day results. Send-out testing, cultures and other tests may take longer. For specific test information, find your desired test in the test menu.

Routine Testing

Test Type	Results Delivery Turnaround Time	Sample Collection and Submission Guidelines For seasonal shipment tips, click here .
Chemistry (e.g., ClinChem 25, Total T4, SDMA)	Next day by 8 AM	<ul style="list-style-type: none"> - Fasted samples are recommended to avoid lipemia - Collect in SST or RTT - If using a RTT, separate serum from the cellular component within 1 hr of collection - Follow centrifugation best practices - Label aliquoted serum sample as SERUM
Hematology (e.g., CBC)	Next day by 8 AM	<ul style="list-style-type: none"> - Fasted samples are recommended to avoid lipemia - Use the largest vein and largest needle that can safely be used - Collect in LTT (CBC) or BTT (most coagulation tests) - Fill tubes to the blood tube's volume fill line. Failure to do so may clinically affect results. - Label aliquoted plasma as CITRATED PLASMA, if relevant - Coagulation testing requires special handling and separation of citrated plasma immediately after collection
Infectious Disease Serology (e.g., Heartworm, Lyme)	Next day by 8 AM	<ul style="list-style-type: none"> - May use whole blood, plasma or serum, as indicated by specific test ordered - Fasted samples are recommended to avoid lipemia - Follow separation/centrifugation best practices - Label aliquoted serum or plasma as such, including the type of plasma (e.g., EDTA PLASMA or CITRATED PLASMA)
Urinalysis	Next day by 8 AM	<ul style="list-style-type: none"> - Recommend submission of ≥ 5 ml of urine for most accurate results - Urine volume submissions < 1 mL may result in incomplete urinalysis testing - Collect in a sterile, single-use yellow top urine tube with no additive or sterile urine collection cup - Label the tube URINE
Fecal Parasite Tests (e.g., Ova & Parasites, Giardia Antigen Test)	Next day by 5 PM	<ul style="list-style-type: none"> - 3-5 grams of fresh feces in sealed fecal cup

For general guidelines and additional shipping requirements, see [Sample Guidelines & Submission](#) for the type of test you are ordering.

[Core Lab Specimens](#) | [Parasitology Specimens](#) | [Microbiology Specimens](#) | [Histopathology Specimens](#) | [Cytology Specimens](#)

Turnaround Times (cont'd)

Pathology

Test Type	Results Delivery Turnaround Time	Sample Collection and Submission Guidelines For seasonal shipment tips, click here .
Cytopathology	1-2 days, Mon-Fri	<ul style="list-style-type: none"> - Use new, clean glass slides with frosted edge - Label slides clearly with patient name and site sampled using a pencil on frosted edge - Collect sample using the appropriate method for the lesion type; FNA is most common - Fluid submissions for cytology should be collected into an EDTA tube and refrigerated until sample pickup. - Glass slides should NOT be refrigerated.
Histopathology	<p>Simple 3-7 days, Mon-Fri</p> <p>Complex Level 1 3-9 days, Mon-Fri</p> <p>Complex Level 2 Varies by nature of sample</p>	<ul style="list-style-type: none"> - Formalin fixed or fresh tissue - History and description of the lesion is essential to a quality interpretation - Indicate the number of containers and specimens being sent - Label containers with the source of specimen - Formalin fixed tissue may be stored at room temperature - Keep fresh tissue refrigerated until sample pickup (do not freeze)

Turnaround Time Policy

Our goal is to meet published turnaround times and ensure results are delivered as quickly as possible.

- Turnaround times for routine testing do not include Sunday.
- Turnaround times for pathology testing do not include Saturday and Sunday.
- Electronic orders are preferred. Manual Paper requisition forms may extend turnaround time.
- If the quality of a sample we receive prevents us from producing quality results within our promised turnaround time, we will notify you. Examples include samples that are damaged, contaminated, incorrectly preserved or packed, and samples that do not meet volume requirements.
- Once a specimen is received at the lab, we assume any risk of loss or damage.
- We reserve the right to reject any specimen likely to pose an unreasonable health or safety risk to our lab staff.

Send-Out Testing


Some tests in our Directory are not performed at Zoetis Reference Laboratories. Turnaround times for these tests are dependent on the partner laboratory that processes these tests.

Zoetis Reference Labs Testing Options: Hematology

Choose the CBC test that's best for your patients.

At Zoetis Reference Laboratories, our CBCs are grounded in quality standard operating procedures developed by a team of board-certified clinical pathologists. Expertly trained technicians produce quality results, collaborating directly with pathologists as needed along the way.

	Basic Test code: 2000	Comprehensive Test code: CBCC
	Recommended for wellness screens in non-geriatric, healthy patients	Recommended for sick or geriatric patients
WBC, HGB, RBC, HCT, MCV, MCHC, CHCM, RDW, 5-Part WBC Differential, Platelet Count, MPV	✓	✓
Performed if indicated based on automated results: Spun PCV, Manual WBC Differential, Assessment for Platelet Clumping +/-, Visual Platelet Estimate, and/or Reticulocyte Count	✓	✓
Comprehensive manual smear review by technician		✓
Pathologist review if indicated by laboratory established criteria		✓



See the full list of Hematology Tests for...

[Canine & Feline](#) | [Large Animals](#) | [Avian, Exotic, & Small Mammal](#)

Test requirements: EDTA whole blood, air-dried slides


Zoetis Reference Labs Testing Options: Fecal/Parasitology

Solutions tailored to your patients' medical needs

- Ova and Parasite (Centrifugal Fecal Flotation) Test | **2320A**
- *Giardia* Antigen Test | **2321**
- *Giardia* Antigen Test, Canine/Feline Add-on | **A2321**
- Ova and Parasite (Centrifugal Fecal Flotation) Test and *Giardia* Antigen Test | **2320**
- *Giardia*-Focused O&P | **25001611**
- Fecal Direct Smear Test | **25001610**
- Fecal Sedimentation Test | **25001612**
- Baermann Test | **BAER**
- O&P (Centrifugal Fecal Flotation) Panel | **25001613**
- O&P (Centrifugal Fecal Flotation) with Direct Smear | **25001614**

See the full list of Fecal/Parasitology Tests for...

[Canine & Feline](#) | [Large Animals](#)



* If *Giardia* is confirmed and treated, follow-up testing 24-48 hours after completion of treatment is recommended.

† *Giardia*-focused O&P is recommended for both asymptomatic and symptomatic patients. If no ova or parasites are seen in symptomatic patients, recommend *Giardia* Ag test

Zoetis Reference Labs Testing Options: Cytopathology

Solutions tailored to your medical need

Our standard turnaround time is 1-2 days, Mon-Fri.



Cytopathology Simple/Single Site | 7120

For the following sample types:

- Single lesion/tissue (FNA, swab, impression smear)
- Single or multiple lymph nodes
- Single joint (slides only)
- Fluid from mass lesion (cytology only)



Cytopathology Complex | 7201

For the following sample types:

- Multiple lesions/tissues (FNA, swab, impression smear) - up to 7 sites
- Joint fluid in LTT (single or multiple)
- Multiple joints (slides only)
- Cavitory effusion (single or multiple)
- CSF
- Bone marrow
- BAL/TTW
- Other fluids/lavage

See the full list of Cytopathology Tests for...


[Canine & Feline](#) | [Large Animals](#)

Learn more about the clinical pathologists supporting you and your patients

[Meet the Experts](#)

Zoetis Reference Labs Testing Options: Histopathology

Is your sample simple or complex? Find solutions tailored to your needs



Multiple sample combinations: Which test code is right for you?

Number and Type of Samples Included	Test Selection
2-5 Simple Sites	Simple Multiple Sites ZN101
6-15 Simple Sites	Complex Level 1 ZN102
>15 Simple Sites	Complex Level 2 ZN103
2+ Complex Level 1 Sites	Complex Level 2 ZN103
Complex Level 1 Site & <5 Simple Sites	Complex Level 1 ZN102
Complex Level 1 Site & >5 Simple Sites	Complex Level 2 ZN103

Specialty Panels

- Dermatopathology Panel | **ZN400**
- Liver Panel | **ZN500**
- Liver Panel & Medical Consult | **25001436**
- Renal Panel | **ZN25000541**

See the full list of Histopathology Tests for...

Canine & Feline | **Large Animals**

Learn more about the anatomic pathologists supporting you and your patients

Meet the Experts

Zoetis Reference Labs Testing Options: Histopathology (cont'd)

Synoptic reporting designed to improve interpretation

We're pleased to tell you about our histopathology report format for malignancies and gastrointestinal biopsies.

This report structure, known as synoptic reporting, breaks down results into easy-to-read tables and checklists. By simplifying data presentation and ensuring consistency of formatting, we're making it easier for you to:






- Locate and evaluate key diagnostic and prognostic features
- Scan and digest report findings
- Communicate effectively with pathologists






Tube Guide

Use this quick reference guide to identify which tubes are appropriate for sample collection, based on sample type and common test usage.

For instructions on how to prepare and submit your sample, please refer to the collection guidelines specific to your sample type.

Tube Type	Additive	Sample Type(s)	Tube Description	Common Uses
Blue-top tube ● BTT 	3.2% sodium citrate anticoagulant	Plasma	Light blue top – Microtube also available for low-volume samples (very small, pediatric, and exotic species)	Coagulation testing (PT, aPTT, coagulation factor assays)
Purple-top or lavender-top tube ● LTT 	Potassium EDTA anticoagulant	Whole blood or plasma (EDTA plasma is rarely required)	Purple or lavender top – Microtube also available for low-volume samples (very small, pediatric, and exotic species)	Most common use: – Whole blood for CBC
Green-top tube ● GRNTT 	Lithium heparin	Whole blood or plasma	Green top – Microtube also available for low-volume samples (very small, pediatric, and exotic species)	Most common use: – Free (ionized) calcium, whole blood for point-of-care chemistry testing Other uses: – Lithium heparin plasma sometimes required for specialized testing
Mint- or light green-top microtube ■ M-GRNTT 	Lithium heparin	Plasma	Mint- or light green-top microtube – Opaque, light yellow gel at bottom of tube – Differs from green-top tube in that it contains a gel barrier that separates RBCs from plasma when centrifuged	Most common use: – Lithium heparin plasma sometimes required for specialized testing
Red-top tube ● RTT 	Clot activator	Serum	Red top with chalky residue (clot activator) in tube Please note: No-additive RTTs are also available. Tubes with clot activator should not be used when a no-additive tube is required.	Most common use: – Serum chemistry testing Other uses: – Immunology – Serology – Many others

Tube Guide (cont'd)

Tube Type	Additive	Sample Type(s)	Tube Description	Common Uses
Serum-separator tube or tiger-top tube ○ SST 	Clot activator (embedded in gel)	Serum	Red top, red and black marble top (tiger top) <ul style="list-style-type: none"> - Pediatric tube with gold-top microtube also available for low-volume samples (very small, pediatric, and exotic species) - Important not to confuse gold-top tube with yellow-top urine tube - All tubes have opaque light yellow gel at bottom of tube 	Most common use: <ul style="list-style-type: none"> - Serum chemistry testing Other uses: <ul style="list-style-type: none"> - Immunology - Serology - Many others
White-top tube ○ WTT 	No additive	Serum or Plasma	White top (no additive)	Used when serum or plasma needs to be separated from RBCs and removed from original collection tube to avoid prolonged contact with RBCs or separator gel <ul style="list-style-type: none"> - Used when special handling of sample required (such as freezing) - Sometimes used for collection of fluids/urine and materials for culture
Yellow-top tube ● YTT 	No additive	Urine	Bright yellow top <ul style="list-style-type: none"> - May have a conical or tapered base - Do not confuse with gold-top serum tubes used occasionally 	Urine for urinalysis and other various tests

Sample Guidelines & Submission

Core Lab Specimens

Core Lab specimens include whole blood, serum, plasma, urine, other liquid samples and blood smears (slides).

- Label all tubes with the patient's first and last name.
- Serum, plasma or sometimes urine cannot be differentiated visibly. Be sure to label the tube with the type of sample it contains. Plasma must be further differentiated as containing EDTA, Citrate or Heparin plasma.
- Refrigerate within 30 minutes of collection and keep refrigerated until sample pickup.
- Do not freeze samples unless expressly indicated in the test menu.
- If shipping refrigerated samples via UPS, include a frozen ice pack in the packaging. Wrap the ice pack in a paper towel or similar covering to avoid direct contact with the sample, which may cause the sample to freeze.
- For CBC submissions, we recommend making a blood smear at the time of sample collection. Store any prepared slides at room temperature. Refrigeration of glass slides causes condensation that can destroy the sample. Submit any prepared slides along with the sample tubes.
- Samples that must be shipped frozen require additional special handling. Consult the [Test Menu](#) or contact [Customer Service](#) for specific shipping recommendations.

See the [full blood sample collection and submission guidelines](#), which includes best practices for whole blood, plasma, serum and urine.

Parasitology Specimens

Parasitology specimens include feces.

- Refrigerate as soon as possible after collection and keep refrigerated until sample pickup.
- Include a frozen ice pack in the packaging.
- Wrap the ice pack in a paper towel or similar covering to avoid direct contact with the sample, which may cause the sample to freeze.

Microbiology Specimens

Microbiology specimens include culture swabs, liquid samples (e.g., urine), tissue and feces.

- Collect specimens as aseptically as possible to prevent overgrowth by normal flora.
- Collect before the administration of microbials, whenever possible.
- Label sample with the patient's first and last name, the source/site, and date of collection.
- Store and ship culture swabs at room temperature.
- Refrigerate urine, feces or tissue being submitted for culture within 30 minutes of collection, and keep refrigerated until sample pickup.
- If shipping refrigerated samples via UPS, include a frozen ice pack in the packaging. Wrap the ice pack in a paper towel or similar covering to avoid direct contact with the sample, which may cause the sample to freeze.
- Do not send needles, razors, broken glass, plastic or anything sharp with your sample or as your sample.
- Please properly discard used needles or other sharps prior to transportation. For the submission of syringes, the needle must be removed and the syringe capped prior to transportation to ensure there is no leakage.

See the [full microbiology sample collection and submission guidelines](#).

Sample Guidelines & Submission (cont'd)

Histopathology Specimens

Histopathology specimens include biopsies — formalin-fixed and fresh tissue.

- Complete the requisition form as indicated for pathology submissions. History and description of the lesion are vital to a quality histopathology report.
- Indicate the number of containers and specimens being sent so the lab can appropriately verify receipt of all samples.
- Label the container(s) with the patient's first and last name, collection date, and source of specimen. Ensure this labeling matches information supplied on the requisition form.
- Fresh tissue must be kept cold with bagged ice or ice packs during shipment to avoid tissue necrosis. See below for additional guidance.
- Place the specimen and any forms in separate, sealed plastic bags. Double bag specimens, whenever possible.
- Ice packs and bags can break or leak during shipping, which is why we stress sealing accompanying forms in separate, sealed plastic bags.

Formalin-fixed tissue

- Always place tissues in 10% neutral buffered formalin to ensure appropriate fixation. A 10:1 formalin to tissue ratio is recommended.
- In extreme cold, add isopropyl or ethyl alcohol to your formalin jars (1 part alcohol to 10 parts formalin) to prevent freezing. Note "alcohol added" on the container or in the submission notes.
- In extreme heat, submit samples in an insulated shipping container with frozen ice packs.

Fresh tissue

- In extreme cold, submit samples in an insulated shipping container and include cold (not frozen) soft gel ice packs.
- In extreme heat, submit samples in an insulated shipping container with frozen ice packs.

See the [full histopathology sample collection and submission guidelines](#).

Cytology Specimens

Cytology specimens include customer-prepared slides and liquid samples.

- Use new, clean glass slides with a frosted edge.
- Using a pencil (ink labeling can wash off), label slides clearly on the frosted edge with the patient name and sites.
- Store any prepared slides at room temperature. Refrigeration of glass slides causes condensation that can destroy the sample.
- For shipping cytology slides, we recommend using hard plastic slide boxes designed for this purpose. Ensure all slides have completely air-dried before packaging. DO NOT enclose slides in slide boxes prior to drying, as this can damage the sample due to humidity.
- Avoid shipping any cytology slides or liquid samples inside packaging that also contains any form of formalin. Fumes can escape from closed formalin containers or liquid formalin can leak en route. Both formalin fumes and liquid formalin can destroy a cytology slide.
- Refrigerate any liquid samples (e.g., cavitory effusion, CSF) until sample pickup.
- If liquid samples are being shipped for cytology, include a frozen ice pack in the packaging. Wrap the ice pack in a paper towel or similar covering to avoid direct contact with the sample, which may cause the sample to freeze.
- If a liquid cytology sample is intended for cytological examination, we recommend using an EDTA tube. If a liquid cytology sample is also intended for culture, submit an additional sterile, no-additive tube aliquot.
- We recommend preparing a direct fluid smear at the time of liquid specimen sample collection. Submit any prepared slides with the sample tubes of liquid samples.

See the [full cytopathology sample collection and submission guidelines](#).

Seasonal Shipment Tips

Tips for Summer

Extreme heat can compromise sample integrity. Follow these guidelines to help maintain sample quality.



Do not place samples in lockboxes or outdoors until close to sample pickup.



Include ice packs for any refrigerated samples and consult the specimen instructions for any frozen sample shipping requirements (dry ice may be required).



Multiple ice packs may be utilized when temperatures are 90°F and above.



Wrap the ice pack in a paper towel or similar covering to avoid direct contact with the sample, which may cause the sample to freeze.

Tips for Winter

Extreme cold can compromise sample integrity. Unintentional freezing of samples can render them useless. Follow these guidelines to help maintain sample quality.



Do not place samples in lockboxes or outdoors until close to sample pickup.



Do not use a frozen ice pack in cold weather. Instead, use a cold (not frozen) soft gel ice pack, placing the pack around the samples. This will help ensure temperature distribution during transit, decreasing the likelihood of sample freezing.



In extremely cold weather consider using a Styrofoam-insulated box for shipping.

Sending Your Sample

Our goal is to offer you the best reference lab experience. To make sure sending samples to us is as convenient as possible, we are covering the cost of transportation fees for qualified customers.*

Next-day Service with UPS®



Please follow the instructions to the right, but feel free to give us a call at **1-888-965-9652** or reach out to your local representatives if you have any questions.

Simple Steps for Sending a Package Via UPS:

1. Request a UPS label

- Call Customer Service at **1-888-965-9652** to create and print a shipping label. We will provide you with UPS shipping labels that are addressed to Zoetis Reference Laboratories.
- You must use these UPS labels on all packages you send to Zoetis Reference Laboratories. This will ensure that your packages are delivered to us as quickly and as safely as possible.

2. Prepare the package

- Put your labeled and packaged samples in the UPS Laboratory Pak and seal.
- We have insulated mailers available for shipping in extreme temperatures. If you are using these, place everything, including ice packs, inside the insulated mailer, place the insulated mailer inside the UPS Laboratory Pak, and seal.
- Peel and affix the UPS label. Multiple samples can go in one UPS Laboratory Pak.
- Keep the separate, smaller “customer copy” label that contains the UPS tracking number for your records. If any issues arise during shipping, having the UPS tracking number is essential when trying to locate your package.
- For seasonal tips, [click here](#).

3. Set up your clinic as a UPS Package Pickup Point

- We recommend setting up your clinic as a UPS Package Pickup Point (UPS-PPP). This is a free service that will be completed for you during the onboarding process.
- If you anticipate sending samples more than 3 times a week, scheduling a daily UPS pickup is recommended.
- Once your package pickup point is activated and a schedule established, the driver will stop by daily within the assigned time frame, with no call or email required. (Note: The UPS driver will pick up any UPS packages at your location.)

* \$50.00 monthly Processing and Handling fee applied to customers spending \$500 or less per month.

To find out if courier service is available for your practice, please contact your Zoetis Diagnostics representative or call us at **1-888-965-9652**.

Access Your Results Anytime, Anywhere, from Any Device

A Zoetis-connected practice means you can access the power of a comprehensive and fully integrated suite of Zoetis products and services from your practice management software, your phone and everywhere in between.

The ZoetisDx portal is the connectivity solution from Zoetis that makes it easy to:

- Order and view test results digitally, eliminating the need for paper transactions
- Access the most up-to-date test menu
- Drive informative discussions to improve pet owner-veterinarian communication
- Connect with a global network of board-certified specialists to support your clinical decisions
- Access Zoetis Reference Laboratories and Zoetis point-of-care results on a single, easy-to-use platform

Ready to get Zoetis-connected?

Visit the [ZoetisDx portal](#) to sign up today.



The ZoetisDx portal keeps you
Zoetis-connected

Customer Service

Dedicated Customer Service

Get help your way!

Login to zoetisdx.com

1-888-965-9652

8:00 AM–9:30 PM ET (Weekdays)

9:00 AM–8:30 PM ET (Saturdays)

referencelabsupport@zoetis.com


On-demand Solutions

Visit zoetisreflabs.com/customers for these self-service options and more:

- Order supplies
- Print shipping labels
- Register for the ZoetisDx portal

Easy Supply Ordering

We are happy to provide routine laboratory supplies at no charge to our active Zoetis Reference Laboratories customers. To order supplies, please call Customer Service at **1-888-965-9652** or access eCommerce at shop.zoetis.com for all your supplies.



Request your complimentary supplies through customer service or online at shop.zoetis.com

Test Menu

Key

TAT: Turnaround Time

○ **SST:** Serum Separator (gel)

● **RTT:** Clot Activator (no gel)

● **LTT:** K EDTA

● **BTT:** Na Citrate

● **GRNTT:** Lithium Heparin

○ **WTT:** No Additive

● **Feces**

● **Urine**

● **GRYTT:** K Oxalate/Na Fluoride

● **RBTT:** Na EDTA/Na Heparin

† Send-out test

🧊 Refrigerate

🧊🚫 Refrigerate, do not freeze


🧊💡 Refrigerate, protect from light

❄️ Freeze

❄️💡 Freeze, protect from light

🏠 Room temperature

For additional information, see the [Tube Guide](#)





Canine/ Feline



Bundle Overview

Our bundle menu offers comprehensive options organized into two categories:



Wellness Bundles

Named for **life stages** and generally intended for wellness testing:

- ZRL™ Junior
- ZRL™ Adult
- ZRL™ Senior



Clinical Bundles

Named for **clinical area of need** and generally intended for evaluation of specific organ systems or investigation of specific disease states.



Upgrade any bundle to a Comprehensive CBC. Use test code **BCCUP**.

Upgrade any bundle with a Urinalysis to a Urinalysis with Reflex UPC. Use test code **25001702**.

Junior and Adult Bundles Overview



Upgrade any bundle to a Comprehensive CBC. Use test code **CBCCUP**.

Upgrade any bundle with a Urinalysis to a Urinalysis with Reflex UPC. Use test code **25001702**.

Junior Bundles

All Junior bundles include a Basic CBC and a ClinChem 18 with SDMA.

Select each code for full test details.

Components	Test Name	Test Code	Price
Basic CBC, ClinChem 18 with SDMA, Urinalysis	ZRL™ Junior Complete	JWP	\$51.01
+ HW	ZRL™ Junior Complete with HW	25001642	\$57.58
Basic CBC, ClinChem 18 with SDMA, Ova and Parasite	ZRL™ Junior with O&P	25001635	\$71.97
+ HW, Tick	ZRL™ Junior with O&P, HW, Tick	25001636	\$80.03
+ Giardia, HW, Tick	ZRL™ Junior with O&P, Giardia, HW, Tick	25001637	\$82.66
Basic CBC, ClinChem 18 with SDMA	ZRL™ Junior	25001634	\$48.65
+ T4	ZRL™ Junior with T4	25001641	\$64.85
+ HW	ZRL™ Junior with HW	25001638	\$54.86
+ HW, Tick	ZRL™ Junior with HW, Tick	25001639	\$66.85
+ PT/aPTT	ZRL™ Junior with PT/aPTT	25001640	\$61.86

Adult Bundles

All Adult bundles include a Basic CBC and a ClinChem 25 with SDMA.

Components	Test Name	Test Code	Price
Basic CBC, ClinChem 25 with SDMA, Urinalysis	ZRL™ Adult Complete	AWP	\$50.74
+ O&P, HW, Tick	ZRL™ Adult Complete with O&P, HW, Tick	25001631	\$75.56
Basic CBC, ClinChem 25 with SDMA, Ova and Parasite	ZRL™ Adult with O&P	25001630	\$66.63
Basic CBC, ClinChem 25 with SDMA	ZRL™ Adult	25000501	\$43.64
+ HW	ZRL™ Adult with HW	25000502	\$50.21

Senior Bundles Overview

Senior Bundles

All Senior bundles include a Basic CBC, a ClinChem 25 with SDMA, and a T4.



Comprehensive CBC upgrade is recommended for all Senior bundles
Use test code **CBCCUP**.

Upgrade any bundle with a Urinalysis to a Urinalysis with Reflex UPC.
Use test code **25001702**.

Select each code for full test details.

Components	Test Name	Test Code	Price
Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis	ZRL™ Senior Complete	SWP	\$74.89
+ FT4	ZRL™ Senior Complete with FT4	25001656	\$87.24
+ HW, TG	ZRL™ Senior Complete with HW, TG	25001657	\$86.13
+ ZRL Premium Lipase (DGGR Lipase), TG	ZRL™ Senior Complete with ZRL Premium Lipase (DGGR Lipase), TG	25001655	\$82.15
+ O&P, FeLV/FIV	ZRL™ Senior Complete with O&P, FeLV/FIV	25001646	\$89.86
+ O&P, FeLV/FIV, HW	ZRL™ Senior Complete with O&P, FeLV/FIV, HW	25001647	\$103.31
+ O&P, Giardia	ZRL™ Senior Complete with O&P, Giardia	25001650	\$86.13
+ O&P, Giardia, HW, Tick	ZRL™ Senior Complete with O&P, Giardia, HW, Tick	25001651	\$103.17
+ O&P, HW, Tick, ZRL Premium Lipase (DGGR Lipase), TG	ZRL™ Senior Complete with O&P, HW, Tick, ZRL Premium Lipase (DGGR Lipase), TG	25001653	\$111.23
+ O&P, HW	ZRL™ Senior Complete with O&P, HW	25001648	\$113.32
+ O&P, HW, Tick	ZRL™ Senior Complete with O&P, HW, Tick	25001652	\$101.25
+ O&P, HW, Tick, UPC	ZRL™ Senior Complete with O&P, HW, Tick, UPC	25001654	\$113.96
+ O&P, fPLI	ZRL™ Senior Complete with O&P, fPLI	25001649	\$140.93
Basic CBC, ClinChem 25 with SDMA, T4, O&P	ZRL™ Senior with O&P	25001643	\$71.50
+ HW	ZRL™ Senior with O&P, HW	25001644	\$78.54
Basic CBC, ClinChem 25 with SDMA, T4	ZRL™ Senior	SBASIC	\$59.46
+ FT4	ZRL™ Senior with FT4	25001645	\$74.31

Clinically Focused Bundles

Our Clinical bundles provide you with the information you need to evaluate organ systems or disease states in your patients.



Upgrade any bundle to a Comprehensive CBC. Use test code **CBCCUP**.

Upgrade any bundle with a Urinalysis to a Urinalysis with Reflex UPC. Use test code **25001702**.

Select each code for full test details.

Components	Test Code	CBC	ClinChem	UA	T4	SDMA	Reflex UPC	Fructosamine	Price
ZRL™ Mature Cat	25000500	Basic	Mature Cat	✓	✓	✓	✓	✓	\$86.16
ZRL™ Critical Care	25000508	Basic	Critical Care	✓					\$46.01
ZRL™ Renal	25000478	Basic	Renal	✓		✓	✓		\$60.37
ZRL™ Metabolic	25000493	Basic	Metabolic	✓	✓	✓	✓	✓	\$77.24
ZRL™ Pre-Anesthetic	25000496	Basic	Pre-Anesthetic						\$25.72
ZRL™ Feline Hyperthyroid	25000490	Basic	Feline Hyperthyroid	✓	✓	✓	✓		\$59.93

GI-Focused Bundles

Components	Test Code	ClinChem	ZRL Premium Lipase (DGGR Lipase)	Cbl	Fol	TLI	Cortisol	Giardia Ag	O&P	Price
ZRL™ Canine GI	25000491	25	✓	✓	✓	✓	✓			\$129.61
ZRL™ Canine GI with O&P, Giardia	25001632	25	✓	✓	✓	✓	✓	✓	✓	\$152.16
ZRL™ Feline GI	25000492	25	✓	✓	✓					\$100.47
ZRL™ Feline GI with O&P, Giardia	25001633	25	✓	✓	✓			✓	✓	\$101.02

See the [ClinChem Panels overview](#) for a list of the components in each panel.

UA = Urinalysis
 Cbl = Cobalamin
 Fol = Folate
 TLI = Trypsin-Like Immunoreactivity

ClinChem Panels

Test Name	Test Code	Components	Recommended Use Case	Price
ClinChem 18 with SDMA	25001624	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, TBili, Chol, SDMA	Health screen or undiagnosed sick patient	\$30.18
ClinChem 25 with SDMA	25001625	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, TBili, Chol, CO2, Anion Gap, Mg, AST, GGT, CK, Osmol(c), SDMA	Comprehensive health screen or undiagnosed sick patient	\$33.11
Feline Hyperthyroid ClinChem	25000507	BUN, Creat, K, Phos, Alb, ALT, SDMA, TT4	Monitoring feline hyperthyroid patients	\$23.26
Metabolic ClinChem	25000512	Glu, BUN, Creat, BUN:Creat Ratio, TP, ALP, ALT, Chol, Na, K, Na:K Ratio, Trig, Fructosamine, TT4	Monitoring diabetic patients	\$34.39
Mini Renal ClinChem	25000513	BUN, Creat, K, Phos, Alb	Monitoring renal disease	\$9.11
Renal ClinChem	25000479	Glu, Na, K, Cl, BUN, Creat, BUN:Creat Ratio, Ca, P, CO2, AG, TP, ALB, Chol, Osmol(c), SDMA	Expanded renal disease monitoring	\$22.53
Mature Cat ClinChem	25000482	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, AST, GGT, CK, TBili, Chol, TT4, Fructosamine	Monitoring senior cat health	\$33.47
Critical Care ClinChem, Canine/Feline	25000488	Glu, BUN, Creat, BUN:Creat Ratio, TP, Alb, Na, K, Na:K Ratio, Cl, CO2, Anion Gap, TBili, Osmol(c)	Monitoring electrolytes, acid base, and renal status	\$10.88
Liver ClinChem	25000510	Glu, BUN, TP, Alb, Glob, Alb:Glob Ratio, AST, ALT, ALP, GGT, TBili, Chol	Monitoring liver disease	\$10.88
NSAID ClinChem	25000514	BUN, Creat, BUN:Creat Ratio, ALT	Monitoring effects of NSAID therapy	\$9.11
Phenobarbital ClinChem	25000515	Glu, Na, K, Creat, Ca, TP, Alb, ALT, ALP, TBili, Phenobarbital	Monitoring phenobarbital therapy	\$36.35
Pre-Anesthetic ClinChem	25000516	Glu, Na, K, Creat, Ca, TP, ALT, TBili	Health screen prior to anesthesia/sedation	\$9.11

Bundles

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ZRL™ Junior Complete	JWP	●●●○ 1 mL serum or 2 mL EDTA whole blood, air-dried slides; 5 mL urine in WTT/sterile container	Basic CBC, ClinChem 18 with SDMA, Urinalysis	8 AM next day	\$51.01
ZRL™ Junior Complete with HW	25001642	●●●○ 2 mL serum or 2 mL EDTA whole blood (LTT), air-dried slides; 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 18 with SDMA, Urinalysis, Heartworm Antigen	8 AM next day	\$57.58
ZRL™ Junior with O&P	25001635	●●●● 1 mL serum or 2 mL EDTA whole blood (LTT), air-dried slides; 3g fresh feces	Basic CBC, ClinChem 18 with SDMA, Ova and Parasite	8 AM next day Fecal 5 PM next day	\$71.97
ZRL™ Junior with O&P, HW, Tick	25001636	●●●● 2 mL serum or 2 mL EDTA whole blood (LTT), air-dried slides; 3g fresh feces	Basic CBC, ClinChem 18 with SDMA, Ova and Parasite, Heartworm Antigen, Tick Panel	8 AM next day Fecal 5 PM next day	\$80.03
ZRL™ Junior with O&P, Giardia, HW, Tick	25001637	●●●● 2 mL serum or 2 mL EDTA whole blood (LTT), air-dried slides; 3g fresh feces	Basic CBC ClinChem 18 with SDMA, Ova and Parasite + and <i>Giardia</i> Ag ELISA, Heartworm Antigen, Tick Panel	8 AM next day Fecal 5 PM next day	\$82.66
ZRL™ Junior	25001634	●●● 1 mL serum or 2 mL EDTA whole blood (LTT), air-dried slide	Basic CBC, ClinChem 18 with SDMA	8 AM next day	\$48.65
ZRL™ Junior with T4	25001641	●●● 1 mL serum or 2 mL EDTA whole blood (LTT), air-dried slide	Basic CBC, ClinChem 18 with SDMA, T4	8 AM next day	\$64.85
ZRL™ Junior with HW	25001638	●●● 2 mL serum or 2 mL EDTA whole blood (LTT), air-dried slide	Basic CBC, ClinChem 18 with SDMA, Heartworm Antigen	8 AM next day	\$54.86
ZRL™ Junior with HW, Tick	25001639	●●● 2 mL serum or 2 mL EDTA whole blood (LTT), air-dried slide	Basic CBC, ClinChem 18 with SDMA, Heartworm Antigen, Tick Panel	8 AM next day	\$66.85

Bundles (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ZRL™ Junior with PT/aPTT	25001640	●●●○ 1 mL serum or 2 mL EDTA whole blood (LTT), air-dried slide, 0.5 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT)	Basic CBC, ClinChem 18 with SDMA, PT/aPTT	8 AM next day	\$61.86
ZRL™ Adult Complete	AWP	●●●○● 1 mL serum or 2 mL EDTA whole blood, air-dried slides; 5 mL urine in WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, Urinalysis	8 AM next day	\$50.74
ZRL™ Adult Complete with O&P, HW, Tick	25001631	●●●○●● 2 mL serum or 2 mL EDTA whole blood, air-dried slides; 5 mL urine in WTT/sterile container, 3-5g fresh feces	Basic CBC, ClinChem 25 with SDMA, Urinalysis, Ova and Parasite, Heartworm Antigen, Tick Panel	8 AM next day Fecal 5 PM next day	\$75.56
ZRL™ Adult with O&P	25001630	●●●○● 1 mL serum or 2 mL EDTA whole blood, air-dried slides; 3-5g fresh feces	Basic CBC, ClinChem 25 with SDMA, Ova and Parasite	8 AM next day Fecal 5 PM next day	\$66.63
ZRL™ Adult	25000501	●●● 1 mL serum or 2 mL EDTA whole blood (LTT), air-dried slides	Basic CBC, ClinChem 25 with SDMA	8 AM next day	\$43.64
ZRL™ Adult with HW	25000502	●●● 2 mL serum or 2 mL EDTA whole blood (LTT), air-dried slides	Basic CBC, ClinChem 25 with SDMA, Heartworm Antigen	8 AM next day	\$50.21
ZRL™ Senior Complete	SWP	●●●○ 1 mL serum or 2 mL EDTA whole blood, air-dried slides, 5 mL urine in WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis	8 AM next day	\$74.89

Bundles (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ZRL™ Senior Complete with FT4	25001656	●●●○ 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 5 mL urine in WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, FT4	8 AM next day	\$87.24
ZRL™ Senior Complete with HW, TG	25001657	●●●○ 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 5 mL urine in WTT/sterile container, 0.5 mL of serum or EDTA plasma or whole blood	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Heartworm Antigen, Triglycerides	8 AM next day	\$86.13
ZRL™ Senior Complete with Zoetis Premium Lipase (DGGR), TG	25001655	●●●○ 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 5 mL urine in WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Zoetis Premium Lipase (DGGR), Triglycerides	8 AM next day 1-2 days	\$82.15
ZRL™ Senior	SBASIC	●●● 1 mL serum or 2 mL EDTA whole blood, air-dried slides	Basic CBC, ClinChem 25 with SDMA, T4	8 AM next day	\$59.46
ZRL™ Senior with O&P	25001643	●●●● 1 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces	Basic CBC, ClinChem 25 with SDMA, T4, O&P	8 AM next day Fecal 5 PM next day	\$71.50
ZRL™ Senior with O&P, HW	25001644	●●●● 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces	Basic CBC, ClinChem25 with SDMA, T4, Ova and Parasite, Heartworm Antigen	8 AM next day Fecal 5 PM next day	\$78.54
ZRL™ Senior with FT4	25001645	●●● 2 mL serum or 2 mL EDTA whole blood, air-dried slides	Basic CBC, ClinChem25 with SDMA, T4, FT4	8 AM next day	\$74.31
ZRL™ Senior Complete O&P with FeLV/FIV	25001646	●●●● 1 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, FeLV/FIV	8 AM next day Fecal 5 PM next day	\$89.86


Bundles (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ZRL™ Senior Complete O&P with FeLV/FIV, HW	25001647	●●●○● 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and parasite, FeLV/FIV, Heartworm Antigen	8 AM next day Fecal 5 PM next day	\$103.31
ZRL™ Senior Complete O&P with HW	25001648	●●●○● 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, Heartworm Antigen	8 AM next day Fecal 5 PM next day	\$113.32
ZRL™ Senior Complete O&P, HW, Tick, UPC	25001654	●●●○● 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, Heartworm Antigen, Tick Panel, Urine Protein:Creatinine Ratio	8 AM next day Fecal 5 PM next day	\$113.96
ZRL™ Senior Complete O&P with fPLI	25001649	●●●○● 1 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, Feline Pancreatic Lipase (PLI) (Texas A&M GI Lab)	8 AM next day Fecal 5 PM next day fPLI 2-5 days	\$140.93
ZRL™ Senior Complete O&P with Giardia	25001650	●●●○● 1 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, and <i>Giardia</i> Ag ELISA	8 AM next day Fecal 5 PM next day	\$86.13
ZRL™ Senior Complete O&P with Giardia, HW, Tick	25001651	●●●○● 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, and <i>Giardia</i> Ag ELISA, Heartworm Antigen, Tick Panel	8 AM next day Fecal 5 PM next day	\$103.17
ZRL™ Senior Complete O&P with HW, Tick	25001652	●●●○● 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, Heartworm Antigen, Tick Panel	8 AM next day Fecal 5 PM next day	\$101.25




Bundles (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ZRL™ Senior Complete O&P with HW, Tick, Zoetis Premium Lipase (DGGR), TG	25001653	●●●○ 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 3-5g fresh feces, 5 mL urine in a WTT/sterile container	Basic CBC, ClinChem 25 with SDMA, T4, Urinalysis, Ova and Parasite, Heartworm Antigen, Tick Panel, Zoetis Premium Lipase (DGGR), Triglycerides	8 AM next day Fecal 5 PM next day	\$111.23
ZRL™ Critical Care	25000508	●●●● 2 mL EDTA whole blood, air-dried slides, 1 mL serum, 5 mL urine in WTT/sterile container	Basic CBC, Critical Care ClinChem, Urinalysis	8 AM next day	\$46.01
ZRL™ Feline Hyperthyroid	25000490	●●●● 2 mL EDTA whole blood, air-dried slides, 2 mL serum, 5 mL urine in WTT/sterile container	Basic CBC, Feline Hyperthyroid ClinChem, Urinalysis, Reflex UPC	8 AM next day	\$59.93
ZRL™ Mature Cat	25000500	●●●● 2 mL EDTA whole blood, air-dried slides, 2 mL serum, 5 mL urine in WTT/sterile container	Basic CBC, Mature Cat ClinChem, SDMA, UA, Reflex UPC	8 AM next day	\$86.16
ZRL™ Canine GI	25000491	●● 2 mL serum in RTT	ClinChem 25, ZRL Premium Lipase (DGGR), Cobalamin, Folate, Cortisol, TLI	1-2 days	\$129.61
ZRL™ Canine GI with O&P, Giardia	25001632	●● ZRL™ Canine GI: 2mL serum in RTT 5g fresh feces	ClinChem 25, ZRL Premium Lipase (DGGR), Cobalamin, Folate, Cortisol, TLI, O&P, and <i>Giardia</i> Ag ELISA	ZRL™ Canine GI: 1-2 days O&P & <i>Giardia</i> : 5 PM next day	\$152.16
ZRL™ Complete Coccidioides	25001711	●●●○ 2 mL serum or 2 mL EDTA whole blood, air-dried slides, 2 mL EDTA whole blood, air-dried slides, 5 mL urine in a WTT/sterile container	ClinChem 25 w/ SDMA, Basic CBC, UA, Coccidioides Ab Screening Test, Ehrlichia canis IFA, Triglycerides, ZRL Premium Lipase (DGGR)	2-3 days	\$130.03
ZRL™ Feline GI	25000492	●● 2 mL serum in SST or RTT	ClinChem 25, ZRL Premium Lipase (DGGR), Cobalamin, Folate	1-2 days	\$100.47

Bundles (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ZRL™ Feline GI with O&P, Giardia	25001633	●●● ZRL™ Feline GI: 2mL serum in SST or RTT 5g fresh feces	ClinChem 25, ZRL Premium Lipase (DGGR), Cobalamin, Folate, O&P, and <i>Giardia</i> Ag ELISA	ZRL™ Feline GI: 1-2 days O&P & <i>Giardia</i> : 5 PM next day	\$101.02
ZRL™ Metabolic	25000493	●●●● 2 mL EDTA whole blood, air-dried slides, 1 mL serum, 5 mL urine in WTT/sterile container	Basic CBC, Metabolic ClinChem, Urinalysis, Reflex UPC	1-3 days	\$77.24
ZRL™ Pre-Anesthetic	25000496	●●● 1 mL serum or 2 mL EDTA whole blood, air-dried slides	Basic CBC, Pre-Anesthetic ClinChem	8 AM next day	\$25.72
ZRL™ Renal	25000478	●●●● 2 mL EDTA whole blood, air-dried slides, 1 mL serum, 5 mL urine in WTT/sterile container	Basic CBC, Renal ClinChem with SDMA, Urinalysis, Reflex UPC	8 AM next day	\$60.37
Urinalysis with Urine Culture & Sensitivity	25001629	○● Urinalysis: 5 mL urine in WTT/sterile container Urine C&S: 1 mL urine in WTT or sterile container	Urinalysis, Urine Culture & Sensitivity	Urinalysis: 8 AM next day Urine C&S: 1-5 days	\$82.17
Canine Diarrhea PCR[†] with O&P, Giardia Ag	25001623	O&P: 5g fresh feces <i>Giardia</i> Ag: 3g fresh feces Diarrhea Panel: feces, fecal swab or GI tissue. Minimum 1 g feces in a clean leakproof container. If submitting fecal swab, please submit two. 	Canine Diarrhea PCR Panel, Ova and Parasite, and <i>Giardia</i> antigen ELISA. Canine Diarrhea PCR Panel: <i>Campylobacter jejuni/coli</i> , <i>Salmonella</i> spp., <i>Lawsonia intracellularis</i> , <i>Clostridium difficile</i> toxins A & B, <i>Clostridium perfringens</i> enterotoxin, <i>Giardia intestinalis</i> , Rotavirus, <i>Cryptosporidium</i> spp., Canine Parvovirus, and Coronavirus	O&P & <i>Giardia</i> Ag: Next day by 5 PM Diarrhea Panel: 5-7 days	\$107.91






Bundles (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Feline Diarrhea PCR[†] with O&P, Giardia Ag	25001626	O&P: 5g fresh feces <i>Giardia</i> Ag: 3g fresh feces Diarrhea Panel: feces, fecal swab or GI tissue. Minimum 1 g feces in a clean leakproof container. If submitting fecal swab, please submit two. Diarrhea Panel: RTT, feces, culture swab feces, fecal swab or GI tissue. Minimum 1 g feces in a clean leakproof container. If submitting fecal swab, please submit two. 	Feline PCR Diarrhea Panel, Ova and Parasite, and <i>Giardia</i> antigen ELISA Feline Diarrhea PCR Panel: <i>Campylobacter jejuni/coli</i> , <i>Salmonella</i> spp., <i>Clostridial difficile</i> toxins A & B, <i>Clostridium perfringens</i> enterotoxin, Feline parvovirus (Panleukopenia), <i>Trichomonas foetus</i> , Coronavirus	O&P & <i>Giardia</i> Ag: Next day by 5 PM Diarrhea Panel: 5-7 days	\$117.76
O&P with HW & Tick	25001627	O&P: 5g fresh feces HW: 0.5 mL of serum or EDTA plasma Tick Panel: 0.25 mL of serum or EDTA plasma or whole blood 	Tick Panel: Lateral flow assay for <i>Ehrlichia</i> , <i>Anaplasma</i> , <i>Borrelia</i> (Lyme) Antibodies HW: DiroCHEK® Ag ELISA, Ova, Parasite	O&P: 5 PM next day HW & Tick: 8 AM next day	\$36.84
Tick with HW	25001628	●●● SST, RTT, LTT 0.5 mL of serum or EDTA plasma 	Tick Panel: Lateral flow assay for <i>Ehrlichia</i> , <i>Anaplasma</i> , <i>Borrelia</i> (Lyme) Antibodies HW: DiroCHEK® Ag ELISA	8 AM next day	\$21.59









See the complete Directory of Services online at zoetisreflabs.com/DOS

Hematology/Coagulation




Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Basic CBC, Canine/Feline (Healthy, Non-geriatric)	2000	<ul style="list-style-type: none"> LTT 2 mL EDTA whole blood, air-dried slides  	RBC, HCT, HGB, MCV, MCHC, CHCM, RDW, WBC, WBC Differential, Platelet Count, MPV. If results indicate, will include a PCV, Manual Differential, and/or Reticulocyte Count.	8 AM next day	\$18.19
Comprehensive CBC, Canine/Feline	CBCC	<ul style="list-style-type: none"> LTT 2 mL EDTA whole blood, air-dried slides  	CBC with blood smear review by technician. If results indicate, will include a PCV, Manual Differential, Saline Agglutination Test, Reticulocyte Count, and/or Pathologist Review.	8 AM next day*	\$29.96
Pathologist Review	PATH		Blood smear review by pathologist when ordered pre-emptively. Only comprehensive CBCs are eligible for pathologist review.	1-2 days	\$80.25
Comprehensive CBC Upgrade (add to Bundle)	CBCCUP	<ul style="list-style-type: none"> LTT 2 mL EDTA whole blood, air-dried slides  	CBC with blood smear review by technician. If results indicate, will include a PCV, Manual Differential, Saline Agglutination Test, Reticulocyte Count, and/or Pathologist Review.	8 AM next day*	\$11.77
Reticulocyte	RETIC	<ul style="list-style-type: none"> LTT 1 mL EDTA whole blood  	RBC, % Reticulocyte, Reticulocyte	8 AM next day	\$18.54
Heartworm Antigen Test	6001	<ul style="list-style-type: none"> SST, RTT, LTT 0.25 mL of serum or EDTA plasma  	DiroCHEK® Ag ELISA	8 AM next day	\$10.52
Heartworm Antigen with Heat Treatment	HWHEAT	<ul style="list-style-type: none"> SST, RTT, LTT 0.5 mL plasma or serum (or 1 mL whole blood) 	DiroCHEK® Ag ELISA	1-2 days	\$25.48

* Pathologist Review may extend turnaround time by 1-2 days for final report.

Hematology/Coagulation (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Heartworm Antigen Test (Alternate Methodology)	25000533	<p>●●● SST, RTT, LTT 0.5 mL serum or EDTA plasma or whole blood</p> 	WITNESS® Heartworm Ag lateral flow assay	8 AM next day	\$12.71
Feline Heartworm Antibody (ELISA)†	6007	<p>●● SST, RTT 1 mL serum</p> 		5-7 days	\$142.70
Antithrombin (AT3)†	ANTTHRO3	<p>●○ BTT, WTT 1 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT)</p> 		3-5 days	\$114.64
Blood Typing	BLTYP	<p>● LTT 0.5 mL EDTA whole blood</p> 		8 AM next day	\$114.61
Crossmatch, One Donor	XMTCH1	<p>●● RTT, LTT 1 mL serum AND 1 mL EDTA whole blood required for recipient and each donor</p> 	Blood pigtailed also accepted as donor specimens.	8 AM next day	\$63.65
Crossmatch, Two Donors	XMTCH2	<p>●● RTT, LTT 1 mL serum AND 1 mL EDTA whole blood required for recipient and each donor</p> 	Blood pigtailed also accepted as donor specimens.	8 AM next day	\$114.61



Hematology/Coagulation (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Crossmatch, Three Donors	XMTCH3	<ul style="list-style-type: none"> ● ● RTT, LTT 1 mL serum AND 1 mL EDTA whole blood required for recipient and each donor  	Blood pigtails also accepted as donor specimens.	8 AM next day	\$152.85
Crossmatch, Additional Donor	XMTCHX	<ul style="list-style-type: none"> ● ● RTT, LTT 1 mL serum AND 1 mL EDTA whole blood required for recipient and each donor  	Blood pigtails also accepted as donor specimens.	8 AM next day	\$38.18
DIC Panel†	6565	<ul style="list-style-type: none"> ● ○ BTT, WTT 2-3 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT)  	aPTT, PT, Fibrinogen, Antithrombin, D-Dimer	2-5 days	\$185.25



See the complete Directory of Services online at zoetisreflabs.com/DOS






Hematology/Coagulation (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Flow Cytometry†	FLOW	<ul style="list-style-type: none"> ● LTT 1 mL sample in EDTA. May be performed on blood, bone marrow, cavity fluid, or tissue aspirate. Peripheral blood Samples must have had CBC within previous 48 hours; please submit results with sample. If a CBC is not submitted with a peripheral blood sample, a Comprehensive CBC will be added and performed by ZRL for an additional charge due to the time sensitive nature of this test. 		3-7 days	\$230.05
Flow Cytometry with CBC†	FLOWCBC	<ul style="list-style-type: none"> ● LTT May be performed on peripheral blood only. ZRL will perform a Comprehensive CBC using additional blood included in submission. For Flow Cytometry: 1 mL EDTA whole blood (LTT). For CBC: 2 mL EDTA whole blood (LTT) 		3-7 days	\$279.27

Hematology/Coagulation (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Hemophilia Panel, Canine/Feline†	HEMO	<ul style="list-style-type: none"> ● ○ BTT, WTT 2-3 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 	aPTT, PT, TCT, Factor VIII:C, Factor IX:C, vWF:Ag	3-5 days	\$154.23
Iron, Serum/Plasma, Canine†	IRON	<ul style="list-style-type: none"> ● ● RTT, GRNTT 0.5 mL separated serum or heparinized plasma 📄 		3-5 days	\$84.14
Iron Profile†	IRONPRO	<ul style="list-style-type: none"> ● ● RTT, GRNTT 0.5 mL serum or heparinized plasma 📄 	Iron, TIBC, and Ferritin	7-10 days	\$142.99
Platelet Surface Associated Antibody†	PBA	<ul style="list-style-type: none"> ● LTT 1 mL EDTA whole blood 📄 		4-7 days	\$147.40
Protein C†	PROC	<ul style="list-style-type: none"> ● ○ BTT, WTT 0.5 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 		2-3 days	\$118.92
PT & aPTT	4101	<ul style="list-style-type: none"> ● ○ BTT, WTT 0.5 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 		8 AM next day	\$25.48

Hematology/Coagulation (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
von Willebrand factor (vWF:Ag) [†]	4106	<ul style="list-style-type: none"> ● ○ BTT, WTT 0.5 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT). Freeze. Ship in insulated container with ice packs for next-day delivery 		3-5 days	\$110.89
Buffy Coat	6295	<ul style="list-style-type: none"> ● LTT 1 mL EDTA whole blood 		1-2 days	\$108.25
D-Dimer	4104	<ul style="list-style-type: none"> ● ○ BTT, WTT 0.5 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) 	Semiquantitative D-Dimer assay	1-2 days	\$40.76
Direct Coombs	2701	<ul style="list-style-type: none"> ● LTT 2 mL EDTA whole blood 		1-2 days	\$63.65
Fibrinogen, Canine/Feline	4103	<ul style="list-style-type: none"> ● BTT Citrated whole blood filled according to tube manufacturer 		2-3 days (3-4 days for feline samples)	\$25.48




Hematology/Coagulation (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Hemostasis Panel, Canine/Feline†	6560	<ul style="list-style-type: none"> ● ○ BTT, WTT 2-3 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 	aPTT, PT, TCT, vWF:Ag	3-5 days	\$125.34
Visual Platelet Estimate	2211	<ul style="list-style-type: none"> 1 mL EDTA or citrated whole blood, air-dried slides 📄 		8 AM next day	\$18.55
Technician Blood Smear Review	SMREV	<ul style="list-style-type: none"> ● LTT 2 unstained, air-dried slides; include 2 mL EDTA whole blood, if possible (for in-lab preparation of smears) 📄 Blood 🏠 Slides 	Morphologic review of RBC, WBC, and platelets; Manual WBC Differential, and Visual Platelet Estimate. If indicated based on morphologic findings, a Pathologist Review will be included.	8 AM next day If pathologist review added, TAT is delayed	\$17.80
Coagulation Panel, Canine/Feline (Cornell)†	COAGFA	<ul style="list-style-type: none"> ● ○ BTT, WTT 1 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 	aPTT, PT, TCT	2-3 days	\$104.47
Clotting Factor FIX†	FACTORIX	<ul style="list-style-type: none"> ● ○ BTT, WTT 2 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 		2-3 days	\$91.10

Hematology/Coagulation (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Clotting Factor FVIII†	FACTORVIII	<ul style="list-style-type: none"> ● ○ BTT, WTT 2 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 		2-3 days	\$91.10
Clotting Factor FX†	FACX	<ul style="list-style-type: none"> ● ○ BTT, WTT 2 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 		2-3 days	\$91.10
Erythropoietin (EPO)†	50227	<ul style="list-style-type: none"> ● ○ SST or plasma separator tube 1-2 mL plasma or serum from SST or plasma separator tube (EDTA plasma not acceptable) 		5-7 days	\$185.25
<i>Babesia</i> Blood Smear†	BABBLDSMR	<ul style="list-style-type: none"> ● LTT 2 mL EDTA whole blood 📄 	Giemsa blood smear examination for <i>Babesia</i> organisms (available for export animals only)	3-5 days	\$78.69
Clotting Factor FXII†	FACTORXII	<ul style="list-style-type: none"> ● ○ BTT, WTT 2 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) ❄ 		2-3 days	\$91.10

Hematology/Coagulation (cont'd)









Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Low Molecular Weight Heparin†	LMWHEPARIN	<ul style="list-style-type: none"> ● BTT 0.5 mL citrated plasma (frozen) 	Anti-Xa activity of LMWH	2-3 days	\$118.38
Coagulation Panel and Fibrinogen, Canine/Feline†	25001390	<ul style="list-style-type: none"> ● ○ BTT, WTT 1 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) and frozen 	aPTT, PT, TCT, and Fibrinogen	2-3 days	\$112.50
ZRL™ Coagulation Panel I	COAG	<ul style="list-style-type: none"> ● ● LTT, BTT 1 mL EDTA whole blood and 0.5 mL separated citrated plasma	PT, aPTT, Platelet Count, Total Plasma Protein	8 AM next day	\$95.23
ZRL™ Coagulation Panel II	25001391	<ul style="list-style-type: none"> ● ○ BTT, WTT 0.5 mL separated citrated plasma. The plasma should be aliquoted into a plain tube with no additive (WTT) 	PT, aPTT, Fibrinogen, D-Dimer	2-3 days	\$91.71



See the complete Directory of Services online at zoetisreflabs.com/DOS







Chemistry

Chemistry Panels

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Feline Hyperthyroid ClinChem	25000507	<ul style="list-style-type: none"> ● SST or RTT 2 mL serum  	BUN, Creat, K, Phos, Alb, ALT, SDMA, TT4	8 AM next day	\$23.26
Metabolic ClinChem	25000512	<ul style="list-style-type: none"> ● SST or RTT 2 mL serum  	Glu, BUN, Creat, BUN:Creat Ratio, TP, ALP, ALT, Chol, Na, K, Na:K Ratio, Trig, Fructosamine, TT4	1-3 days	\$34.39
Renal ClinChem	25000479	<ul style="list-style-type: none"> ● SST or RTT 2 mL serum  	Glu, Na, K, Cl, BUN, Creat, BUN:Creat Ratio, Ca, P, FeLV CO2, AG, TP, ALB, Chol, Osmol(c), SDMA	8 AM next day	\$22.53
Mature Cat ClinChem	25000482	<ul style="list-style-type: none"> ● SST or RTT 1 mL serum  	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, AST, GGT, CK, TBili, Chol, TT4, Fructosamine	1-2 days	\$33.47
Critical Care ClinChem, Canine/Feline	25000488	<ul style="list-style-type: none"> ● SST or RTT 1 mL serum  	Glu, BUN, Creat, BUN:Creat Ratio, TP, Alb, Na, K, Na:K Ratio, Cl, CO2, Anion Gap, TBili, Osmol(c)	8 AM next day	\$10.88
Liver ClinChem	25000510	<ul style="list-style-type: none"> ● SST or RTT 1 mL serum  	Glu, BUN, TP, Alb, Glob, Alb:Glob Ratio, AST, ALT, ALP, GGT, TBili, Chol	8 AM next day	\$10.88
Mini Renal ClinChem	25000513	<ul style="list-style-type: none"> ● SST or RTT 1 mL serum  	BUN, Creat, K, Phos, Alb	8 AM next day	\$9.11
NSAID ClinChem	25000514	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  	BUN, Creat, BUN:Creat Ratio, ALT	8 AM next day	\$9.11








Chemistry (cont'd)

Chemistry Panels (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Phenobarbital ClinChem	25000515	<ul style="list-style-type: none"> ● RTT 2 mL serum, separated 	Glu, Na, K, Creat, Ca, TP, Alb, ALT, ALP, TBili, Phenobarbital RTT serum preferred for phenobarbital testing. SST serum may affect phenobarbital results.	1-2 days	\$36.35
Pre-Anesthetic ClinChem	25000516	<ul style="list-style-type: none"> ● ● SST or RTT 1 mL serum 	Glu, Na, K, Creat, Ca, TP, ALT, TBili	8 AM next day	\$9.11
ClinChem 18	CHEM18	<ul style="list-style-type: none"> ● ● SST or RTT 0.5 mL serum 	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, TBili, Chol	8 AM next day	\$20.33
ClinChem 18 with SDMA	25001624	<ul style="list-style-type: none"> ● ● SST or RTT 0.5 mL serum 	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, TBili, Chol, SDMA	8 AM next day	\$30.18
ClinChem 25	1000	<ul style="list-style-type: none"> ● ● SST or RTT 1 mL serum 	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, TBili, Chol, CO ₂ , Anion Gap, Mg, AST, GGT, CK, Osmol(c)	8 AM next day	\$26.75
ClinChem 25 with SDMA	25001625	<ul style="list-style-type: none"> ● ● SST or RTT 1 mL serum 	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, TP, Alb, Glob, Alb:Glob Ratio, ALT, ALP, TBili, Chol, CO ₂ , Anion Gap, Mg, AST, GGT, CK, Osmol(c), SDMA	8 AM next day	\$33.11








Chemistry (cont'd)

Individual Chemistry Tests

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Triglycerides	1023	●● SST or RTT 0.5 mL serum 		8 AM next day	\$8.93
SDMA	7037	●● SST or RTT 0.5 mL serum 		8 AM next day	\$18.00
Bile Acids Pre/Resting	6180	●● SST or RTT 0.5 mL serum 		1-2 days	\$27.99
Bile Acids Post	6181	●● SST or RTT 0.5 mL serum 		1-2 days	\$27.99
Free (Ionized) Calcium	IONISTAT	● GRNTT 0.5 mL lithium heparinized whole blood 	Ensure that the sample is treated anaerobically and that the correct ratio of heparin anticoagulant to blood is achieved (do not underfill GRNTT).	1-2 days	\$50.95
Free (Ionized) Calcium (Mich. State Univ.)†	6304	● RTT 0.75 mL serum, separated 		3-6 days	\$79.86
LDH†	LDH	●● RTT or GRNTT 0.5 mL separated serum or heparinized plasma (effusion fluid and urine also accepted) 		3-5 days	\$87.35









Chemistry (cont'd)

Individual Chemistry Tests (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Osmolality†	OSMOLAL	●● RTT or GRNTT 1 mL separated heparinized or EDTA plasma 		2-3 days	\$95.91
Albumin	1001	●● SST or RTT 0.5 mL serum 		8 AM next day	\$7.10
Canine C-reactive Protein (CRP)†	CRP	●●● RRT, GRNTT or LTT 0.5 mL serum or plasma (heparinized or EDTA) 		5-7 days	\$142.70
High-Sensitivity Cardiac Troponin I†	TROPONIN	●● SST or RTT 0.5 mL non-hemolyzed serum (frozen) 		4-5 days	\$136.57
Corticosteroid-Induced Alk Phos (CIALP)†	CIALP	●●● SST, RTT or GRNTT 0.5 mL serum or heparinized plasma 		2-3 days	\$148.72
NT-proBNP†	PROBNP	●●● LTT, SST or RTT 1 mL EDTA plasma (canine, feline) or 1 mL serum (feline only) 		3-7 days	\$170.65
ALT	1002	●● SST, RTT 0.5 mL serum 		8 AM next day	\$8.75









Chemistry (cont'd)

Individual Chemistry Tests (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ALP (Alk Phos)	1003	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
AST	1005	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.49
BUN	1006	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
Chloride	1007	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.49
Calcium (Total Calcium)	1008	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
Cholesterol	1009	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
Creatinine	1010	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
Creatinine Kinase (CK)	1011	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.10



Chemistry (cont'd)

Individual Chemistry Tests (cont'd)



Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
GGT	1013	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.49
Glucose	1014	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
Magnesium	1016	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.49
Phosphorus	1017	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
Potassium	1018	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$8.75
Sodium	1019	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.49
CO2	1020	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.49
Total Bilirubin	1021	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum  		8 AM next day	\$7.49

Chemistry (cont'd)





Individual Chemistry Tests (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Total Protein	1022	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum 		8 AM next day	\$7.49
Glomerular Filtration Rate (Iohexol Clearance)†	GFRIC	<ul style="list-style-type: none"> ● SST or RTT 0.5 mL serum (3 samples required) 	Serum samples should be taken two, three, and four hours following a single injection of iohexol (300 mg iodine equivalents/kg body weight). Each sample should yield at least 0.5 mL serum.	6-8 days	\$195.42

Urine

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Urinalysis	3000	<ul style="list-style-type: none"> ● Urine 5 mL urine in YTT/sterile container 		8 AM next day	\$16.95
Urinalysis with Urine Culture & Sensitivity	25001629	Urinalysis: 5 mL urine in YTT/sterile container Urine C&S: 1 mL urine in YTT or sterile container		Urinalysis: 8 AM next day Urine C&S: 1-5 days	\$82.17
Urine Protein:Creatinine Ratio	1555	<ul style="list-style-type: none"> ● Urine 5 mL urine in YTT/sterile container 	Urine Creatinine, Urine Protein, Urine Protein:Creatinine Ratio	8 AM next day	\$16.95








Urine (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Urinalysis with Reflex UPC	25000521	<ul style="list-style-type: none"> ● Urine 5 mL urine in YTT/sterile container 	An automatic Urine Protein:Creatinine Ratio will be performed when proteinuria is present.	8 AM next day	\$23.25
Urinalysis with Reflex UPC upgrade	25001702	<ul style="list-style-type: none"> ○ ● WTT, Urine 5 mL urine in sterile container 	Add to ZRL Bundles including Urinalysis		\$6.30
Urine Electrolyte Panel (for Fractional Excretion) [†]	FREXCRE	<ul style="list-style-type: none"> ● Urine Sterile Leakproof Container YTT/Sterile Leakproof Container 2 mL urine in a YTT/Sterile Leakproof Container AND 2 mL separated serum OR heparinized plasma 	Cl, Creat, K, Na Concurrent measurement of serum creatinine and electrolytes is needed for determination of fractional excretion.	3-5 days	\$118.38
Stone Analysis [†]	7032	Dry stones in sterile container with no additive 		4-9 days	\$109.82



See the complete Directory of Services online at zoetisreflabs.com/DOS








Endocrinology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Total T4	5102	●● SST or RTT 0.5 mL serum 		8 AM next day	\$16.95
Free T4	5103	●● SST or RTT 0.5 mL serum 		8 AM next day	\$19.80
TSH	5106	●● SST or RTT 0.5 mL serum 		8 AM next day (1-2 days for feline samples)	\$16.08
Thyroid Panel	25000484	●● SST or RTT 1 mL serum 	TT4, FT4, TSH	8 AM next day (1-2 days for feline samples)	\$29.43
Free T4 by ED [†]	6330	● RTT 0.5 mL serum, separated and frozen 	Ship frozen on ice.	5-7 days	\$120.52
Thyroglobulin Autoantibody [†]	6437	● RTT 0.5 mL serum, separated		3-5 days	\$94.84
Total T3 (TT3), Canine/Feline [†]	1396	● RTT 1.5 mL serum, separated or heparinized plasma		2-5 days	\$79.28
Premium Thyroid Panel, Canine [†]	6545	● RTT 2 mL serum, separated and frozen 	TT4, TT3, FT4 by ED, TSH, TgAA, T4AA, and T3AA Ship frozen on ice.	5-7 days	\$150.48
Premium Thyroid Panel + TSH, Feline [†]	6550	● RTT 2 mL serum, separated and frozen 	TT4, TT3, FT4 by ED, TSH Ship frozen on ice.	5-7 days	\$147.40







Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
OFA Thyroid Panel, Canine†	K9THYROID	<ul style="list-style-type: none"> ● RTT 2 mL serum, separated ❄️ 	TT4, TT3, FT4 by ED, TSH, TgAA, T4AA, and T3AA Ship frozen on ice. Animal must not have received thyroid therapy for 3 months prior to sampling. A completed OFA Thyroid Registry form must accompany the sample. This may be faxed or emailed to your clinic on request to the laboratory, or may be downloaded from the OFA website at the following link: www.ofa.org . Please note that OFA Registry requires prepayment by the animal owner of \$15.00* for the registry service. This fee is not included in the price shown to the right.	2-3 weeks	\$148.34
Fructosamine	1024	<ul style="list-style-type: none"> ● ● SST or RTT 0.5 mL serum 📄 		1-3 days	\$22.55
Insulin:Glucose Ratio†	MSUIG	<ul style="list-style-type: none"> ● RTT 1 mL serum, separated. Freeze SEPARATED serum prior to shipping. 📄 	Ship on ice. Customer service may assist with frozen package pickup.	3-7 days	\$142.70
Malignancy (Hypercalcemia) Panel†	6555	<ul style="list-style-type: none"> ● ● RRT and LTT 1 mL serum AND 0.5 mL EDTA plasma, both separated and frozen. ❄️ 	PTH, PTHrp, Free (Ionized) Calcium Sample must arrive frozen to the laboratory. Shipping on dry ice is recommended. Customer service may assist with frozen package pickup.	5-7 days	\$156.90
Parathyroid Panel†	6306	<ul style="list-style-type: none"> ● RTT 1 mL serum, separated and frozen ❄️ 	PTH, Free (Ionized) Calcium Sample must arrive frozen to the laboratory. Shipping on dry ice is recommended. Customer service may assist with frozen package pickup.	5-7 days	\$122.16
Parathyroid Hormone-Related Protein (PTHrP)†	PTHRP	<ul style="list-style-type: none"> ● LTT 1 mL EDTA plasma, separated and frozen ❄️ 	Sample must arrive frozen to the laboratory. Shipping on dry ice is recommended. Customer service may assist with frozen package pickup.	5-7 days	\$126.94

Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Vitamin D (25-hydroxyvitamin D3)[†]	VITD	● RTT 0.75 mL serum, separated and frozen 	Sample must arrive frozen to the laboratory. Shipping on dry ice is recommended. Customer service may assist with frozen package pickup.	4-10 days	\$130.54
Vitamin D Profile[†]	VITP	● RTT 1.5 mL serum, separated and frozen 	PTH, Free (Ionized) Calcium, and 25-hydroxyvitamin D Sample must arrive frozen to the laboratory. Shipping on dry ice is recommended. Customer service may assist with frozen package pickup.	3-9 days	\$140.53
Cortisol (Baseline/Resting)	5101	● ● SST or RTT 0.5 mL serum 		8 AM next day (1-2 days for feline/equine samples)	\$16.08
ACTH Stimulation (2 samples)	6100	● ● SST or RTT 0.5 mL serum each sample (2 samples required) 		8 AM next day (1-2 days for feline samples)	\$30.57
ACTH Stimulation (3 samples)	6101	● ● SST or RTT 0.5 mL serum each sample (3 samples required) 		8 AM next day (1-2 days for feline samples)	\$45.85
Cortisol (Single Sample, Post-ACTH)	CORTPOST	● ● SST or RTT 0.5 mL serum 	This test code is intended to support submission of a single, post-ACTH cortisol serum sample when other samples for dynamic testing have been submitted on another occasion.	8 AM next day (1-2 days for feline samples)	\$22.95
Dexamethasone Suppression (2 samples)	6103	● ● SST or RTT 0.5 mL serum each sample (2 samples required) 		8 AM next day (1-2 days for feline samples)	\$30.57







Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Dexamethasone Suppression (3 samples)	6104	<ul style="list-style-type: none"> ● ● SST or RTT 0.5 mL serum each sample (3 samples required)  		8 AM next day (1-2 days for feline samples)	\$50.95
Urine Cortisol:Creatinine Ratio	3310	<ul style="list-style-type: none"> ● Urine 1 mL in YTT/sterile container  		8 AM next day (1-2 days for feline samples)	\$70.06
Adrenal Panel (with ACTH Stimulation), Canine/Feline†	ACTHADR	<ul style="list-style-type: none"> ● RTT 2 mL serum each, separated (2 samples)  	Baseline and post-ACTH Cortisol, Estradiol, Androstenedione, Progesterone, 17-Hydroxyprogesterone, Testosterone Requires two samples. Testing may be performed on canines and neutered felines.	12 days	\$288.90
17-Hydroxyprogesterone, Pre & Post-(ACTH Stim)†	17HPROG	<ul style="list-style-type: none"> ● RTT 0.5 mL serum, separated  	Requires two samples. Testing may be performed on canines and neutered felines.	7-10 days	\$91.63
Aldosterone Baseline†	ALDOS	<ul style="list-style-type: none"> ● ● ● LTT or SST, RTT 0.5 mL EDTA plasma (preferred); 0.5 mL serum also accepted. Do not send whole blood.  	Centrifuge EDTA sample immediately and separate EDTA plasma within 30 minutes of collection. Transfer to sterile, no additive tube and label tube EDTA plasma. Freeze sample and ship on ice.	7-10 days	\$100.19
Aldosterone Post†	ALDOSPOST	<ul style="list-style-type: none"> ● ● ● LTT or SST, RTT 0.5 mL EDTA plasma (preferred); 0.5 mL serum also accepted. Do not send whole blood.  	Centrifuge EDTA sample immediately and separate EDTA plasma within 30 minutes of collection. Transfer to sterile, no additive tube and label tube EDTA plasma. Freeze sample and ship on ice.	7-10 days	\$96.92

Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Endogenous ACTH, Canine†	25001384	<ul style="list-style-type: none"> ● LTT 0.5 mL EDTA plasma (separated and frozen) ❄️ 	<p>Centrifuge immediately and pipette EDTA plasma into a sterile, no-additive, plastic tube. Glass tubes should not be used.</p> <p>Freeze sample. Sample must arrive frozen to the laboratory. Shipping on dry ice is recommended. Customer service may assist with frozen package pickup.</p>	2-5 days	\$105.54
Reproductive Hormone Panel†	RHP	<ul style="list-style-type: none"> ● ● LTT, RTT (no SST) 2 mL serum or plasma (no SST) 	Estrogen, Progesterone, Testosterone	5-10 days	\$153.71
Estradiol, Baseline†	ESTRADIOL	<ul style="list-style-type: none"> ● RTT 0.5 mL serum, separated 🧪 		7-10 days	\$99.12
Testosterone, Baseline†	TESTOS	<ul style="list-style-type: none"> ● RTT 2 mL serum or EDTA 🧪 		2-9 days	\$96.98
Testosterone Post-HCG/GnRH†	TESTPOST	<ul style="list-style-type: none"> ● ● RTT, LTT 1 mL serum or EDTA plasma 🧪 		3-8 days	\$89.49
Progesterone	5105	<ul style="list-style-type: none"> ● RTT 0.5 mL serum 🧪 	Serum separator gel may affect results.	8 AM next day	\$44.60
Anti-Müllerian Hormone†	5555	<ul style="list-style-type: none"> ● RTT 0.5 mL non-hemolyzed serum, separated and frozen ❄️ 	Ship frozen and on ice. Sample should not be stored in a refrigerator for more than 2 days. Wait 30 days after spay or remnant removal surgery before testing.	5-8 days	\$148.34
Ovarian Remnant Syndrome†	OVR5	<ul style="list-style-type: none"> ● RTT 0.5 mL serum, separated and frozen ❄️ 	<p>Anti-Müllerian Hormone, Progesterone</p> <p>Ship frozen on ice. Wait 30 days after spay or remnant removal surgery before testing.</p>	8-10 days	\$160.00












Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Luteinizing Hormone (LH), Canine/Feline	LH	<ul style="list-style-type: none"> ●● SST, RTT 0.5 mL serum  same day 		8 AM next day	\$70.06
Pregnancy Test (Relaxin), Canine/Feline	6009	<ul style="list-style-type: none"> ●●● SST, RTT or LTT 0.5 mL serum or EDTA plasma  	Relaxin degrades rapidly at room temperature. Keep refrigerated until shipment. Ship on ice packs.	1-2 days	\$57.30
Insulin-like Growth Factor†	ILGROFAC	<ul style="list-style-type: none"> ●● RTT, LTT 0.5 mL serum or EDTA plasma  		5-10 days	\$139.78
Gastrin†	GASTRIN	<ul style="list-style-type: none"> ● RTT 0.5 mL serum, separated and frozen  	Fast patient for 12 hours. Freeze sample. Sample must arrive frozen to the laboratory. Shipping on dry ice is recommended. Alternatively, ship on ice packs. Customer service may assist with frozen package pickup.	3-12 days	\$113.03
Metanephrines, Plasma-free†	METAPLAS	<ul style="list-style-type: none"> ●● LTT, GRNTT 2 mL EDTA or heparinized plasma, separated and frozen  	Centrifuge and separate plasma within 1 hr of collection. Freeze immediately. Ship on dry ice. Customer service may assist with frozen package pickup.	5-7 days	\$216.78
Metanephrines, Urine Fractionated†	METAUR	<ul style="list-style-type: none"> ● Urine 2.5-4 mL urine in YTT/sterile container  	Freeze immediately. Ship on dry ice. Customer service may assist with frozen package pickup.	5-7 days	\$115.86











See the complete Directory of Services online at zoetisreflabs.com/DOS



GI/Pancreatic

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ZRL Premium Lipase (DGGR Lipase)	1030	● ● SST, RTT 0.5 mL serum 		8 AM next day	\$12.53
ZRL Premium (DGGR Lipase), TLI, Cobalamin, and Folate, Canine	7250	● ● SST, RTT 2 mL serum  		1-2 days	\$165.55
Pancreatic Lipase Immunoreactivity (PLI), Canine (TX GI)†	PLITX	● ● SST, RTT 0.5 mL separated serum 		3-5 days	\$125.87
Pancreatic Lipase Immunoreactivity (PLI), Feline (TX GI)†	7240	● ● SST, RTT 0.5 mL separated serum 		2-5 days	\$125.87
Amylase (Mich. State Univ.)†	1004	● SST 0.5 mL serum 		2-3 days	\$88.32
PLI, Cobalamin, and Folate, Canine/Feline (TX GI)†	FEPLICOFOL	● ● SST, RTT 1 mL separated non-hemolyzed serum  		2-5 days	\$125.87
PLI, TLI, Cobalamin, and Folate, Canine (TX GI)†	K9PTCFTX	● ● SST, RTT 1 mL separated non-hemolyzed serum  		2-5 days	\$160.00
PLI and TLI, Feline (TX GI)†	7253	● ● SST, RTT 1 mL separated serum 		2-5 days	\$128.45






GI/Pancreatic (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
PLI, TLI, Cobalamin, and Folate, Feline (TX GI)†	7251	● ● SST, RTT 2 mL separated non-hemolyzed serum 		2-5 days	\$160.00
TLI, Canine	6206	● ● SST, RTT 		1-2 days	\$53.52
Trypsin-Like Immunoreactivity (TLI), Feline (TX GI)†	6203	● ● SST, RTT 0.2 mL separated serum 		2-5 days	\$98.05
TLI, Cobalamin, and Folate, Canine	7252	● ● SST, RTT 2 mL serum 		1-2 days	\$108.25
TLI, Cobalamin, and Folate, Feline (TX GI)†	F7251	● ● SST, RTT 1 mL separated non-hemolyzed serum 		2-5 days	\$125.87
Cobalamin (Vitamin B12)	6207	● ● SST, RTT 1 mL separated serum 		1-2 days	\$38.18
Folate	6208	● ● SST, RTT 1 mL serum 		1-2 days	\$38.18
Cobalamin and Folate	KFCOFO	● ● SST, RTT 2 mL serum 		1-2 days	\$63.65










GI/Pancreatic (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Microbiota Dysbiosis Index, Canine†	DYSINDEX	● Feces 1 gm feces 	Quantifies abundance of <i>Faecalibacterium</i> , <i>Turicibacter</i> , <i>Streptococcus</i> , <i>E coli</i> , <i>Blautia</i> , and <i>C hiranonis</i>	2-3 days	\$119.65
Fecal Occult Blood	7255	● Feces 1-3 gm fresh feces 		1-2 days	\$38.18






Fecal/Parasitology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Ova and Parasite (Centrifugal Fecal Flotation) Test	2320A	● Feces 3-5g fresh feces 	Fecal floatation by sugar solution. Recommended for most comprehensive fecal parasite assessment.	Next day by 5 PM	\$16.80
<i>Giardia</i> Focused O&P	25001611	● Feces 3-5g fresh feces 	Fecal floatation by zinc sulfate. Recommended for diagnosis of Giardiasis.	Next day by 5 PM	\$10.70
O&P (Centrifugal Fecal Flotation) Panel	25001613	● Feces 10 g fresh feces 	Two fecal flotations by both sugar and zinc sulfate for most comprehensive assessment of parasitic infection. Quantification not included.	Next day by 5 PM	\$21.41
<i>Giardia</i> Antigen Test	2321	● Feces 3 gm fresh feces 	<i>Giardia</i> Ag ELISA	Next day by 5 PM	\$13.16
Ova and Parasite (Centrifugal Fecal Flotation) Test and <i>Giardia</i> Antigen Test	2320	● Feces 5 gm fresh feces 	Fecal floatation by sugar solution and <i>Giardia</i> Ag ELISA	Next day by 5 PM	\$19.48

Fecal/Parasitology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Fecal Direct Smear Test	25001610	● Feces 1-3 g fresh feces 	Feces in sealed, leakproof container or plastic bag with excess air expressed	1-3 days	\$8.01
<i>Giardia</i> & <i>Cryptosporidium</i> IFA, Canine/Feline	2404	● Feces 3-5 gm fresh feces 	<i>Giardia</i> and <i>Cryptosporidium</i> spp. Ag by IFA. Not recommended for use in birds/reptiles.	1-4 days	\$39.50
O&P (Centrifugal Fecal Flotation) with Direct Smear	25001614	● Feces 5g fresh feces 	Feces in sealed, leakproof container or plastic bag with excess air expressed	1-3 days	\$18.71
Fecal Sedimentation Test	25001612	● Feces 10 g fresh feces 	Feces in sealed, leakproof container or plastic bag with excess air expressed	1-4 days	\$53.47
McMaster (Fecal Egg Counting) Test	2328	● Feces 5 gm fresh feces 	Quantitative fecal egg count	1-2 days	\$25.48
Centrifugal Fecal Flotation Test (O&P) & McMaster (Fecal Egg Counting) Test	25000487	● Feces 8-10 gm fresh feces 	Fecal floatation by sugar solution and quantitative fecal egg count.	1-2 days	\$25.48
Baermann Test (Lungworm)	BAER	● Feces 10 gm fresh feces 		1-4 days	\$53.47
Heartworm Antigen Test	6001	● ● ● SST, RTT, LTT 0.5 mL serum or EDTA plasma 	DiroCHEK® HW Ag ELISA	8 AM next day	\$10.52
Heartworm Antigen with Heat Treatment	HWHEAT	● ● ● SST, RTT, LTT 0.5 mL serum or EDTA plasma 	DiroCHEK® HW Ag ELISA	1-2 days	\$25.48



Fecal/Parasitology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Heartworm Antigen Test (Alternate Methodology)	25000533	<ul style="list-style-type: none"> ● ● ● SST, RTT, LTT 0.5 mL serum or EDTA plasma or whole blood  	Zoetis Witness® Heartworm Ag lateral flow assay. Intended for use as alternative methodology, confirmatory testing after initial HW Ag positive test by different methodology per American Heartworm Society guidelines.	8 AM next day	\$12.71
O&P with HW & Tick	25001627	<ul style="list-style-type: none"> O&P: 3g fresh feces HW: 0.25 mL of serum or EDTA plasma Tick Panel: 0.25 mL of serum or EDTA plasma or whole blood  	<ul style="list-style-type: none"> O&P: Fecal flotation by sugar solution Tick Panel: Lateral flow assay for <i>Ehrlichia</i>, <i>Anaplasma</i>, <i>Borrelia</i> (Lyme) Antibodies HW: DiroCHEK® Ag ELISA, Ova and Parasite 	<ul style="list-style-type: none"> O&P: Next day by 5 PM HW & Tick: 8 AM next day 	\$36.84
Modified Knott's Test for Microfilariae (Heartworm)	6225	<ul style="list-style-type: none"> ● LTT 2 mL EDTA whole blood  	Quantification not included.	1-2 days	\$38.18
Modified Knott's Test with Quantification	MODKNOTTSQ	<ul style="list-style-type: none"> ● LTT 2 mL EDTA whole blood  		1-2 days	\$37.10
Feline Heartworm Antibody (ELISA)†	6007	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum  	HW Ab ELISA	5-7 days	\$142.70
Gross Parasite Identification	25001435	Whole parasite (adult worm, larvae, ectoparasite, etc.) in sealed, leakproof container		1-5 days	\$34.35




See the complete Directory of Services online at zoetisreflabs.com/DOS

Microbiology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Aerobic Culture	7156	○ Culture swab or fluid/ tissue		2-5 days	\$69.55
Anaerobic Culture	7157	○ Culture swab or fluid/ tissue		2-7 days	\$38.18
Aerobic and Anaerobic Culture	7135	○ Culture swab or fluid/ tissue		2-7 days	\$107.00
Urine Culture and Sensitivity	7140	○ ● YTT/sterile tube, Urine 1 mL urine in WTT or sterile container 		1-5 days	\$69.55
Urinalysis with Urine Culture & Sensitivity	25001629	● Urinalysis: 5 mL urine in YTT/sterile container Urine C&S: 1 mL urine in WTT or sterile container 		Urinalysis: 8 AM next day Urine C&S: 1-5 days	\$82.17
Hold Urine Culture for 48 hours	25000534		Bacterial culture will be held for an additional 24 hours (48 hours total). This test is intended for use in working up complicated/chronic urinary tract infections and must be ordered at the time the urine culture is initially requested.	2 days	\$0.01
Fungal Culture	7145	○ Culture swab or fluid/ tissue		Up to 26 days for final report	\$70.06
Fecal Culture	7155	● Feces 5 gm fresh feces		2-7 days	\$70.06
Dermatophyte Culture	7160	○ WTT/sterile tube Plucked hair, nail clippings, skin scrapings or brushings		12 days for final report	\$70.06

Microbiology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Mycoplasma</i> Culture†	MYCOPLASCU	Amies transport media swab, fluid or tissue in WTT/sterile container.		Up to 21 days	\$93.77
Blood Culture	7150	Two blood culture bottles, each inoculated with 3 mLs of blood 	Submit two blood samples, drawn at least one hour apart and from different, disinfected venipuncture sites. Store at room temperature. DO NOT REFRIGERATE.	3-9 days	\$140.07
<i>Brucella canis</i> Culture†	BRUCLCULT	As indicated: Amies transport media swab, whole blood inoculated into blood culture bottle, 5 g lung tissue, 5 g placenta, 5 mL semen, 1 g stomach contents, 1 g testicle; submit in sterile, leakproof container(s).	Blood submitted in any tube other than blood culture bottle will be rejected.	7-10 days	\$105.54
<i>Mycobacterium</i> Culture†	MYCO	Intestine, lymph node, rectal, or skin biopsy in sterile, leakproof container		60 days	\$104.47
Acid-Fast Stain, Non-Avian/Reptile	ACIDFAST	○ Culture swab or fluid/tissue		1-4 days	\$43.28







Microbiology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Fungal Susceptibility 1 Drug†	FUNGALSUS1	ZRL fungal culture	Add on fungal susceptibility testing for one antifungal drug.	1-2 weeks	\$136.87
Fungal Susceptibility 2 Drugs†	FUNGALSUS2	ZRL fungal culture	Add on fungal susceptibility testing for two antifungal drugs.	1-2 weeks	\$216.50
Fungal Susceptibility 3 Drugs†	FUNGALSUS3	ZRL fungal culture	Add on fungal susceptibility testing for three antifungal drugs.	1-2 weeks	\$296.13
Skin Culture and Sensitivity (Dermatology)	25001434	○ Culture swab or fluid/ tissue	Culture and AMS designed for skin lesions, including ear infections.	2-5 days	\$63.66











See the complete Directory of Services
online at zoetisreflabs.com/DOS







Serology/Infectious Diseases

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Tick Panel Screening Test	6002	<p>● ● ● SST, RTT or LTT 0.5 mL of serum or EDTA plasma or whole blood</p> 	Lateral flow assay for <i>Ehrlichia</i> , <i>Anaplasma</i> , <i>Borrelia</i> (Lyme) Antibodies Canines only.	8 AM next day	\$14.96
FeLV/FIV Rapid Test	6005	<p>● ● ● SST, RTT or LTT 0.5 mL of serum or EDTA plasma or whole blood</p> 	Lateral flow assay for FeLV Ag and FIV Ab with confirmation of positive FeLV results by ViraCHEK® ELISA	8 AM next day	\$16.08
O&P with HW & Tick	25001627	<p>O&P: 3g fresh feces HW: 0.25 mL of serum or EDTA plasma Tick Panel: 0.25 mL of serum or EDTA plasma or whole blood</p> 	<p>O&P: Fecal flotation by sugar solution Tick Panel: Lateral flow assay for <i>Ehrlichia</i>, <i>Anaplasma</i>, <i>Borrelia</i> (Lyme) Antibodies HW: DiroCHEK® Ag ELISA, Ova and Parasite</p>	<p>O&P: Next day by 5 PM HW & Tick: 8 AM next day</p>	\$36.84
<i>Cryptococcus</i> Antigen Titer by Latex Agglutination	7216	<p>● ● SST, RTT 0.5 mL serum or CSF</p> 		1-2 days	\$95.54
<i>Coccidioides</i> Antibody Screening Test	25001690	<p>● ● SST or RTT 1 mL serum</p>	Lateral flow assay assessing IgG and IgM Abs combined.	2-3 days	\$48.15
<i>Coccidioides</i> Titer AGID [†]	6129	<p>● ● SST or RTT 1 mL serum</p> 	IgG, IgM	3-5 days	\$107.68
<i>Blastomyces</i> Quantitative Antigen EIA [†]	BLASTANT	<p>● ● ● SST, RTT or LTT 1.2 mL serum or plasma, separated; or 0.8 mL CSF; or 0.5 mL urine or BAL fluid</p> 		2-5 days	\$145.13









Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Fungal Titer Panel AGID[†]	6130	●● SST or RTT 1 mL serum 	<i>Coccidioides</i> (IgG/IgM), <i>Blastomyces</i> (IgG), and <i>Histoplasma</i> (IgG). If indicated, can substitute <i>Aspergillus</i> (IgG) and/or <i>Cryptococcus</i>	3-6 days	\$160.50
Vector Serology Panel, Canine[†]	EXTENTICK	●● SST or RTT 2 mL serum 	<i>Babesia canis</i> , <i>B gibsoni</i> , <i>Ehrlichia canis</i> , <i>Rickettsia rickettsii</i> , <i>Bartonella henselae</i> , <i>B vinsoni</i> , <i>B koehlerae</i> (by IFA); <i>Borrelia</i> (Lyme), <i>Anaplasma</i> , <i>Ehrlichia</i> , and heartworm (by rapid test)	10-14 days	\$187.25
Tick Serology Panel	7185	●● SST, RTT 3 mL serum 	<i>Ehrlichia canis</i> (IFA), <i>Anaplasma phagocytophilum</i> (IFA), <i>Borrelia burgdorferi</i> (Cornell Multiplex Antibody Test) [†]	2-7 days	\$152.85
Toxoplasma gondii ELISA[†]	TOXAB189	●● SST, RTT 0.5 mL serum 	IgG and IgM Ab to <i>Toxoplasma gondii</i> . CSF and aqueous humor also accepted. Do not freeze.	5-7 days	\$93.77
Tick with HW	25001628	●●● SST, RTT, LTT 0.5 mL of serum or EDTA plasma 	Tick Panel: Lateral flow assay for <i>Ehrlichia</i> , <i>Anaplasma</i> , <i>Borrelia</i> (Lyme) Heartworm: DiROCHEK [®] ELISA	8 AM next day	\$21.59
Vaccine Titer (Qualitative), Canine	6145	●● SST, RTT 1 mL serum 	Qualitative CDV and parvovirus titer by ELISA	1-3 days	\$47.62
Feline Vaccine Panel 2[†]	7200	● RTT 1 mL serum 	Herpesvirus, calicivirus, parvovirus (panleukopenia) Ab	6-10 days	\$131.22
Rabies Titer for Export[†]	7205E	● RTT 2 mL serum 	FAVN-OIE completed form must accompany sample submission.	6-7 weeks	\$168.70









Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Vaccine Titer (Quantitative), Canine(Cornell)	DPVAC	●● SST, RTT 1 mL serum 	FAVN-OIE completed form must accompany sample submission.	7-10 days	\$98.17
Core Vaccine Titer Panel (includes Rabies, non-export), Canine†	6160	● RTT 2 mL serum 	Ab titers for canine distemper, canine adenovirus-2, canine parvovirus, and rabies	5-6 weeks	\$185.87
Feline Infectious Peritonitis Virus (FIP) IFAT	25001570	● RTT 1 mL Feline effusions (ascites, peritoneal fluid) in red-top or preservative-free tube 	Indirect fluorescent antibody	3-4 days	\$96.17
Leptospira MAT, 5 Serovar†	6164	● RTT 0.5 mL serum 	Detection of antibodies to <i>Leptospira interrogans</i> serovars <i>pomona</i> , <i>hardjo</i> , <i>icterohaemorrhagiae/copenhageni</i> , <i>grippityphosa</i> , and <i>canicola</i>	Up to 10 days	\$109.57
Leptospira MAT, 5 Serovar & PCR†	6163	●● RTT, LTT For MAT: 0.5 mL serum in RTT; For PCR: 5 mL Urine; OR 1 mL EDTA whole blood (LTT); OR 1 gm tissue (e.g., kidney and/or liver); OR 1 mL aqueous humor; OR 0.5 mL CSF 	Detection of antibodies to <i>Leptospira interrogans</i> serovars <i>pomona</i> , <i>hardjo</i> , <i>icterohaemorrhagiae/copenhageni</i> , <i>grippityphosa</i> , and <i>canicola</i>	Up to 10 days	\$147.40
Leptospira Antibodies (IgM), Canine	6008	●● SST, RTT, LLT 0.5 mL serum or EDTA plasma or whole blood 	Detects IgM for the following canine serovars: <i>canicola</i> , <i>icterohaemorrhagiae</i> , <i>grippityphosa</i> , <i>pomona</i>	1-3 days	\$35.72










Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Feline Coronavirus (IFA)	CORONAIFA	●● SST, RTT 1 mL serum 	Detects IgG antibodies against Feline coronavirus. Positive results reported as a titer.	2-3 days	\$95.54
<i>Brucella canis</i> Titer IFA	BRUCIFA	●● SST, RTT 1 mL serum 		1-4 days	\$95.54
Neurology Infectious Disease Panel†	NEUROID	●● SST or RTT 2 mL serum 	<i>E canis</i> , <i>R rickettsii</i> , <i>B burgdorferi</i> (IFA); <i>Blastomyces</i> IgG, <i>Histoplasma</i> IgG, <i>Coccidioides</i> IgG/IgM; <i>Toxoplasma</i> IgG/IgM, <i>Cryptococcus neoformans</i> antigen, <i>Neosporium caninum</i> IgG/IgM, <i>A phagocytophilum</i> , <i>Aspergillus</i> IgG	3-5 days	\$542.49
Lyme Multiplex Test, Canine†	LYMEMPLX	● RTT 2 mL serum, separated 		5-7 days	\$123.73
<i>Anaplasma phagocytophilum</i> (IFA)	ANATITER	●● SST or RTT 1 mL serum 		2-3 days	\$50.95
<i>Aspergillus</i> Titer AGID†	ASP	●● SST or RTT 1 mL serum 	IgG	3-5 days	\$111.96
Australia Export Completion Panel†	AUSECPN	●● SST, RTT 2 mL serum 	<i>Brucella canis</i> , <i>Leishmania infantum</i> , <i>Leptospira interrogans</i> serovar <i>canicola</i> and <i>Ehrlichia canisthen</i> Appropriate FAVN forms must accompany submission. See online DOS for more comprehensive instructions.	9-12 days	\$167.07
<i>Brucella canis</i> Slide AGGL/AGID II†	6315	● RTT 0.5 mL serum, separated 		3-5 days	\$120.52
<i>Borrelia</i> (Lyme) (Western Blot)	LYMEWB	●● SST, RTT 1 mL serum		1-7 days (test performed on Mondays)	\$57.30









Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Neospora caninum</i> IFA†	NEO	●● SST or RTT 1 mL serum IgG, IgM 		3-5 days	\$118.38
Rabies Screen Titer (not for export)†	7205	● RTT 1 mL serum 		4-5 weeks	\$124.80
<i>Trypanosoma cruzi</i> (Chagas Disease) IFA Titert	TCRUIFA	●● SST, RTT 1 mL serum 		3-5 days	\$103.27
<i>Aspergillus</i> Galactomannan Antigen EIA†	ASPGUR	●● SST or RTT 0.8 mL serum (separated), CSF, urine, or BAL fluid 		2-5 days	\$153.71
<i>Bartonella</i> IFA†	6025	● RTT 2 mL serum 	Detects antibodies against <i>Bartonella henselae</i> , <i>B vinsonii</i> subsp <i>berhoffii</i> , and <i>B koehlerae</i>	Up to 4 weeks	\$147.40
Comprehensive (PCR and Serology) Vector Panel, Canine†	K9FCOMP	●● RTT, LTT 2 mL serum in RTT AND 2 mL whole blood in LTT 	<i>Babesia canis</i> , <i>B gibsoni</i> , <i>E canis</i> , <i>R rickettsii</i> (RMSF), <i>Bartonella henselae</i> , <i>B vinsonii</i> , <i>B koehlerae</i> (by IFA); Lyme, <i>Anaplasma</i> , <i>Ehrlichia</i> , and heartworm (by rapid assay); PCR for <i>Babesia</i> , <i>Bartonella</i> , <i>Anaplasma</i> , <i>Ehrlichia</i> , <i>Rickettsia</i> , and hemotropic <i>Mycoplasma</i> organisms	10 days	\$297.33
<i>Toxoplasma gondii</i> ELISA (Canine) & <i>Neospora caninum</i> IFA†	K9TOXOPL	●● SST, RTT 1 mL serum 	IgG and IgM Ab to <i>Toxoplasma gondii</i> . CSF and aqueous humor also accepted. Do not freeze.	7-10 days	\$125.19
Canine Distemper Virus Titer (Cornell)†	DISTEMP	●● SST or RTT 1 mL serum 	Testing also performed on CSF	7-10 days	\$109.57





Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Ehrlichia canis</i> (IFA)	EHRCH	●● SST, RTT 1 mL serum 	Detects IgG antibodies against <i>Ehrlichia canis</i> . Positive results reported as a titer.	2-3 days	\$50.95
FIV by Western Blot & ELISA†	FIVWB	●● RTT, LTT 1 mL serum or EDTA plasma 		4-7 days	\$148.34
<i>Histoplasma</i> Quantitative Antigen EIA†	HISTOANT	●●● SST, RTT, or LTT 1.2 mL serum or plasma, separated; or 0.8 mL CSF; or 0.5 mL urine or BAL fluid 		2-5 days	\$145.13
Canine Adenovirus Type 1 Antibody (VN)†	K9ADNE1	●● SST, RTT 0.5 mL serum 		4 days from run date	\$85.48
<i>Leishmania</i> Titer IFA†	LEISH	●● SST, RTT 1 mL serum 		2-5 days	\$115.17
Lyme (<i>Borrelia</i>) Titer IFA†	LYMEIFA	●● SST, RTT 1 mL serum 		2-3 days	\$95.91
<i>Pythium</i> Serology†	OOPYTHIUM	● RTT 0.5 mL serum 		3-5 days	\$147.27
Rabies Antibody Endpoint (Not for Export)†	7206	● RTT 1 mL serum 		4-5 weeks	\$155.25
<i>Rickettsia rickettsii</i> (RMSF) IFA†	RMSFIFA	● RTT 2 mL serum 		10 days	\$96.92





Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Feline Leukemia Virus (FeLV) (ELISA)	25000531	<ul style="list-style-type: none"> ● ● ● SST, RTT, LTT 0.5 mL serum or EDTA plasma or whole blood  		Next day by 8 AM	\$19.12
<i>Babesia</i> IFA†	BABESIAIFA	<ul style="list-style-type: none"> ● ● SST or RTT 2 mL serum  	<i>Babesia canis, B gibsoni</i>	10 days	\$113.57
Comprehensive (PCR and Serology) Vector Panel, Felinet	FELVECTPNL	<ul style="list-style-type: none"> ● ● RTT, LTT 2 mL serum AND 2 mL EDTA whole blood  	<i>Bartonella henselae, B vinsonii, B koehlerae</i> (by IFA); <i>Borrelia</i> (Lyme), <i>Anaplasma, Ehrlichia</i> , and heartworm (by rapid assay); PCR for <i>Babesia, Bartonella, Cytauxzoon, Anaplasma, Ehrlichia, Rickettsia</i> , and hemotropic <i>Mycoplasma</i> organisms	7-10 days	\$272.15
Canine Influenza Virus Titer†	K9FLUAB	<ul style="list-style-type: none"> ● RTT 0.5 mL serum  	HAI (antibody detection) specific for H3N8	7-10 days	\$99.12
Tick Panel Screening Test + Heartworm Antigen Test	25001389	<ul style="list-style-type: none"> 0.5 mL of serum or EDTA plasma or whole blood  	<i>Ehrlichia, Anaplasma, Borrelia</i> (Lyme) by rapid assay; Heartworm Ag by DiroCHEK® ELISA	8 AM next day	\$21.12
New Zealand Export Panel I†	NZECPN	<ul style="list-style-type: none"> ● ● SST, RTT 2 mL serum  	<i>Leptospira canicola</i> , occult heartworm	8-12 days	\$112.77
New Zealand Export Panel II†	NZECPN2	<ul style="list-style-type: none"> ● ● SST, RTT 2 mL serum  	<i>Babesia gibsoni, Brucella canis</i>	3-5 days	\$195.42
Canine Parvovirus Fecal Test	60002	<ul style="list-style-type: none"> ● Feces 1 gm fresh feces  		Next day by 5 PM	\$57.30

Clinical Pharmacology/Toxicology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Cyclosporine†	6335	<ul style="list-style-type: none"> LTT 0.5 mL EDTA whole blood 		3-5 days	\$128.45
Leflunomide Level†	LEFLUNO	<ul style="list-style-type: none"> RTT (no SST), LTT 0.5 mL serum or plasma (no SST) 		7-14 days	\$172.66
Levetiracetam (Kepra)†	6321	<ul style="list-style-type: none"> RTT 0.5 mL serum, separated (no plasma) 		3-5 days	\$153.71
Quantitative Copper (Blood or Tissue)†	6171	<ul style="list-style-type: none"> RTT, LTT, GRNTT 2 mL serum or whole blood (EDTA or heparin) or 50 mg wet tissue (fresh, frozen, fixed, or embedded in a paraffin block) 		7-10 days	\$194.61
Zinc, Serum/Plasma†	ZINC	<ul style="list-style-type: none"> RBTT (EDTA) 2 mL EDTA plasma, separated 	If possible utilize K2EDTA in RBTT. These tubes are designed for trace mineral testing and are free of zinc. Other tube types may contain zinc.	3-5 days	\$108.75
Anticoagulant (Rodenticide)†	ANTICOAG	<ul style="list-style-type: none"> LTT 3 mL EDTA whole blood 	Includes brodifacoum, bromadiolone, chlorophacinone, difenacoum, difethialone, diphacinone, and warfarin	6-12 days	\$197.90
Bromide	1025	<ul style="list-style-type: none"> SST, RTT 0.5 mL serum 		1-2 days	\$142.70
Digoxin†	6320	<ul style="list-style-type: none"> LTT, RTT (no SST) 0.5 mL serum (no SST) or EDTA plasma 		3-5 days	\$122.16

Clinical Pharmacology/Toxicology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Drugs of Abuse†	DOA	●● SST, RTT 10 mL serum or urine (10 gm liver or kidney also accepted)	Amphetamine, cocaine, benzoylecgonine, norcocaine, ephedrine, lysergic acid diethylamide, (LSD), methamphetamine, MDMA, nicotine, phentermine, psilocin, delta-9-THC, THC-OH, THC-COOH, methadone, morphine, HU-210, JWH-018	10-14 days	\$292.53
LC/Mass Spectrometry Toxicology Screen†	LCMS	● RTT 20 mL fluid; 50 gm tissue; 50 gm feed; OR 50 gm source material 	Toxicology screen for unknown exposure	2-3 weeks	\$248.33
Blood Lead†	6240	●● LTT, GRNTT 0.5 mL EDTA or heparinized whole blood 		5-7 days	\$114.10
Zonisamide†	6322	●● LTT, RTT (no SST) 0.5 mL serum (no SST) or EDTA plasma		6-10 days	\$153.71
Keppra Levels Peak & Trough Samples†	KEPPRA	● RTT (no SST) 2 samples of 0.5 mL serum, separated (no plasma) each		6-10 days	\$235.73
Phenobarbital	1031	● RTT 0.5 mL serum, separated 		1-2 days	\$44.60
Gabapentin†	GAB	●● LTT, RTT (no SST) 0.5 mL EDTA plasma or serum (no SST)		3-5 days	\$142.99
Toxic Heavy Metal Panel†	25001395	●●● LTT, GRNTT, BTT 5 mL EDTA, heparinized, or citrate whole blood; 3 gm tissue (liver or kidney) 	Quantitative analysis of arsenic, cadmium, lead, thallium, mercury, selenium	6-12 days	\$218.96

Pathology

Histopathology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Histopathology Simple Single Site	ZN100	Formalin-fixed tissue	<ul style="list-style-type: none"> - Single site - No more than 10 cm in any dimension - Not for dermatologic disease or with other criteria of Complex Level 1 or 2 	3-7 days	\$90.95
Histopathology Simple Multiple Sites (Up to 5 Sites)	ZN101	Formalin-fixed tissue	<ul style="list-style-type: none"> - Simple samples from up to 1-5 unique sites - No more than 10 cm in any dimension 	3-7 days	\$133.75
Histopathology Complex Level 1	ZN102	Fresh or formalin-fixed tissue	<ul style="list-style-type: none"> - Skin biopsies for dermatologic disease (without dermatologist consult) - Multiple gastrointestinal biopsies - Bone marrow - Whole organs/major components of organs: eye, spleen, testicles, uterus/ovaries, mammary chains (2 or more connected glands), amputated digit/nail or tail, liver lobe, or lung lobe - Large samples (>10cm in any dimension) - Individual teeth or other small samples containing bone - Simple samples from 6-15 unique sites - Formalin-fixed whole bodies <300 g 	3-9 days	\$168.25
Histopathology Complex Level 2	ZN103	Fresh or formalin-fixed tissue	<ul style="list-style-type: none"> - Mandible, maxilla (all jaws, regardless of size) - Amputated limb - Other large/atypical samples containing bone - Multiple Complex Level 1 samples - Simple samples from >15 sites. 	Varies by nature of sample	\$213.93



See the complete Directory of Services online at zoetisreflabs.com/DOS

Pathology (cont'd)

Specialty Histopathology Tests



Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Dermatopathology Panel	ZN400	<ul style="list-style-type: none"> - Formalin-fixed skin biopsies - Detailed history and clinical images are required - Images can be sent to: referencelabpathology@zoetis.com 	<ul style="list-style-type: none"> - Includes Histopathology & Dermatologist Consult - Addendum with case management guidance from boarded dermatologist 	Dermatopathology: 3-7 days. Consult provided after dermatopathology.	\$184.67
Liver Panel	ZN500	<ul style="list-style-type: none"> - At least 3 x 1 cm 14G needle biopsies, 1+ wedge biopsy (ideally multiple), or other larger liver sample - Samples for histologic review should be formalin-fixed - Samples for copper quantification may be fresh or formalin-fixed - Detailed history is required - Images can be sent to: referencelabpathology@zoetis.com 	<ul style="list-style-type: none"> - Biopsy report by liver specializing pathologist - H&E, connective tissue, and copper stains - Copper quantification 	Histopathology: 3-7 days. Quantitative copper analysis: Additional 3-7 days	\$299.28
Liver Panel & Medical Consult	25001436	<ul style="list-style-type: none"> - Biopsy report by liver specializing pathologist - H&E, connective tissue, and copper stains - Copper quantification 	<ul style="list-style-type: none"> - Biopsy report by liver specializing pathologist - H&E, connective tissue, and copper stains - Copper quantification - Addendum with case management guidance from boarded internist with liver expertise 	Histopathology: 3-7 days Copper quantification and consult: Up to 21 days	\$413.77

Pathology (cont'd)

Specialty Histopathology Tests (cont'd)


Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Nephropathology (Renal) Panel	ZN25000541	<ul style="list-style-type: none"> - 1 or more 14G to 18G core needle biopsy, a wedge biopsy, or other larger renal sample (eg, nephrectomy) - Detailed history is required - Images can be sent to: referencelabpathology@zoetis.com 	<ul style="list-style-type: none"> - Biopsy report by kidney specializing pathologist - H&E, PAS, trichrome, Jones methenamine silver stains (Congo Red also included, if indicated) 	3-7 days	\$230.05

Other Histopathology Tests

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
C-Kit Mutation Assay†	6541	Paraffin block (preferred); Formalin-fixed tissue; fresh tissue; or minimum 6 unstained positively charged slides. Must also include the histopathology report. 	C-Kit PCR Exons 8 & 11 This test must be ordered on an existing histopathology accession by calling customer service at 1-888-965-9652 . See specific recommendations made by the pathologist in the existing pathology report. Please include histopathology report.	10-14 days	\$206.12
Melanoma Prognostic Panel, Canine†	6426	Fixed tissue, paraffin tissue block, or unstained slides 	IHC for Ki-67 and prognostic interpretation This test must be ordered on an existing histopathology accession. See specific recommendations made by the pathologist in the existing pathology report.	7-10 days	\$117.70

Pathology (cont'd)

Other Histopathology Tests (cont'd)









Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Complete Muscle Profile†	CM	0.5 x 0.5 x 1.0 cm biopsy from affected, but not end-stage muscle from open biopsy procedure. 	Wrap specimen in a saline moistened (not dripping) gauze sponge, place in dry water-tight container, keep refrigerated, and ship on cold packs. Collect a second smaller biopsy and immersion fix in 10% neutral buffered formalin. For optimal results, keep specimens chilled and ship immediately.	7-10 days	\$221.10
Lymphoid Neoplasia PCR/ PARR (T&B cell) from Histopathology Sample†	6544	Biopsy sample	PARR for lymphoma and lymphoid leukemia This test must be ordered on an existing histopathology accession by calling customer service at 1-888-965-9652 . See specific recommendations made by the pathologist in the existing pathology report.	10-14 days	\$279.87
Mast Cell Tumor Prognostic Panel, Canine†	6540	Paraffin block, fresh or formalin-fixed tissue, or unstained, non-coverslipped, positively charged slides (minimum 10 slides)	IHC for C-Kit and Ki-67, PCR for C-Kit Exons 8 and 11, and AgNOR special staining This test must be ordered on an existing histopathology accession by calling customer service at 1-888-965-9652 . See specific recommendations made by the pathologist in the existing pathology report. Please include a copy of the original biopsy report.	10-14 days	\$348.00
Melanoma Diagnostic Panel, Canine†	6425	Fixed tissue, paraffin tissue block, or unstained slides	IHC for melan A, PNL-2, TRP-1, and TRP-2 This test must be ordered on an existing histopathology accession by calling customer service at 1-888-965-9652 . See specific recommendations made by the pathologist in the existing pathology report. Please include a copy of the original biopsy report. Numeric results will be complete within 3-5 days. An interpretation will be provided an additional 1-10 days post numeric result.	3-5 days	\$153.71
Copper, Zinc, Iron Levels on Hepatic Tissue†	COPZNCIRON	Liver may be fresh (minimum 50 mg wet tissue), frozen, fixed, or embedded in a paraffin block		1-2 weeks	\$134.81

Pathology (cont'd)

Cytopathology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Cytopathology Simple/ Single Site	7120	<ul style="list-style-type: none"> LTT Unstained, air-dried slides; fluid from a mass in LTT <ul style="list-style-type: none"> Slides Liquid 	For the following sample types: <ul style="list-style-type: none"> Single lesion/tissue (FNA, swab, impression smear) Single or multiple lymph nodes Single joint (slides only) Fluid from mass lesion (cytology only) 	1-2 days	\$96.30
Cytopathology Complex	7201	<ul style="list-style-type: none"> LTT Unstained, air-dried slides and/or fluid in LTT <ul style="list-style-type: none"> Slides Liquid 	For the following sample types: <ul style="list-style-type: none"> Multiple lesions/tissues (FNA, swab, impression smear) Joint fluid in LTT (single or multiple) Multiple joints (slides only) Cavitary effusion (single or multiple) CSF Bone marrow BAL/TTW Other fluids/lavage 	1-2 days	\$130.52
Lymphoid Neoplasia PCR/PARR (T&B Cell) from Blood, Fluid, Cytopathology Sample†	6542	<ul style="list-style-type: none"> LTT May be performed on blood, bone marrow, cavity fluid, CSF, or tissue aspirate. For blood, bone marrow, or other fluids: 1 mL in EDTA (LTT). For tissue aspirate: aspirate into EDTA or submit cytology preparation (slides).		7-12 days	\$279.87

Clinical Immunology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
2M Antibody Masticatory Muscle Myositis†	6403	● RTT 1 mL serum 		7-12 days	\$204.37
Acetylcholine Receptor Antibody/Myasthenia Gravis†	6402	● RTT 1 mL serum 		7-12 days	\$204.37
Complete Amino Acid Profile†	COMPAA	● GRNTT 1 mL separated heparinized plasma 		7-9 days	\$235.73
IgA, Canine	IGA	● ● SST, RTT 2 mL serum 		3-5 days	\$57.30
Anti-Nuclear Antibody (IFA)	ANA	● ● SST, RTT 0.5 mL serum 	Detects antinuclear IgG antibodies. Positive results reported as a titer.	2-5 days	\$50.95
Protein Electrophoresis	6345	● ● SST, RTT, sterile leakproof container 0.5 mL serum, >10 mL urine 		3-6 days	\$70.06
IgG, IgA, IgM, Canine†	IGGAM	● RTT 1 mL separated serum 		7-10 days	\$161.18
IgG, Canine†	IGG	● RTT 1 mL separated serum 		7-10 days	\$94.70
Allergy Panel, Canine/Feline†	25001385	● ● SST, RTT 3 mL serum		10-14 days	\$332.88









Clinical Immunology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Immunofixation Free Lt Chain Add-on†	25001431	<ul style="list-style-type: none"> ○ SST, sterile leakproof container 1 mL serum, 5-10 mL urine 	Must be ordered in addition to Protein Electrophoresis.	3-4 days	\$152.27
Protein Electrophoresis Immunofix†	25001430	<ul style="list-style-type: none"> ● ○ RTT, WTT sterile leakproof container 0.5 mL serum, >=10 mL urine in WTT/Sterile leakproof container 		5-7 days	\$203.79



See the complete Directory of Services online at zoetisreflabs.com/DOS





Genetic

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Degenerative Myelopathy†	DEGMYE	<ul style="list-style-type: none"> ● LTT 3-7 mL EDTA whole blood  Blood  Tissue 	If submitting tissue, do not place in formalin. Place tissue in sterile container/bag, double bag, freeze, and ship with ice pack. UMC CANINE DM DNA TESTING & RESEARCH FORM required.	2 weeks	\$141.11
Fanconi Syndrome Test†	FANCONI	<ul style="list-style-type: none"> ● Urine 3 mL urine (frozen)  		1-4 weeks	\$141.92
Centronuclear Myopathy†	CNM	Interdental/GUM brushes supplied by owner or Cytology Brush supplied by VGL at no additional charge		10-14 days	\$185.25
Methylmalonic Aciduria (Cobalamin Malabsorption)†	COBMAL	<ul style="list-style-type: none"> ● Urine 3 mL urine (frozen)  		5-7 days	\$134.81
Cystinuria Nitroprusside†	CYSTINURIA	<ul style="list-style-type: none"> ● Urine 3 mL urine (frozen)  		7-10 days	\$125.87
Congenital (CKCS) Macrothrombocytopenia†	CKCSMACRO	<ul style="list-style-type: none"> ● LTT 1 mL EDTA whole blood AND 3-5 unstained blood smears  Blood  Slides 	β1-tubulin	7-10 days	\$210.37
Pyruvate Kinase Deficiency DNA Test†	PYRUDNA	<ul style="list-style-type: none"> ● LTT 1-2 mL EDTA whole blood or cheek swab  		1-4 weeks	\$149.12

Nutrition

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Selenium†	SELENIUM	<ul style="list-style-type: none"> ● ● LTT, RTT 0.5 mL EDTA whole blood preferred; 0.5 mL serum or EDTA plasma also accepted ❄ Serum/Plasma 📄 Whole Blood 		3-5 days	\$100.19
Trace Element Screen†	TMINPAN	<ul style="list-style-type: none"> ● RBTT 2 mL serum from RBTT ❄ 	Ca, Cu, Fe (unbound), Mg, Na, K, Phos, Zn	5-7 days	\$88.16
Taurine (Plasma and Whole Blood)†	TAUR	<ul style="list-style-type: none"> ● GRNTT 1 mL separated heparinized plasma AND 2 mL heparinized whole blood 📄 		10 days	\$195.42
Taurine:Whole Blood†	TAURWB	<ul style="list-style-type: none"> ● GRNTT 2 mL heparinized whole blood 📄 		10 days	\$147.27
Taurine:Plasma†	TAURPLAS	<ul style="list-style-type: none"> ● GRNTT 1 mL separated heparinized plasma 📄 		10 days	\$147.27
Thiamine (Vitamin B1)†	THIAMINE	<ul style="list-style-type: none"> ● ● GRNTT, LTT 3 mL EDTA or heparinized whole blood ❄💡 		5-7 days	\$122.16

Nutrition (cont'd)







Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Vitamin E [†]	2331	<p>●●●● SST, RTT, LTT, GRNTT</p> <p>1.5 mL serum or heparinized or EDTA plasma, separated</p> 		3-5 days	\$100.19
Vitamin E and Mineral Panel [†]	VITEMINE	<p>●●● RTT, LTT, RBTT</p> <p>1 mL EDTA whole blood, 1 mL RBTT serum (separated), 1 mL serum no-additive RTT (separated and protected from light)</p> 	Vitamin E, Ca, Cu, Fe, Mg, Phos, Zn, Se	8-10 days	\$144.89
Carnitine [†]	CARNITINE	<p>●○ GRNTT, WTT</p> <p>0.5 mL heparinized plasma (separated) OR 5 mL urine (frozen)</p> 		15-21 days	\$196.49
Vitamin A [†]	VITA	<p>●●●● SST, RTT, LTT, GRNTT</p> <p>1.5 mL serum or EDTA or heparinized plasma, separated</p> 		3-5 days	\$100.19











See the complete Directory of Services online at zoetisreflabs.com/DOS

Add-Ons








You can quickly and easily add any of the tests below to an existing order, provided sufficient sample remains at the lab and is within stability guidelines. If you order electronically, you can place add-on orders through your PIMS or the ZoetisDx portal. If you use paper requisition forms, call customer service at **1-888-965-9652**.

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Anaerobic ID Add-on	A7157	○ WTT/sterile tube Culture swab or fluid/tissue in WTT, anaerobic transport medium, or sterile tube		2-7 days	\$32.45
Bile Acid Pre/Resting Add-on	A6180	● ● SST, RTT 0.5 mL serum 		1-2 days	\$23.79
ZRL Premium Lipase (DGGR Lipase) Add-on	A1030	● ● SST, RTT 0.5 mL serum 		8 AM next day	\$10.65
FELV/FIV Rapid Test Add-on	A6005	● ● ● SST, RTT, LTT 0.25 mL of serum or EDTA plasma or whole blood 	Lateral flow assay with confirmation of positive FeLV results by ViraCHEK® ELISA	8 AM next day	\$13.67
Fructosamine Add-on	A1024	● ● SST, RTT 0.5 mL serum 		1-2 days	\$19.16
Heartworm Antigen Test Add-on	A6001	● ● ● SST, RTT, LTT 0.25 mL of serum or EDTA plasma 	DiroCHEK® ELISA	8 AM next day	\$9.47
SDMA Add-on	AN7037	● ● SST, RTT 0.5 mL serum 		8 AM next day	\$17.10





Add-Ons (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Tick Panel Screening Test Add-on	A6002	<ul style="list-style-type: none"> ● ● ● SST, RTT, LTT 0.25 mL of serum or EDTA plasma or whole blood 	Lateral flow assay for <i>Ehrlichia</i> , <i>Anaplasma</i> , <i>Borrelia</i> (Lyme)	8 AM next day	\$13.47
Triglycerides Add-on	A1023	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL of serum or EDTA plasma or whole blood 		8 AM next day	\$7.59
TLI Canine Add-on	A6206	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum 		1-2 days	\$45.49
Cortisol Add-on	A5101	<ul style="list-style-type: none"> ● RTT 0.5 mL serum 		8 AM next day (1-2 days for feline/equine samples)	\$13.67
Direct Coombs Add-on	A2701	<ul style="list-style-type: none"> ● LTT 2 mL EDTA whole blood 		1-2 days	\$54.11
Feline Heartworm Antibody (ELISA) Add-on†	A6007	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum 		5-7 days	\$114.16
Free T4 Add-on	A5103	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum 		8 AM next day	\$16.83
<i>Giardia</i> Antigen Test, Canine/Feline Add-on	A2321	<ul style="list-style-type: none"> ● Feces 3 gm fresh feces 		Next day by 5 PM	\$11.18

Add-Ons (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Leptospira Antibodies (IgM), Canine Add-on	A6008	<ul style="list-style-type: none"> ● ● ● SST, RTT, LTT 0.5 mL serum or EDTA plasma or whole blood  		1-3 days	\$30.36
Phenobarbital Add-on	A1031	<ul style="list-style-type: none"> ● RTT 0.5 mL serum  		1-2 days	\$37.91
Total T4 Add-on	A5102	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum  		8 AM next day	\$14.41
TSH Add-on	A5106	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum  		8 AM next day (1-2 days for feline samples)	\$13.67
Urine Culture and Sensitivity Add-on	A7140	<ul style="list-style-type: none"> ○ ● WTT/sterile tube, Urine 1 mL urine in WTT or sterile container  		1-5 days	\$59.12
Urine Protein:Creatinine Ratio Add-on	A1555	<ul style="list-style-type: none"> ○ ● WTT, Urine 1 mL urine in sterile container  	Urine Creatinine, Urine Protein, Urine Protein:Creatinine Ratio	8 AM next day	\$14.41
Urine Cortisol:Creatinine Ratio Add-on†	A3310	<ul style="list-style-type: none"> ○ ● WTT, Urine 1 mL urine in sterile container  		8 AM next day (1-2 days for feline samples)	\$59.55
Coccidioides Titer AGID Add-on	A6129	<ul style="list-style-type: none"> ● SST 1 mL serum 	IgG, IgM Titers	3-5 days	\$86.14

Add-Ons (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Urinalysis Add-on	25001700	<ul style="list-style-type: none"> ● Urine 5 mL urine in YTT/sterile container 			\$14.41
Urinalysis w/ reflex UPC Add-on	25001701	<ul style="list-style-type: none"> ● Urine 5 mL urine in YTT/sterile container 	An automatic Urine Protein:Creatinine Ratio will be performed when proteinuria is present.		\$19.77
Free (Ionized) Calcium (Mich. State Univ.) Add-on†	A6304	<ul style="list-style-type: none"> ● RTT 0.75 mL serum, separated 	Overnight fasting recommended for best results. Heparinized plasma is suitable for analysis, but must be indicated on form.	3-6 days	\$70.62
Reticulocyte Count Add-on	25000525	<ul style="list-style-type: none"> ● LTT 0.5 mL EDTA whole blood 		8 AM next day	\$15.76



See the complete Directory of Services online at zoetisreflabs.com/DOS



Large Animal
















Bundles at a glance

Our bundles equip you with the information you need to help you provide the best medicine for your patients.

Components	Test Code	CBC	ClinChem	Endocrine	UA	Reflex UPC	Price
Wellness	EQWP	Basic	ClinChem 24: Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Mg, TP, Alb, Glob, Alb:Glob Ratio, AST, GGT, TBili, CK, CO2, Anion Gap, Phos, GLDH, ALP, Chol				\$52.06
Exercise	25000489		Exercise ClinChem: Glu, Na, K, Cl, BUN, Creat, BUN:Creat Ratio, AST, CK		✓		\$27.50
Metabolic	25000494		Glu, BUN, Creat, BUN:Creat Ratio, TP, Chol, Trig	Endogenous ACTH, Insulin			\$97.03
Renal I	25000497	Comp	Renal ClinChem: Glu, Na, K, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, CO2, Anion Gap, TP, Alb		✓	✓	\$63.83
Renal II	25000498		Renal ClinChem: Glu, Na, K, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, CO2, Anion Gap, TP, Alb		✓	✓	\$34.88

Comp = Comprehensive

Bundles

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Large Animal Wellness Bundle	EQWP	 1 mL serum; 2 mL EDTA whole blood, air-dried slides  Blood  Slides	Basic CBC, Large Animal ClinChem 24	1-2 days	\$52.06
Large Animal Exercise Bundle	25000489	 SST, RTT, Urine 1 mL serum; 5 mL urine in WTT/sterile container 	Exercise ClinChem + UA	1-2 days	\$27.50
Large Animal Metabolic Bundle	25000494	 SST, RTT 2 mL serum 	Glu, BUN, Creat, BUN:Creat Ratio, TP, Chol, Trig, Endogenous ACTH, Insulin	1-3 days	\$97.03
Large Animal Pre-Anesthetic Bundle	25000495	 SST, RTT, LTT 1 mL serum; 2 mL EDTA whole blood 	Basic CBC + Pre-Anesthetic ClinChem	1-2 days	\$29.22
Large Animal Renal Bundle I	25000497	 SST, RTT, LTT, Urine 2 mL EDTA whole blood; 1 mL serum; 5 mL urine in WTT/sterile container 	Basic CBC + Renal ClinChem, UA, Reflex UPC	1-2 days	\$63.83
Large Animal Renal Bundle II	25000498	 SST, RTT, Urine 2 mL EDTA whole blood; 1 mL serum; 5 mL urine in WTT/sterile container 	Renal ClinChem, UA, Reflex UPC	1-2 days	\$34.88

Hematology

Choose the CBC option that's best for your patients.







At Zoetis Reference Laboratories, our CBCs are grounded in quality standard operating procedures developed by a team of board-certified clinical pathologists. Expert technicians review all results manually and run additional tests if needed.

	Basic Test code: 7000	Comprehensive Test code: 8000
	Recommended for wellness screens in non-geriatric, healthy patients	Recommended for sick or geriatric patients
WBC, HGB, RBC, HCT, MCV, MCHC, CHCM, RDW, 5-Part WBC Differential, Platelet Count, MPV, Plasma Protein, Fibrinogen, Plasma Protein:Fibrinogen Ratio	✓	✓
Perform if indicated based on automated results: Spun PCV, Manual WBC Differential, Visual Platelet Estimate, and/or Reticulocyte Count	✓	✓
Complete manual smear review by technician		✓
Perform if indicated based on automated results: complete pathologist review		✓







Tests listed above are intended for horses and other domestic hoofstock, including: cattle, sheep, goats, pigs, camels, llamas, and alpacas.

Test requirements: 2 mL EDTA whole blood, air-dried slides

Hematology (cont'd)







Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Basic CBC, Large Animal (Healthy Animals)	7000	<ul style="list-style-type: none"> LTT 2 mL EDTA whole blood, air-dried slides 	RBC, HGB, HCT, MCV, MCHC, CHCM, RDW, WBC, WBC Differential, PLT, MPV, Plasma Protein, Fibrinogen, Plasma Protein:Fibrinogen Ratio. If results indicate, will include spun PCV, Manual Differential, and/or Reticulocyte Count.	1-2 days	\$26.79
Comprehensive CBC, Large Animal	8000	<ul style="list-style-type: none"> LTT 2 mL EDTA whole blood, air-dried slides 	CBC with blood smear review by technician. If results indicate, will include a PCV, Manual Differential, Reticulocyte Count, and/or Pathologist Review.	1-2 days. Pathologist Review may extend turnaround time to end of day.	\$38.56
Comprehensive CBC Upgrade, Large Animal	25001590	<ul style="list-style-type: none"> LTT 2 mL EDTA whole blood, air-dried slides		1-2 Days. Pathologist Review may extend turnaround time to end of day.	\$11.77
Pathologist Review	PATH	<ul style="list-style-type: none"> LTT 2 mL EDTA whole blood, air-dried slides	Blood smear review by Pathologist. May be added on to Comprehensive CBC only.	1-2 days	\$80.25
Technician Blood Smear Review	SMREV	<ul style="list-style-type: none"> LTT 2 unstained, air-dried slides; include 2 mL EDTA whole blood, if possible (for in-lab preparation of smears)  Blood  Slides	Morphologic review of RBC, WBC, and PLT, Manual WBC Differential, and visual PLT estimate. If indicated based on morphologic findings, a Pathologist Review will be included.	8 AM next day	\$17.80
Reticulocyte	RETIC	<ul style="list-style-type: none"> LTT 1 mL EDTA whole blood 	% Reticulocyte, RBC, Reticulocyte	8 AM next day	\$18.54
Buffy Coat	6295	<ul style="list-style-type: none"> LTT 1 mL EDTA whole blood 		1-2 days	\$108.25

Hematology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Iron Profile†	IRONPRO	<ul style="list-style-type: none"> ● RTT, GRNTT 0.5 mL serum or heparinized plasma  	Fe, TIBC, and Ferritin	7-10 days	\$142.99
PT & aPTT	4101	<ul style="list-style-type: none"> ● BTT 0.5 mL separated citrated plasma  		8 AM next day	\$25.48
Antithrombin (AT3)†	ANTTHRO3	<ul style="list-style-type: none"> ● BTT 1 mL separated citrated plasma (frozen)  		3-5 days	\$114.64
Protein C†	PROC	<ul style="list-style-type: none"> ● BTT 1 mL separated citrated plasma (frozen)  		2-3 days	\$118.92
Low Molecular Weight Heparin†	LMWHEPARIN	<ul style="list-style-type: none"> ● BTT 0.5 mL citrated plasma (frozen)  	Anti-Xa activity of LMWH	2-3 days	\$118.38
Erythropoietin (EPO)†	50227	<ul style="list-style-type: none"> ● SST or plasma separator tube 1-2 mL plasma or serum from SST or plasma separator tube (EDTA plasma not acceptable)  		5-7 days	\$185.25






Chemistry

ClinChem Panels

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
ClinChem 18, Large Animal	CHEM18EQ	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum  	Glu, Na, K, Na:K Ratio, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Mg, TP, Alb, Alb:Glob Ratio, AST, GGT, TBili, CK	1-2 days	\$16.08
Large Animal ClinChem 24	CHEM24EQ	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum  	ClinChem 18 + CO2, Anion Gap, Phos, GLDH, ALP, Chol	1-2 Days	\$22.27
Exercise ClinChem, Large Animal	25000509	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum  	Glu, Na, K, Cl, BUN, Creat, BUN:Creat Ratio, AST, CK	1-2 days	\$10.88
Large Animal Liver ClinChem	25000511	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum  	Glu, BUN, TP, Alb, Glob, Alb:Glob Ratio, AST, GLDH, GGT, TBili, Chol	1-2 days	\$10.88
Large Animal Renal ClinChem	25000485	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum  	Glu, Na, K, Cl, BUN, Creat, BUN:Creat Ratio, Ca, Phos, CO2, Anion Gap, TP, Alb	1-2 days	\$12.54
Large Animal Pre-Anesthetic ClinChem	25000486	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum  	Glu, Na, K, Creat, Ca, TP, GLDH, TBili	1-2 days	\$10.88









Chemistry (cont'd)

Individual Chemistry Tests

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Free (Ionized) Calcium	IONISTAT	<ul style="list-style-type: none"> ● GRNTT Lithium heparinized whole blood. Ensure that the sample is treated anaerobically and that the correct ratio of heparin anticoagulant to blood is achieved (do not underfill GRNTT). 		1-2 days	\$50.95
Free (Ionized) Calcium (Mich. State Univ.)†	6304	<ul style="list-style-type: none"> ● RTT 0.5 mL serum, separated Ensure that the sample is treated anaerobically. 		3-6 days	\$79.86
LDH†	LDH	<ul style="list-style-type: none"> ● ● RTT, GRNTT 0.5 mL separated serum or heparinized plasma (fluid and urine also accepted)		3-5 days	\$87.35
Osmolality†	OSMOLAL	<ul style="list-style-type: none"> ● ● RTT, GRNTT 1 mL separated heparinized or EDTA plasma 		2-3 days	\$95.91
Serum Amyloid A (SAA), Equinet†	SAA	<ul style="list-style-type: none"> ● RTT 3 mL serum 		2-3 days	\$117.31
Albumin	1001	<ul style="list-style-type: none"> ○ ● SST, RTT 0.5 mL serum 		1-2 days	\$7.10










Chemistry (cont'd)

Individual Chemistry Tests (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
High-Sensitivity Cardiac Troponin I†	TROPONIN	● ● SST, RTT 0.5 mL non-hemolyzed serum (frozen) 		4-5 days	\$136.57
ALP (Alk Phos)	1003	● ● SST, RTT 0.5 mL serum 		1-2 days	\$8.75
AST	1005	● ● SST, RTT 0.5 mL serum 		1-2 days	\$7.49
BUN	1006	● ● SST, RTT 0.5 mL serum 		1-2 days	\$8.75
Chloride	1007	● ● SST, RTT 0.5 mL serum 		1-2 days	\$7.49
Calcium (Total Calcium)	1008	● ● SST, RTT 0.5 mL serum 		1-2 days	\$8.75
Cholesterol	1009	● ● SST, RTT 0.5 mL serum 		1-2 days	\$8.75
Creatinine	1010	● ● SST, RTT 0.5 mL serum 		1-2 days	\$8.75






Chemistry (cont'd)

Individual Chemistry Tests (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Creatinine Kinase (CK)	1011	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum   		1-2 days	\$7.10
GLDH	1040	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum  		1-2 days	\$7.49
GGT	1013	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum  		1-2 days	\$7.49
Glucose	1014	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum  		1-2 days	\$8.75
Magnesium	1016	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum  		1-2 days	\$7.49
Phosphorus	1017	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum  		1-2 days	\$8.75
Potassium	1018	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum  		1-2 days	\$8.75
Sodium	1019	<ul style="list-style-type: none"> ● SST, RTT 0.5 mL serum  		1-2 days	\$7.49

Chemistry (cont'd)






Individual Chemistry Tests (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
CO2	1020	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum  		1-2 days	\$7.49
Total Bilirubin	1021	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum   		1-2 days	\$7.49
Total Protein	1022	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum  		1-2 days	\$7.49
Lipid Profile†	LIPIDPROF	<ul style="list-style-type: none"> ● ● ● SST, RTT, GRNTT 1 mL serum OR 1 mL heparinized plasma 	Chol, Trig, Chol HDL, Chol LDL, Non-HDL Chol, HDL:Chol Ratio	2-4 days	\$125.39
Triglycerides	1023	<ul style="list-style-type: none"> ● ● SST, RTT 0.5 mL serum  		1-2 days	\$8.93











See the complete Directory of Services online at zoetisreflabs.com/DOS

Urine

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Urinalysis	3000	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 5 mL urine in a sterile container 		8 AM next day	\$16.95
Urine Protein:Creatinine Ratio	1555	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 1 mL urine in sterile container 		8 AM next day	\$16.95
Urinalysis with Reflex UPC	25000521	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 5 mL urine in sterile container 	An automatic Urine Protein:Creatinine Ratio will be performed when proteinuria is present.	8 AM next day	\$23.25
Urinalysis with Reflex UPC upgrade	25001702	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 5 mL urine in sterile container 	Add to ZRL Bundles including Urinalysis		\$6.30
Urine Electrolyte Panel (for Fractional Excretion) [†]	FREXCRE	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 2 mL urine in a no-additive tube or leakproof plastic container 	Cl, Creat, K, Na	3-5 days	\$118.38








Endocrinology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Cortisol (Baseline/Resting)	5101	●● SST, RTT 0.5 mL serum 		8 AM next day (1-2 days for feline/ equine samples)	\$16.08
Cortisol AM/PM, Equine	AMPMEC	●● SST, RTT 2 mL serum (2 samples required) 		1-2 days	\$30.57
Progesterone	5105	● RTT 0.5 mL serum 		8 AM next day	\$44.60
Total T4	5102	●● SST, RTT 0.5 mL serum 		8 AM next day	\$16.95
Free T4	5103	●● SST, RTT 0.5 mL serum 		8 AM next day	\$19.80
Free T4 by ED†	6330	●● SST, RTT 0.5 mL serum 		5-7 days	\$120.52
Total T3 (TT3), Equine	TT3	●● SST, RTT 1 mL serum 		1-2 days	\$38.18
Thyroid Panel	EQTHYROID	●● SST, RTT 1.5 mL serum 	TT4, TT3, free T4	1-2 days	\$45.85








Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Estrone sulfate, Equine†	25001386	<ul style="list-style-type: none"> ●● SST, RTT PREFERRED: 1 mL serum; 1 mL separated EDTA or heparinized plasma also accepted 		3-10 days	\$107.68
Endogenous ACTH	ENDOACTH	<ul style="list-style-type: none"> ● LTT 1 mL EDTA plasma; separated within 4 hours and shipped frozen ❄ 		1-3 days	\$52.27
Equine Metabolic Panel (ACTH, Insulin, Glucose)	EQCUSHING	<ul style="list-style-type: none"> ●●●● SST, RTT, LTT, GRYTT 1 mL EDTA plasma (separated within 4 hrs and frozen) AND 1 mL serum (separated within 1 hr) ❄ Plasma 📄 Serum 	Endogenous ACTH, Insulin, Glu If unable to separate serum within 1 hr: 1 mL EDTA plasma (separated within 4 hrs and frozen), 1 mL serum, AND 2 mL GRYTT whole blood	1-3 days	\$91.71
Metabolic Syndrome Panel I, Equine†	CUSHING	<ul style="list-style-type: none"> ●● LTT AND RTT 1 mL separated EDTA plasma, frozen, AND 1 mL separated serum ❄ 	Endogenous ACTH, Insulin, Glu, T4	1-3 days	\$191.14
TRH Stimulation Test for Equine Cushing's/PPID	2222	<ul style="list-style-type: none"> ● LTT 1 mL EDTA plasma for each sample 📄 	Baseline and post-TRH ACTH	1-4 days	\$101.90
Insulin-like Growth Factor†	ILGROFAC	<ul style="list-style-type: none"> ●●● SST, RTT, LTT 0.5 mL serum or EDTA plasma 📄 		5-10 days	\$139.78

Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Vitamin D (25-hydroxyvitamin D3)†	VITD	<ul style="list-style-type: none"> ● ● ● RTT, LTT, GRNTT 0.5 mL serum or EDTA or heparinized plasma, separated   	Overnight fasting recommended for best results. Centrifuge samples and separate serum/plasma within an hour of collection. Protect sample from light.	4-10 days	\$130.54
GCT Panel II, Equine†	EQGCTPII	<ul style="list-style-type: none"> ○ ● SST, RTT 3 mL serum  	Anti-Müllerian Hormone, Inhibin, Testosterone	5-7 days	\$216.50
Anti-Müllerian Hormone, Equine†	ANTIMULL	<ul style="list-style-type: none"> ○ ● SST, RTT 1 mL serum  		3-14 days	\$134.82
Pregnancy Panel, Equine†	25001387	<ul style="list-style-type: none"> ○ ● SST, RTT 3 mL serum (separated)  	Estrone Sulfate, Pregnant Mare Serum Gonadotropin, Progesterone	3-10 days	\$145.46
Pregnant Mare Serum Gonadotropin (PMSG)†	25001388	<ul style="list-style-type: none"> ○ ● SST, RTT 1 mL serum  	Horses only - samples should be taken 50-100 days post-breeding.	3-10 days	\$104.47
Cryptorchid Panel, Equine†	EQCRYPTO	<ul style="list-style-type: none"> ● ● RTT, LTT PREFERRED: 2 mL serum; 2 mL separated EDTA or heparinized plasma also accepted 	Testosterone, Estrone Sulfate	7-10 days	\$123.73
Gastrin†	GASTRIN	<ul style="list-style-type: none"> ● RTT 1 mL serum, frozen  		3-12 days	\$113.03
Reproductive Hormone Panel†	RHP	<ul style="list-style-type: none"> ● ● LTT, RTT (no SST) 2 mL serum or plasma (no SST) 	Estrogen, Progesterone, Testosterone	5-10 days	\$153.71

Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Aldosterone Baseline†	ALDOS	<ul style="list-style-type: none"> ● ● LTT or SST, RTT 0.5 mL EDTA plasma (preferred); 0.5 mL serum also accepted. Do not send whole blood.  		7-10 days	\$100.19
Insulin:Glucose Ratio, Equine	EQINSGLU	<ul style="list-style-type: none"> ● ● RTT, GRYTT 1 mL serum  		1-4 days	\$40.76
Estrogens, Total†	EQUEST	<ul style="list-style-type: none"> ● ● LTT, RTT (no SST) 2 mL serum or plasma 		3-5 days	\$93.77
Testosterone, Baseline†	TESTOS	<ul style="list-style-type: none"> ● ● RTT, LTT 2 mL serum or EDTA or heparinized plasma  		2-9 days	\$96.98
Oral Sugar Test (Insulin Testing), Equine†	EQGLUCOSE	<ul style="list-style-type: none"> ● ● RTT, LTT PREFERRED: 1 mL serum. 1 mL separated EDTA or heparinized plasma also accepted  		3-5 days	\$94.84
Insulin, Equine	5104	<ul style="list-style-type: none"> ● RTT 0.5 mL serum  		1-4 days	\$37.10
Testosterone Post-HCG/GnRH†	TESTPOST	<ul style="list-style-type: none"> ● ● RTT, LTT 1 mL serum or EDTA plasma  		3-8 days	\$89.49
Metanephrines, Plasma-free†	METAPLAS	<ul style="list-style-type: none"> ● ● LTT, GRNTT 1-2 mL EDTA or heparinized plasma, frozen  		5-7 days	\$216.78






Endocrinology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Metanephrines, Urine Fractionated†	METAUR	○● WTT, Urine 2.5-4 mL urine		5-7 days	\$115.86
Parathyroid Hormone-Related Protein (PTHrP)†	PTHRP	● LTT 0.5 mL separated EDTA plasma (frozen). Do not send whole blood. ❄️		5-7 days	\$126.94
Vitamin D Profile†	VITP	●● RTT, GRNTT 1.5 mL serum, separated ❄️💡	PTH, Free (Ionized) Calcium, and 25-hydroxyvitamin D Overnight fasting recommended for best results. Centrifuge samples and separate serum/plasma within an hour of collection. Protect sample from light. Numeric results will be complete within 3-9 days. An interpretation will be provided an additional 1-10 days post numeric result.	3-9 days	\$140.53



GI/Pancreatic

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Amylase (Mich. State Univ.)†	1004	○ SST 0.5 mL serum 📄		2-3 days	\$88.32

Fecal/Parasitology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Ova and Parasite (Centrifugal Fecal Flotation) Test	2320A	● Feces 3-5g fresh feces 		Next day by 5 PM (1-2 days for equine samples)	\$16.80
McMaster (Fecal Egg Counting) Test	2328	● Feces 5 gm fresh feces 	Quantitative analysis	1-2 days	\$25.48
Centrifugal Fecal Flotation Test (O&P) & McMaster (Fecal Egg Counting) Test	25000487	● Feces 8-10 gm fresh feces 		1-2 days	\$25.48
Baermann Test (Lungworm)	BAER	● Feces 10 gm fresh feces 		3-5 days	\$53.47
Fecal Occult Blood	7255	● Feces 1-3 gm fresh feces 		1-2 days	\$38.18
Gross Parasite Identification	25001435	Sterile leakproof container Whole parasite (adult worm, larvae, ectoparasite, etc.) in sealed, leakproof container.	Comment, Gross Parasite ID	1-5 days	\$34.35

Microbiology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Equine Uterine Culture	7154	Culture swab	Bacterial culture and identification; Sensitivity performed when appropriate; varies by host species, source, and organism	2-5 days	\$63.65
Aerobic Culture	7156	○ WTT/sterile tube Culture swab or fluid/tissue in WTT or sterile tube	Bacterial culture and identification; Sensitivity performed when appropriate; varies by host species, source, and organism	2-5 days	\$69.55
Anaerobic Culture	7157	○ WTT/sterile tube Culture swab or fluid/tissue in WTT, anaerobic transport medium, or sterile tube	Anaerobic bacterial culture and identification	2-7 days	\$38.18
Aerobic and Anaerobic Culture	7135	○ WTT/sterile tube Culture swab or fluid/tissue in WTT, anaerobic transport medium, or sterile tube	Bacterial culture and identification; Sensitivity varies by host species, site, and organism	2-7 days	\$107.00
Mycoplasma Culture†	MYCOPLASCU	Swab/wash (Amies with charcoal) or tissue in a sterile, leakproof container 		Up to 21 days	\$93.77
Urine Culture and Sensitivity	7140	○ ● WTT/sterile tube, Urine 1 mL urine in WTT or sterile container 	Bacterial culture and identification; Sensitivity performed when appropriate; varies by host species, source, and organism	1-5 days	\$69.55
Hold Urine Culture for 48 hours	25000534		Bacterial culture will be held for an additional 24 hours (48 hours total). This test is intended for use in working up complicated/chronic urinary tract infections and must be ordered at the time the urine culture is initially requested.	2 days	\$0.01
Reflex Aerobic Culture (Contingent)	25000536	○ WTT/sterile tube Culture swab, fluid in WTT or sterile tube, or fresh tissue in WTT or sterile tube	An aerobic bacterial culture will be automatically added, if indicated (bacteria seen); can be applied to urinalysis or cytology samples	2-5 days	\$59.99

Microbiology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Blood Culture	7150	Two blood culture bottles should be submitted. Blood should be drawn from two sites (i.e. left jugular, right saphenous)	Bacterial culture and identification; Sensitivity performed when appropriate; varies by host species, source, and organism	3-9 days	\$140.07
Dermatophyte Culture	7160	○ WTT/sterile tube Plucked hair, nail clippings, skin scrapings or brushings	Dermatophyte fungal culture and identification	12 days for final report	\$70.06
Fecal Culture	7155	● Feces 5 gm fresh feces, submitted in a fecal container	Bacterial culture and identification; Sensitivity performed when appropriate; varies by host species, source, and organism	2-7 days	\$70.06
Fungal Culture	7145	○ WTT/sterile tube Culture swab or fluid/tissue in WTT or sterile tube	Fungal culture and identification	Up to 26 days for final report	\$70.06
<i>Mycobacterium</i> Culture†	MYCO	Intestine, lymph node, rectal, or skin biopsy in sterile, leakproof container		60 days	\$104.47
Acid-Fast Stain, Non-Avian/Reptile	ACIDFAST	○ WTT/sterile tube Culture swab or fluid/tissue in WTT or sterile tube	For detection of acid-fast positive bacilli	1-4 days	\$43.28
Fungal Susceptibility 1 Drug†	FUNGALSUS1		Add-on susceptibility to fungal plate grown at Zoetis Reference Laboratories.	1-2 weeks	\$136.87
Fungal Susceptibility 2 Drugs†	FUNGALSUS2		Add-on susceptibility to fungal plate grown at Zoetis Reference Laboratories.	1-2 weeks	\$216.50







See the complete Directory of Services online at zoetisreflabs.com/DOS







Microbiology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Fungal Susceptibility 3 Drugs[†]	FUNGALSUS3		Add-on susceptibility to fungal plate grown at Zoetis Reference Laboratories.	1-2 weeks	\$296.13
Skin Culture and Sensitivity (Dermatology)	25001434	○ WTT, Culture swab Culture swab or fluid/tissue in WTT/sterile tube	Aerobic bacterial culture of SKIN specimens. Bacterial identification; antimicrobial susceptibility testing performed when appropriate; varies by organism, source, and host species.	2-5 days	\$63.66









Serology/Infectious Diseases

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Coggins AGID, Paper Submission	COGAGID	● ● SST, RTT 1 mL serum 		2-5 days	\$31.83
Coggins ELISA, Paper Submission	COGELISA	● ● SST, RTT 1 mL serum 		1-4 days	\$31.83
Coggins AGID, Electronic Submission	E1560	● ● SST, RTT 1 mL serum 		2-5 days	\$19.12
Coggins ELISA, Electronic Submission	E1561	● ● SST, RTT 1 mL serum 		1-4 days	\$19.12









Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Aspergillus Titer AGID[†]	ASP	<ul style="list-style-type: none"> ● SST 1 mL serum  	IgG	3-5 days	\$111.96
Blastomyces Quantitative Antigen EIA[†]	BLASTANT	<ul style="list-style-type: none"> ● ● ● ● SST, LTT, GRNTT, BTT 1.2 mL serum or plasma, separated; or 0.8 mL CSF; or 0.5 mL urine or BAL fluid  		2-5 days	\$145.13
Coccidioides Titer AGID[†]	6129	<ul style="list-style-type: none"> ● SST 1 mL serum  	IgG, IgM	3-5 days	\$107.68
Equine Viral Arteritis[†]	EVA	<ul style="list-style-type: none"> ● ● SST, RTT 1 mL serum  		5-9 days	\$142.70
Leptospira MAT, 5 Serovar[†]	6164	<ul style="list-style-type: none"> ● RTT 0.5 mL serum  	Detection of antibodies to <i>Leptospira interrogans</i> serovars <i>pomona</i> , <i>hardjo</i> , <i>icterohaemorrhagiae/copenhageni</i> , <i>grippityphosa</i> , and <i>canicola</i>	Up to 10 days	\$109.57
Leptospira MAT, 5 Serovar & PCR[†]	6163	<ul style="list-style-type: none"> ● ● RTT, LTT 0.5 mL serum in RTT for MAT; For PCR: 1. 5 mL urine; 2. 1 mL EDTA whole blood; 3. Urine and EDTA whole blood; 4. 1 gm tissue (eg, kidney and/or liver); 5. 1 mL ocular fluids; or 6. Bacterial culture  	MAT includes detection of antibodies to <i>Leptospira interrogans</i> serovars <i>pomona</i> , <i>hardjo</i> , <i>icterohaemorrhagiae/copenhageni</i> , <i>grippityphosa</i> , and <i>canicola</i>	Up to 10 days	\$147.40







Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Sarcocystis neurona</i> by Western Blot†	SNEURONA	●● SST, RTT 1 mL serum or 1 mL CSF (RTT) 		5-7 days	\$159.04
<i>Aspergillus</i> Galactomannan Antigen EIA†	ASPGUR	● RTT 0.8 mL serum (separated), CSF, urine, or BAL fluid 		2-5 days	\$153.71
<i>Cryptococcus</i> Antigen Titer by Latex Agglutination	7216	●● SST, RTT 0.5 mL serum or CSF 		1-2 days	\$95.54
Equine Protozoal Myelitis IFA Titers (SarcoFluor™ & NeoFluor™)†	EPMIFAT	● RTT 2 mL separated serum OR 2 mL CSF  Serum/Plasma  CSF	SarcoFluor detects <i>Sarcocystis neurona</i> antibodies; NeoFluor detects <i>Neospora hughesi</i> antibodies.	7-12 days	\$221.49
<i>Histoplasma</i> Quantitative Antigen EIA†	HISTOANT	●●●● SST, RTT, LTT, GRNTT, BTT 1.2 mL serum or plasma, separated; or 0.8 mL CSF; or 0.5 mL urine or BAL fluid 		2-5 days	\$145.13
Lyme Multiplex Test, Equine†	LYMEMPLXEQ	● RTT 2 mL serum; 1 mL CSF 		5-7 days	\$123.73
<i>Pythium</i> Serology†	OOPYTHIUM	● RTT 0.5 mL serum 		3-5 days	\$147.27




Serology/Infectious Diseases (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Streptococcus equi</i> ELISA†	STREPEQUI	● RTT 1 mL serum 		5-7 days	\$148.34
Equine Protozoal Myelitis IFA Titer (SarcoFluor™)†	SARCOIFAT	● RTT 2 mL separated serum OR 2 mL CSF ❄ Serum/Plasma  CSF	SarcoFluor detects <i>Sarcocystis neurona</i> antibodies.	7-12 days	\$146.59
West Nile Virus IgM ELISA†	WESTNILE	●● SST, RTT 2 mL serum 		3-5 days	\$86.28
Caseous Lymphadenitis (<i>Corynebacterium pseudotuberculosis</i>) Titer†	CALYMTTR	●● SST, RTT 2 mL serum 		3-5 days	\$84.14
Equine Influenza (HAI)†	EQINFLUENZA	● RTT 2 mL serum 		7-10 days	\$92.70
Equine Protozoal Myelitis IFA Titer (NeoFluor™)†	NEOIFAT	● RTT 2 mL separated serum OR 2 mL CSF 	NeoFluor detects <i>Neospora hughesi</i> antibodies.	7-12 days	\$147.27
Potomac Horse Fever (<i>Neorickettsia risticii</i>) IFA†	POTOMAC	○ SST 1 mL serum 	IgG	5-7 days	\$115.17
<i>Rhodococcus equi</i> ELISA†	RHODOCOCCUS	●● SST, RTT 1 mL serum 		3-5 days	\$68.48

Clinical Pharmacology/Toxicology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Drug Screen, Firocoxib†	DGSCRFIO	<p>● ● ● SST, RTT, LTT 10 mL serum or EDTA plasma</p> 		2-4 days	\$355.92
Drug Screen with Tranquilizers, NSAIDs†	DSNSAIDTQ	<p>● ● ● SST, RTT, LTT 10 mL serum or EDTA plasma</p> 	Multiple drugs tested, updated periodically. Call Customer Service at 1-888-965-9652 for more information.	7-10 days	\$216.82
Quantitative Copper (Blood or Tissue)†	6171	<p>● ● ● RTT, LTT, GRNTT 2 mL serum or whole blood (EDTA or heparin) or 50 mg wet tissue (fresh, frozen, fixed, or embedded in a paraffin block)</p>		7-10 days	\$194.61
Selenium and Trace Elements Screen†	SELETREL	<p>● ● LTT, RBTT 2 mL serum (RBTT) AND 1 mL EDTA whole blood (LTT)</p> <p> Whole blood  Serum/Plasma</p>	Se, Ca, Cu, Fe, Mg, Phos, Zn. This test requires both EDTA whole blood and serum. For serum, collect whole blood into a RBTT with clot activator (no anticoagulant) and centrifuge to separate serum. Aliquot serum into a new RBTT with clot activator. Hemolysis may affect calcium, copper, iron, magnesium, phosphorus, selenium, and zinc levels. Whole blood shipped chilled with an ice pack, serum shipped chilled or frozen with an ice pack.	7-10 days	\$134.81
Zinc, Serum/Plasma†	ZINC	<p>● ● RBTT, RTT 2 mL Na EDTA plasma or separated serum (no SST)</p> 		3-5 days	\$108.75
Anticoagulant (Rodenticide)†	ANTICOAG	<p>● LTT 3 mL EDTA whole blood</p> 	Includes brodifacoum, bromadiolone, chlorophacinone, difenacoum, difethialone, diphacinone, and warfarin	6-12 days	\$197.90

Clinical Pharmacology/Toxicology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Drug Screen, NSAIDs, Equine†	DSNSAIDS	<p>● ● ● SST, RTT, LTT 10 mL serum or EDTA plasma</p> 	Multiple drugs tested, updated periodically. Call Customer Service at 1-888-965-9652 for more information.	2-4 days	\$141.92
Drugs of Abuse†	DOA	<p>● ● SST, RTT 10 mL serum or urine (10 gm liver or kidney also accepted)</p>	Amphetamine, cocaine, benzoylecgonine, norcocaine, ephedrine, lysergic acid diethylamide, (LSD), methamphetamine, MDMA, nicotine, phentermine, psilocin, delta-9-THC, THC-OH, THC-COOH, methadone, morphine, HU-210, JWH-018	10-14 days	\$292.53
LC/Mass Spectrometry Toxicology Screen†	LCMS	<p>● RTT 20 mL fluid; 50 gm tissue; 50 gm feed; OR 50 gm source material</p> 	Toxicology screen for unknown exposure	2-3 weeks	\$248.33
Blood Lead†	6240	<p>● ● LTT, GRNTT 0.5 mL EDTA or heparinized whole blood</p>		5-7 days	\$114.10
Gabapentin†	GAB	<p>● ● LTT, RTT (no SST) 0.5 mL EDTA plasma or serum (no SST)</p>		3-5 days	\$142.99
Toxic Heavy Metal Panel†	25001395	<p>● ● ● LTT, GRNTT, BTT 5 mL EDTA, heparinized, or citrated whole blood; 3 gm tissue (liver or kidney)</p> 	Quantitative analysis of arsenic, cadmium, lead, thallium, mercury, selenium	6-12 days	\$218.96



See the complete Directory of Services online at zoetisreflabs.com/DOS

Pathology

Histopathology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Histopathology Simple, Single Site	ZN100	Formalin-fixed tissue	<ul style="list-style-type: none"> - Single site - No more than 10 cm in any dimension - Not for dermatologic disease or with other criteria of Complex Level 1 or 2 	3-7 days	\$90.95
Histopathology Simple, Multiple Sites (Up to 5 Sites)	ZN101	Formalin-fixed tissue	<ul style="list-style-type: none"> - Simple samples from up to 5 unique sites - No more than 10 cm in any dimension 	3-7 days	\$133.75
Histopathology Complex, Level 1	ZN102	Fresh or formalin-fixed tissue	<ul style="list-style-type: none"> - Skin biopsies for dermatologic disease (without dermatologist consult) - Multiple gastrointestinal biopsies - Bone marrow - Whole organs/major components of organs: eye, spleen, testicles, uterus/ovaries, mammary chains (2 or more connected glands), amputated digit/nail or tail, liver lobe, or lung lobe - Large samples (>10cm in any dimension) - Individual teeth or other small samples containing bone - Simple samples from 6-15 unique sites - Formalin-fixed whole bodies <300 g 	3-9 days	\$168.25
Histopathology Complex, Level 2	ZN103	Fresh or formalin-fixed tissue	<ul style="list-style-type: none"> - Mandible, maxilla (all jaws, regardless of size) - Amputated limb - Other large/atypical samples containing bone - Multiple Complex Level 1 samples - Simple samples from >15 sites. 	Varies by nature of sample	\$213.93

Pathology (cont'd)

Specialty Histopathology Tests






Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Dermatopathology Panel	ZN400	<ul style="list-style-type: none"> - Formalin-fixed skin biopsies - Detailed history and clinical images are required - Images can be sent to: referencelabpathology@zoetis.com 	<ul style="list-style-type: none"> - Biopsy report by dermatopathologist - Addendum with case management guidance from boarded dermatologist 	Dermatopathology: 3-7 days. Consult provided after dermatopathology.	\$184.67
Nephropathology (Renal) Panel	ZN25000541	<ul style="list-style-type: none"> - 1 or more 14G to 18G core needle biopsy, a wedge biopsy, or other larger renal sample (eg, nephrectomy) - Detailed history is required - Images can be sent to: referencelabpathology@zoetis.com 	<ul style="list-style-type: none"> - Biopsy report by kidney specializing pathologist - H&E, PAS, trichrome, Jones methenamine silver stains (Congo Red also included, if indicated) 	3-7 days	\$230.05

Pathology (cont'd)

Cytopathology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Cytopathology Simple/ Single Site	7120	<ul style="list-style-type: none"> LTT Unstained, air-dried slides; fluid from a mass in LTT <ul style="list-style-type: none"> Slides Liquid 	For the following sample types: <ul style="list-style-type: none"> Single lesion/tissue (FNA, swab, impression smear) Single or multiple lymph nodes Single joint (slides only) Fluid from mass lesion (cytology only) 	1-2 days	\$96.30
Cytopathology Complex	7201	<ul style="list-style-type: none"> LTT Unstained, air-dried slides and/or fluid in LTT <ul style="list-style-type: none"> Slides Liquid 	For the following sample types: <ul style="list-style-type: none"> Multiple lesions/tissues (FNA, swab, impression smear) Joint fluid in LTT (single or multiple) Multiple joints (slides only) Cavitary effusion (single or multiple) CSF Bone marrow BAL/TTW Other fluids/lavage 	1-2 days	\$130.52


Clinical Immunology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Complete Amino Acid Profile†	COMPAA	● GRNTT 1 mL separated heparinized plasma 		7-9 days	\$235.73
Allergy Panel, Equine†	6307	● ● SST, RTT 3 mL serum		10-14 days	\$361.27
Protein Electrophoresis	6345	● ● SST, RTT, sterile leakproof container 0.5 mL serum, >10 mL urine 		3-6 days	\$70.06
IgM, Equine	25000540	● ● SST, RTT 2 mL serum 		1-3 days. Preliminary reported at 4 hours after setup. Final results at 24 hours after setup.	\$39.56
IgG, Equine	FOALIGG	● ● SST, RTT 2 mL serum 		1-3 days. Preliminary reported at 4 hours after setup. Final results at 24 hours after setup.	\$39.56
IgG, Equine: Rapid, Semi-Quantitative	FOALIGGRAP	● ● ● SST, RTT, LTT 1 mL fresh whole blood or serum 		1-2 days	\$49.46

Nutrition

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Selenium†	SELENIUM	<ul style="list-style-type: none"> ●● LTT, RTT (no SST) 0.5 mL EDTA whole blood preferred; 0.5 mL serum or EDTA plasma also accepted ❄️ Serum/Plasma 📄 Whole blood 		3-5 days	\$100.19
Trace Element Screen†	TMINPAN	<ul style="list-style-type: none"> ● RBTT 2 mL serum from RBTT ❄️ 	Ca, Cu, Fe (unbound), Mg, Na, K, Phos, Zn	5-7 days	\$88.16
Thiamine (Vitamin B1)†	THIAMINE	<ul style="list-style-type: none"> ●● GRNTT, LTT 3 mL EDTA or heparinized whole blood ❄️💡 		5-7 days	\$122.16
Vitamin E†	2331	<ul style="list-style-type: none"> ●●●● SST, RTT, LTT, GRNTT 1.5 mL serum or heparinized or EDTA plasma, separated ❄️💡 		3-5 days	\$100.19
Vitamin E and Mineral Panel†	VITEMINE	<ul style="list-style-type: none"> ●●● RTT, LTT, RBTT 1 mL EDTA whole blood, 1 mL RBTT serum (separated), 1 mL serum no-additive RTT (separated and protected from light) ❄️💡 	Vit E, Ca, Cu, Fe, Mg, Phos, Zn, Se This test requires EDTA whole blood, serum from a royal blue top tube (RBTT) and serum from a no-additive RTT. For RBTT serum, collect blood into RBTT with clot activator (no anticoagulant), allow clot to form, centrifuge and separate serum. Separated serum to be placed in new RBTT for shipment to lab. Label tube serum for no-additive RTT, collect blood into no-additive RTT, allow clot to form, centrifuge and separate serum. Separated serum to be placed in new no-additive RTT and wrap with foil to protect from light. Samples should be kept refrigerated and shipped chilled on an ice pack.	8-10 days	\$144.89

Nutrition (cont'd)







Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Vitamin A†	VITA	●●●● SST, RTT, LTT, GRNTT 1.5 mL serum or EDTA or heparinized plasma, separated 		3-5 days	\$100.19






See the complete Directory of Services online at zoetisreflabs.com/DOS

Add-Ons

You can quickly and easily add any of the tests below to an existing order, provided sufficient sample remains at the lab and is within stability guidelines. If you order electronically, you can place add-on orders through your PIMS or the ZoetisDx portal. If you use paper requisition forms, call customer service at **1-888-965-9652**.

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Anaerobic ID Add-on	A7157	○ WTT/sterile tube Culture swab or fluid/tissue in WTT, anaerobic transport medium, or sterile tube.		2-7 days	\$32.45
Cortisol Add-on	A5101	● RTT 0.5 mL serum 		8 AM next day (1-2 days for feline/equine samples)	\$13.67
Free (Ionized) Calcium (Mich. State Univ.) Add-on†	A6304	● RTT 0.5 mL serum, separated Ensure that the sample is treated anaerobically. 		3-6 days	\$70.62
Free T4 Add-on	A5103	●● SST, RTT 0.5 mL serum 		8 AM next day	\$16.83
Reticulocyte Count Add-on	25000525	● LTT 0.5 mL EDTA whole blood 		8 AM next day	\$15.76
Total T4 Add-on	A5102	●● SST, RTT 0.5 mL serum 		8 AM next day	\$14.41
Triglycerides Add-on	A1023	●● SST, RTT 0.5 mL serum 		8 AM next day	\$7.59

Add-Ons (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Urine Protein:Creatinine Ratio Add-on	A1555	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 1 mL urine in sterile container 		8 AM next day	\$14.41
Urinalysis Add-on	25001700	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 5 mL urine in a sterile container 		8 AM next day	\$14.41
Urinalysis with Reflex UPC Add-on	25001701	<input type="radio"/> <input checked="" type="radio"/> WTT, Urine 5 mL urine in sterile container 			\$19.77



See the complete Directory of Services online at zoetisreflabs.com/DOS



Avian, Exotic & Small Mammal

|

|





|

|








|

|





Bundles

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Well Bird Panel Package†	25001468	<p>●●● RTT, LTT, GRNTT, slides 0.1 mL EDTA whole blood, slides; 0.3 mL plasma or serum</p> 	CBC, Basic Avian Chem Panel, Bile acids, EPH	4-5 days	\$160.29
Advanced Reptile Panel†	25001448	<p>●● RTT, GRNTT, slides 0.1 mL heparinized whole blood, slides, 0.3 mL plasma or serum</p> 	CBC, Advanced Reptile Chem Panel, EPH	4-5 days	\$167.16
Basic Reptile Panel†	25001454	<p>●● RTT, GRNTT, slides 0.1 mL heparinized whole blood, slides, 0.25 mL serum or plasma</p> 	CBC, Basic Reptile Chem Panel, EPH	4-5 days	\$156.85
Basic Mammal Panel†	25001452	<p>●●● RTT, LTT, GRNTT, slides 0.1 mL EDTA whole blood, slides, 0.2 mL serum or plasma</p> 	CBC, Basic Mammal Chem Panel, EPH	4-5 days	\$210.35

Hematology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Fibrinogen Heat Precipitation Method†	25001457	<ul style="list-style-type: none"> LTT 0.1 mL EDTA whole blood 		3-5 days	\$80.14
Mini Avian Panel †	25001459	<ul style="list-style-type: none"> LTT, GRNTT 0.1 mL EDTA whole blood, slides, 0.1 mL heparinized plasma 	CBC, Electrophoresis	4-5 days	\$120.21
Elasmobranch CBC and Diff†	25001455	<ul style="list-style-type: none"> GRNTT, slides 0.1 mL heparinized whole blood, slides 		3-5 days	\$101.90
Mini Elasmobranch Panel†	25001460	<ul style="list-style-type: none"> GRNTT, 2 slides 0.1 mL heparinized whole blood, 2 slides, 0.1 mL heparinized plasma 	CBC, Electrophoresis	4-5 days	\$120.21
CBC Hemocytometer Cnt Man Diff Reptile†	25001444	<ul style="list-style-type: none"> GRNTT, slides 0.1 mL heparinized whole blood, slides 	Automated count and manual diff.	3-5 days	\$101.90
Reptile CBC and Diff†	25001447	<ul style="list-style-type: none"> GRNTT, slides 0.1 mL heparinized whole blood, slides 	By gold standard hemocytometer/Natt & Herrick's method	3-5 days	\$101.90
Mini Reptile Panel†	25001462	<ul style="list-style-type: none"> GRNTT, 2 slides 0.1 mL heparinized whole blood, 2 slides, 0.1 mL heparinized plasma 		4-5 days	\$120.21









Hematology (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Mini Mammal Panel†	25001461	<ul style="list-style-type: none"> LTT, GRNTT, slides 0.1 mL EDTA whole blood, slides, 0.1 mL heparinized plasma 	CBC, Electrophoresis	4-5 days	\$125.94
Mini Teleost Panel†	25001463	<ul style="list-style-type: none"> GRNTT, 2 slides 0.1 mL heparinized whole blood, 2 slides, 0.1 mL heparinized plasma 	CBC, Electrophoresis	4-5 days	\$120.21
Teleost CBC and Diff†	25001465	<ul style="list-style-type: none"> GRNTT, slides 0.1 mL heparinized whole blood, slides 		3-5 days	\$101.90
CBC Automated Count Manual Diff Mammal†	25001443	<ul style="list-style-type: none"> LTT, slides 0.1 mL EDTA whole blood, slides 	Automated count and manual diff.	3-5 days	\$101.90

Chemistry

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Lipoprotein Chemistry Panel†	25001458	<ul style="list-style-type: none"> GRNTT 0.1 mL heparinized plasma 		4-5 days	\$105.33
Adv Avian Chem Panel†	25001437	<ul style="list-style-type: none"> RTT, GTT 0.3 mL heparinized plasma 	Amylase, AST, BUN, Creat, Ca, Chol, Cl, CO2, CPK, GGT, Glu, Lipase, Mg, Phos, Potassium, Sodium, TP, Trig, Uric Acid, GLDH	4-5 days	\$139.68

Chemistry (cont'd)

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Adv Mammal Chem Panel†	25001438	● ● RTT, GRNTT 0.3 mL serum or heparinized plasma 	ALT, AlkPhos, BUN, Crea, Glu, Ca, Na, Phos, Potassium, Cl, Mg, TP, Alb, AST, Amylase, CPK, TBili, GGT	4-5 days	\$139.68
Adv Reptile Chem Panel†	25001439	● GRNTT 0.3 mL heparinized plasma 	AlkPhos, Amy, Anion Gap, AST, BUN, CPK, Ca, Chol, Cl, CO2, GGT, Glu, K, Lipa, Mg, Na, Osmolality, Phos, TP, Trig, Uric Acid, GLDH	4-5 days	\$139.68
Bile Acids †	25001442	● GRNTT 0.1 mL heparinized plasma 		3-5 days	\$91.59
Basic Avian Chemistry Panel†	25001450	● GRNTT 0.2 mL heparinized plasma 	AST, Ca, Glu, TP, Uric Acid, CPK, Phos	4-5 days	\$104.19
Elasmobranch Chemistry Panel†	25001456	● GRNTT 0.3 mL heparinized plasma 	Anion gap, AST, BUN, Ca, Chol, Cl, CO2, CPK, GGT, Glu, K, Mg, Na, Osmolality, Phos, TP, Trig, Uric Acid	4-5 days	\$133.95
Basic Reptile Chemistry Panel†	25001453	● GRNTT 0.25 mL heparinized plasma 	AlkPhos, AST, CPK, Ca, Cl, GGT, Glu, K, Na, Phos, TP, Trig, Uric Acid	4-5 days	\$123.65
Basic Mammal Chemistry Panel†	25001451	● ● RTT, GRNTT 0.2 mL serum or heparinized plasma 	ALT, BUN, Crea, Glu, Phos, Ca, TP	4-5 days	\$104.19
Teleost Chemistry Panel†	25001466	● GRNTT 0.3 mL heparinized plasma 	Anion gap, AST, BUN, Ca, Chol, Cl, CO2, CPK, GGT, Glu, K, Mg, Na, Osmolality, Phos, TP, Trig, Uric Acid	3-5 days	\$133.95



Endocrinology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
T4, Total (Mammals Only)†	25001464	<ul style="list-style-type: none"> GRNTT 0.1 mL heparinized plasma 		3-5 days	\$99.61

Microbiology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Sex Determination (Psittacine)†	25001539	<ul style="list-style-type: none"> GRNTT 0.1 mL whole heparinized whole blood 		7-10 days	\$80.14


Serology/Infectious Diseases

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Aspergillus Panel PLUS †	25001440	<ul style="list-style-type: none"> GRNTT 0.25 mL heparinized plasma 	Antibody (both assays), EPH + Galactomannan, Gliotoxin	3-5 days	\$217.53
Aspergillus Panel – PenguinPLUS †	25001449	<ul style="list-style-type: none"> GRNTT 0.25 mL heparinized plasma 	Hydroxybutyrate, EPH, Gliotoxin	4-5 days	\$166.01

Clinical Immunology

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Protein Electrophoresis (EPH) †	25001446	<ul style="list-style-type: none"> ● GRNTT 0.1 mL heparinized plasma 		3-5 days	\$95.03

Nutrition

Test Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Vitamin D (25-Hydroxy) Mammals Only†	25001467	<ul style="list-style-type: none"> ● ● RTT, GRNTT 0.1 mL heparinized plasma or serum 		3-5 days	\$99.61



See the complete Directory of Services online at zoetisreflabs.com/DOS



PCR

Comprehensive PCR test menu
provided by the University of Georgia
College of Veterinary Medicine.

Tests

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Adult Equine Diarrhea PCR Panel†	25001469	● Feces Minimum 2 g feces in a clean leakproof container	<i>Salmonella</i> spp, <i>Clostridium difficile</i> enterotoxin A, <i>Clostridium difficile</i> cytotoxin B, <i>Clostridium perfringens</i> type A, <i>Clostridium perfringens</i> type C, Equine coronavirus	5-7 days	\$212.93
Anaplasma/Ehrlichia PCR†	25001470	● LTT EDTA whole blood (minimum 0.5 mL)	All species of <i>Anaplasma</i> and <i>Ehrlichia</i>	3-5 days	\$100.75
Anemia PCR Panel, Canine†	25001471	● LTT EDTA whole blood (minimum 0.5 mL)	<i>Anaplasma/Ehrlichia</i> combo, <i>Babesia</i> spp, <i>Hepatozoon</i> spp, <i>Leptospira</i> spp, <i>Rickettsia</i> spp	3-5 days	\$187.25
Anemia PCR Panel, Feline†	25001472	● LTT EDTA whole blood (minimum 0.5 mL)	<i>Anaplasma</i> , <i>Bartonella</i> , <i>Cytauxzoon felis</i> , <i>Ehrlichia</i> , feline hemotropic <i>Mycoplasma</i> , FeLV, FIV	3-5 days	\$181.90
Babesia spp (genus) PCR†	25001473	● LTT EDTA whole blood (minimum 0.5 mL)		3-5 days	\$111.82
Bartonella spp (genus) PCR†	25001474	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), synovial fluid (minimum 0.25 mL)		3-5 days	\$111.82
Blood Bank Disease PCR Panel, Canine†	25001475	● LTT EDTA whole blood (minimum 0.5 mL)	<i>Anaplasma/Ehrlichia</i> combo, <i>Babesia</i> spp, <i>Hepatozoon</i> spp, <i>Rickettsia</i> spp, <i>Bartonella</i> spp, <i>Leishmania</i> spp, <i>Borrelia burgdorferi</i>	7-12 days	\$214.00
Bordetella bronchiseptica PCR†	25001476	● RTT Nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage, respiratory tissue		3-5 days	\$100.75

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Borrelia burgdorferi</i> (Lyme) PCR†	25001477	<ul style="list-style-type: none"> ● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), synovial fluid (minimum 0.25 mL)		3-5 days	\$111.82
<i>Brucella spp</i> (genus) PCR†	25001478	<ul style="list-style-type: none"> ● ● RTT, LTT, YTT Fresh affected tissue, body fluids, swabs, EDTA whole blood (0.5 mL minimum), urine		3-5 days	\$118.74
<i>Campylobacteri coli</i> PCR†	25001480	<ul style="list-style-type: none"> ● Feces 1 g feces, fecal swab (two required) or GI tissue	Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	3-5 days	\$100.75
<i>Campylobacteri jejuni</i> PCR†	25001479	<ul style="list-style-type: none"> ● Feces 1 g feces, fecal swab (two required) or GI tissue.	Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	3-5 days	\$100.75
Canine Adenovirus PCR†	25001481	<ul style="list-style-type: none"> ● RTT, Culture swab Nasopharyngeal or conjunctival swab, transtracheal wash, feces, fresh affected tissue		3-5 days	\$100.75
Canine Puppy Diarrhea PCR w O&P, Giardia	25001623	O&P: 3g fresh feces Giardia: 5g fresh feces Diarrhea Panel: <ul style="list-style-type: none"> ● Feces 1 g feces, fecal swab (two required) or GI tissue.	Diarrhea Panel: Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	O&P & <i>Giardia</i> : Next day by 5 PM Diarrhea Panel: 5-7 days	\$107.91



See the complete Directory of Services online at zoetisreflabs.com/DOS



Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Adult Canine Diarrhea PCR Panel†	25001618	Sterile leakproof container PREFERRED: Submit at least 1 g feces in a clean leak-proof container. Accepted: Fecal swab, please submit two.	<i>Campylobacter jejuni/coli</i> , <i>Salmonella spp.</i> , <i>Lawsonia intracellularis</i> , <i>Clostridium difficile</i> enterotoxin A, <i>Clostridium difficile</i> cytotoxin B, <i>Clostridium perfringens</i> enterotoxin, <i>Giardia intestinalis</i>	5-7 days	\$198.49
Canine Distemper Virus (CDV) PCR†	25001482	● ● RTT, LTT, Culture swab 0.5 mL EDTA whole blood (preferred), fresh tissue, CSF, nasopharyngeal swab, bronchoalveolar lavage or transtracheal wash		3-5 days	\$111.82
Canine Herpesvirus 1 (CHV-1) PCR†	25001483	● RTT, Culture swab Fresh affected tissue (e.g. liver), nasopharyngeal or corneal swab		3-5 days	\$100.75
Canine Influenza PCR†	25001484	● RTT, Culture swab Nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage, respiratory tissue	Influenza A subtypes	3-5 days	\$100.75
Canine Mycoplasma PCR†	25001485	● RTT, Culture swab Nasopharyngeal swabs, nasal swab, transtracheal wash, bronchiolar lavage, respiratory tissue	<i>Mycoplasma canis</i> and <i>Mycoplasma cynos</i>	3-5 days	\$100.75
Canine Neurological Disease, Central Nervous System PCR Panel†	25001486	○ WTT/Sterile Container CSF- 500uL minimum or fresh tissue	<i>Toxoplasma gondii</i> , <i>Cryptococcus neoformans</i> , <i>Neospora caninum</i> , <i>Borrelia burgdorferi</i> (Lyme disease), Canine distemper virus (CDV)	3-5 days	\$168.88

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Canine Parvovirus (CPV) PCR†	25001487	● Feces 1 g feces, fecal swab (two required) or GI tissue	<p>Detects all CPV-2 subtypes. Additional testing is available to distinguish subtypes. Please call the laboratory for more information on CPV subtyping. This PCR is available as an individual test or as part of the Canine Diarrhea Panel. For parvovirus testing on cats, please see Feline panleukopenia PCR.</p> <p>Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.</p>	3-5 days	\$111.82
Canine Puppy Diarrhea PCR Panel†	25001488	● Feces 1 g feces, fecal swab (two required) or GI tissue	<p><i>Campylobacter jejuni/coli</i>, <i>Salmonella</i> spp, <i>Lawsonia intracellularis</i>, <i>Clostridium difficile</i> toxins A & B, <i>Clostridium perfringens</i> enterotoxin, <i>Giardia intestinalis</i>, Rotavirus, <i>Cryptosporidium</i> spp, Canine parvovirus, and coronavirus</p> <p>Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.</p>	5-7 days	\$204.37
Canine Reproductive/ Abortion PCR Panel†	25001489	● RTT Placenta, fetal tissue	<i>Brucella</i> spp, <i>Neospora caninum</i> , <i>Toxoplasma gondii</i> , Canine herpesvirus, <i>Leptospira</i> spp, <i>Mycoplasma</i> spp, <i>Ureaplasma</i> spp	7-14 days	\$156.85
Canine Respiratory PCR Panel†	25001490	● RTT, Culture swab Affected tissue, deep nasal swab, needle aspirates, bronchoalveolar lavage, pleural fluid, CSF	<i>Bordetella bronchiseptica</i> , <i>Mycoplasma canis</i> , <i>Mycoplasma cynos</i> , Canine adenovirus, Influenza A (all subtypes), Canine distemper virus, coronavirus, Canine parainfluenza virus, <i>Streptococcus zooepidemicus</i>	5-7 days	\$198.49
Chlamydia spp (genus) PCR†	25001491	Culture swab For CATS: corneal scraping, eye/corneal swab. For AVIAN species: choanal/cloacal swab. For CATTLE: vaginal scraping		3-5 days	\$100.75
Chlamydophila felis PCR†	25001492	● RTT, Culture swab Corneal scraping, corneal or conjunctival swab, fresh affected tissue		3-5 days	\$100.75

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Comprehensive (PCR and Serology) Vector Panel, Canine†	K9FCOMP	<ul style="list-style-type: none"> ● ● RTT, LTT 2 mL serum in RTT AND 2 mL whole blood in LTT 	<i>Babesia canis</i> , <i>B gibsoni</i> , <i>E canis</i> , <i>R rickettsia</i> (RMSF), <i>Bartonella henselae</i> , <i>B vinsonii</i> , <i>B koehlerae</i> (by IFA); Lyme, <i>Anaplasma</i> , <i>Ehrlichia</i> , and heartworm (by rapid assay); PCR for <i>Babesia</i> , <i>Bartonella</i> , <i>Anaplasma</i> , <i>Ehrlichia</i> , <i>Rickettsia</i> , and hemotropic <i>Mycoplasma</i> organisms	10 days	\$297.33
Comprehensive (PCR and Serology) Vector Panel, Feline†	FELVECTPNL	<ul style="list-style-type: none"> ● ● RTT, LTT 2 mL serum AND mL EDTA whole blood 	<i>Bartonella henselae</i> , <i>B vinsonii</i> , <i>B koehlerae</i> (by IFA); <i>Borrelia</i> (Lyme), <i>Anaplasma</i> , <i>Ehrlichia</i> , and heartworm (by rapid assay); PCR for <i>Babesia</i> , <i>Bartonella</i> , <i>Cytauxzoon</i> , <i>Anaplasma</i> , <i>Ehrlichia</i> , <i>Rickettsia</i> , and hemotropic <i>Mycoplasma</i> organisms	7-10 days	\$297.33
Coronavirus – Pan/Consensus PCR†	25001493	<ul style="list-style-type: none"> ● Feces Feces, fecal swab, GI tissue, effusion, EDTA whole blood (0.5 mL minimum), or airway specimens	Universal coronavirus detection	3-5 days	\$101.65
<i>Cryptococcus neoformans</i> PCR†	25001494	Culture swab Swab of affected tissue, needle aspirates, airway wash fluid, effusion, or CSF	Swabs of affected tissue, deep nasal passages, or needle aspirates applied to swabs should be submitted in a sterile, no additive tube.	3-5 days	\$100.75
<i>Cryptosporidium</i> spp (genus) PCR†	25001495	<ul style="list-style-type: none"> ● Feces 1 g feces, fecal swab (two required) or GI tissue.	Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	3-5 days	\$111.82
<i>Cytauxzoon felis</i> PCR †	25001496	<ul style="list-style-type: none"> ● LTT EDTA whole blood (minimum 0.5 mL)		3-5 days	\$100.75
<i>Ehrlichia</i> spp (genus) PCR†	25001498	<ul style="list-style-type: none"> ● LTT EDTA whole blood (minimum 0.5 mL), spleen		3-5 days	\$111.82
<i>Encephalitozoon cuniculi</i> PCR†	25001499	<ul style="list-style-type: none"> ● RTT, YTT Kidney (fresh or frozen), urine, phacoemulsified lens material		3-5 days	\$100.75

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Equine Abortion/ Reproduction PCR Panel†	25001500	● RTT, Culture swab Preputial wash, vaginal swab, vaginal wash, vaginal tissue, fetal tissue, abomasum, placenta	<i>Bartonella</i> spp, <i>Candida albicans</i> , Equine arteritis virus, EHV1, EHV4, <i>Leptospira</i> spp, <i>Salmonella enterica</i> , <i>Streptococcus equi</i> subsp <i>zooepidemicus</i>	3-5 days	\$166.01
Equine Foal Diarrhea PCR Panel†	25001501	● Feces Feces, GI tissue. Minimum 2 g feces in a clean leakproof container	<i>Salmonella</i> spp, <i>Clostridium difficile</i> enterotoxin A, <i>Clostridium difficile</i> cytotoxin B, <i>Clostridium perfringens</i> type A, <i>Clostridium perfringens</i> Type C, rotavirus, coronavirus	5-7 days	\$194.07
Equine Herpes Virus type 1 PCR†	25001502	● ● RTT, LTT, Culture swab EDTA whole blood, nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage, respiratory tissue		3-5 days	\$100.75
Equine Herpes Virus type 4 PCR†	25001503	● ● RTT, LTT, Culture swab EDTA whole blood, nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage, respiratory tissue		3-5 days	\$100.75
Equine Influenza PCR†	25001504	● RTT, Culture swab Nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage, respiratory tissue	Influenza type A	3-5 days	\$101.65
Equine Respiratory PCR Panel I†	25001505	● RTT, Culture swab Nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage, respiratory tissue	EHV 1 (Equine herpesvirus Type I), Influenza Type A, and <i>Streptococcus equi</i> subsp <i>equi</i>	5-7 days	\$184.04


Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Equine Respiratory PCR Panel II†	25001506	● RTT, Culture swab Nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage, respiratory tissue	EHV 1, EHV 4, Influenza Type A, <i>Streptococcus equi</i> subsp <i>equi</i> , Equine arteritis virus	5-7 days	\$184.04
Equine Tickborne Disease PCR Panel†	25001507	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), tissue	<i>Anaplasma</i> spp, <i>Ehrlichia</i> spp, <i>Borrelia burgdorferi</i> , <i>Rickettsia</i> spp, <i>Bartonella</i> spp, <i>Hepatozoon</i> spp If Equine Piroplasmiasis (<i>Babesia caballi</i> , <i>Theileria equi</i>) is suspected specimens must be submitted to USDA NVSL for testing.	3-5 days	\$163.71
Eubacteria Universal DNA – 16S rRNA PCR†	25001508	● RTT, Culture swab Affected tissue or body fluids, bacterial isolate	Universal Bacterial PCR	3-5 days	\$100.75
Feline Calicivirus PCR†	25001509	Culture swab Conjunctival swab, oropharyngeal swab, nasopharyngeal swab		3-5 days	\$100.75
Feline Coronavirus PCR†	25001510	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL)	Test for FIP diagnosis, pan-coronavirus PCR	3-5 days	\$101.65
Feline Diarrhea PCR Panel†	25001511	● Feces 1 g feces, fecal swab (two required) or GI tissue	<i>Campylobacter jejuni/ coli</i> , <i>Salmonella</i> spp, <i>Clostridium difficile</i> enterotoxin A, <i>Clostridium difficile</i> cytotoxin B, <i>Clostridium perfringens</i> enterotoxin, Feline parvovirus (Panleukopenia), <i>Trichostrongylus axei</i> , coronavirus Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	5-7 days	\$181.37



See the complete Directory of Services online at zoetisreflabs.com/DOS

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Feline Diarrhea PCR† with O&P, <i>Giardia</i>	25001626	O&P: 3g fresh feces <i>Giardia</i> : 5g fresh feces Diarrhea Panel: ● Feces 1 g feces, fecal swab (two required) or GI tissue 	Feline PCR Diarrhea Panel, Ova and Parasite, and <i>Giardia</i> antigen ELISA Diarrhea Panel: Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	O&P & <i>Giardia</i> : Next day by 5 PM Diarrhea Panel: 5-7 days	\$117.76
Feline Hemotropic (blood) <i>Mycoplasma</i> PCR†	25001512	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), spleen	<i>Mycoplasma haemofelis</i> and <i>Mycoplasma haemominutum</i>	3-5 days	\$101.65
Feline Herpes and <i>Chlamydomphila felis</i> PCR Combo†	25001513	Culture swab Corneal or conjunctival scraping, conjunctival swab, oronasal swab		3-5 days	\$105.93
Feline Herpesvirus PCR†	25001514	Culture swab Corneal scraping, oropharyngeal, nasopharyngeal, or conjunctival swab		3-5 days	\$111.82
Feline Immunodeficiency Virus FIV PCR†	25001515	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), bone marrow		3-5 days	\$111.82
Feline Leukemia Virus (FeLV) PCR†	25001516	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), bone marrow		3-5 days	\$111.82
Feline Neurological Disease PCR Panel†	25001517	● RTT CSF (0.5 mL minimum) or fresh tissue	<i>Toxoplasma gondii</i> , <i>Cryptococcus neoformans</i> , Feline infectious peritonitis virus (FIP), Feline panleukopenia virus (FPV), Feline immunodeficiency virus (FIV)	3-5 days	\$192.60

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
Feline Ocular PCR Panel†	25001518	Culture swab Corneal or conjunctival scraping, conjunctival swab, oronasal swab	Feline herpesvirus, <i>Chlamydophila felis</i> , Feline calicivirus	3-5 days	\$126.52
Feline Panleukopenia (Parvovirus) PCR†	25001519	Feces, fecal swab, GI tissue, or EDTA whole blood (0.5 mL minimum)		3-5 days	\$111.82
Feline Respiratory PCR Panel 1†	25001520	● RTT, Culture swab Nasal swab, oropharyngeal swab, bronchoalveolar lavage, transtracheal wash, fresh respiratory tissue	Feline herpesvirus, <i>Chlamydophila felis</i> , <i>Bordetella bronchiseptica</i> , <i>Mycoplasma</i> spp, Influenza A, Feline calicivirus	5-7 days	\$191.53
Fever of Unknown Origin PCR Panel, Canine†	25001521	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), affected tissues	<i>Anaplasma/Ehrlichia</i> combo, <i>Babesia</i> spp, <i>Bartonella</i> spp, <i>Brucella</i> spp, <i>Hepatozoon</i> spp, <i>Leptospira</i> spp, <i>Neospora caninum</i> , <i>Rickettsia</i> spp, <i>Toxoplasma gondii</i> , <i>Trypanosoma cruzi</i> , <i>Leishmania</i> spp	3-5 days	\$204.37
Fever of Unknown Origin PCR Panel, Feline†	25001522	● ● RTT, LTT, Fecal container, culture swab UGA requires all of the following three samples for this panel: 1. EDTA Whole Blood (minimum 0.5 mL) or affected tissues 2. Oral/Nasal/Ocular Swab 3. Feces or Fecal Swab	<i>Anaplasma/Ehrlichia</i> combo, <i>Bartonella</i> spp, <i>Cryptococcus neoformans</i> , <i>Cytauxzoon felis</i> , Feline calicivirus, Feline pan-coronavirus, <i>Mycoplasma haemofelis</i> , <i>Mycoplasma haemominutum</i> , Feline panleukopenia virus, FeLV, FIV, <i>Salmonella enterica</i> , <i>Toxoplasma gondii</i>	4-7 days	\$240.75
<i>Francisella tularensis</i> (Tularemia) PCR†	25001523	● ● RTT, LTT Fresh tissue (liver, spleen, lymph node), EDTA whole blood (0.5 mL minimum), respiratory secretions		3-5 days	\$101.65
Fungal Sequencing and PCR on Tissue Block†	FUNGSEQPCR	One paraffin block and two unstained slides		21 days	\$374.55

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Giardia intestinalis</i> PCR†	25001524	● Feces 1 g feces, fecal swab (two required) or GI tissue	Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	3-5 days	\$100.75
<i>Helicobacter</i> spp (genus) PCR†	25001525	● Feces 1 g feces, fecal swab (two required) or GI tissue	Submit feces in a clean, leakproof container. Use sterile swab for fecal swab. Submit fecal swab or GI tissue in sterile, no additive tube.	3-5 days	\$111.82
<i>Leishmania</i> spp (genus) PCR†	25001526	● ● RTT, LTT Preferred samples include skin lesion, bone marrow, lymph node aspirates, conjunctival swabs, spleen (fresh tissue), and EDTA whole blood (0.5 mL minimum).		3-5 days	\$100.75
<i>Leptospira</i> MAT, 5 Serovar & PCR†	6163	● ● RTT, LTT For MAT 0.5 mL serum in RTT; For PCR: 5 mL Urine; OR 1 mL EDTA whole blood (LTT); OR 1 gm tissue (e.g., kidney and/or liver); OR 1 mL aqueous humor; OR 0.5 mL CSF	MAT includes detection of antibodies to <i>Leptospira interrogans</i> serovars <i>pomona</i> , <i>hardjo</i> , <i>icterohaemorrhagiae/copenhageni</i> , <i>grippityphosa</i> , and <i>canicola</i>	Up to 10 days	\$147.40
<i>Leptospira</i> spp (genus) PCR†	25001527	● ● RTT, LTT, YTT Urine, EDTA whole blood (minimum 0.5 mL) Fresh kidney or liver, placenta, fetal tissue. Keep cold and ship overnight on ice. For increased diagnostic sensitivity please send both Blood and Urine.		3-5 days	\$101.65
<i>Listeria</i> spp (genus) PCR†	25001528	● RTT, Culture swab Brain (fresh), CNS swab, CSF (minimum 0.5 mL)		3-5 days	\$100.75

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Mycobacterium</i> spp (genus) PCR†	25001529	● RTT, Culture swab Affected tissue preferred, body fluids, swabs, feces, bacterial isolate		3-5 days	\$111.63
<i>Mycoplasma</i> spp (genus) PCR†	25001530	● LTT, Culture swab Affected tissue, body fluids, or swabs from affected sites		3-5 days	\$100.75
<i>Neospora caninum</i> PCR†	25001531	● ● ● RTT, LTT, Feces Fresh tissue (brain, fetal tissue), CSF (minimum 0.5 mL), feces		3-5 days	\$100.75
<i>Nocardia</i> spp (genus) PCR†	25001532	● RTT, Culture swab Fresh tissue, swabs from affected sites		3-5 days	\$100.75
Pan-Fungal PCR†	25001533	● RTT, Culture swab Affected fresh tissue, fluids, swab	Detects presence of fungal or yeast organism but does not speciate. Further identification is achieved through Sanger sequencing of the PCR amplicon.	3-5 days	\$101.65
PCR Formalin-Fixed Paraffin Embed Tissue†	25001534	● RTT Formalin fixed paraffin embedded tissue. Please submit 2-3 sliced scrolls in 1.5 mL micro-centrifuge tube.		5-7 days	\$88.74
<i>Rhodococcus equi</i> PCR†	25001535	● ● ● RTT, LTT, Feces, Culture swab Fresh tissue, nasal or oropharyngeal swab, transtracheal wash, bronchiolar lavage, feces, EDTA whole blood (minimum 0.5 mL)		3-5 days	\$111.82

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Rickettsia</i> spp (genus) PCR†	25001536	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), fresh tissue		3-5 days	\$100.75
Ringworm, PCR, Canine/Feline†	25001616	Sterile leakproof container Sterile collection tube Lesion biopsy, lesion swab, or plucked hairs that include roots and attached epithelial cells	<i>Microsporum canis</i> , <i>Trichophyton mentagrophytes</i> (<i>Dermatophytosis</i>)	3-5 days	\$153.42
Ringworm, PCR, Equine†	25001617	Sterile leakproof container Sterile collection tube Lesion biopsy, lesion swab, or plucked hairs that include roots and attached epithelial cells	<i>Trichophyton mentagrophytes</i> (<i>Dermatophytosis</i>)	3-5 days	\$105.33
<i>Salmonella enterica</i> PCR†	25001537	● Feces Minimum 2 g feces in a clean leakproof container.		3-5 days	\$101.65
SARS-COV-2 PCR†	25001538	● Feces, Culture swab Nasal swab (preferred), feces		5-7 days	\$187.20
<i>Strep equi</i> subsp. <i>zooepidemicus</i> PCR†	25001541	● RTT, Culture swab Nasal swab, lung tissue, lymph node		3-5 days	\$100.75
<i>Streptococcus equi</i> subsp. <i>equi</i> PCR†	25001540	● RTT, Culture swab Nasal swab, lung tissue, lymph node		3-5 days	\$101.65
Tickborne Disease PCR Panel Non-equine†	25001542	● ● RTT, LTT EDTA whole blood (minimum 0.5 mL), Fresh tissue	<i>Anaplasma</i> spp, <i>Ehrlichia</i> spp, <i>Borrelia burgdorferi</i> , <i>Rickettsia</i> spp, <i>Babesia</i> spp, <i>Bartonella</i> spp, <i>Hepatozoon</i> spp, <i>Theileria</i> spp	5-10 days	\$174.95

Tests (cont'd)

Name	Test Code	Sample Requirements	Components / Notes	TAT	Price
<i>Toxoplasma gondii</i> PCR†	25001543	●●● RTT, LTT, Feces Fresh tissue, urine, feces, EDTA whole blood (minimum 0.5 mL)		3-5 days	\$101.65
<i>Tritrichomonas foetus</i> PCR†	25001544	● Feces, Culture swab Preputial wash, vaginal scraping, cat feces		3-5 days	\$101.65
<i>Trueperella pyogenes</i> PCR†	25001545	● RTT, Culture swab Nasopharyngeal swab, transtracheal wash, bronchoalveolar lavage (BAL), respiratory tissue, abscess, other affected tissues		3-5 days	\$100.75
<i>Trypanosoma cruzi</i> PCR†	25001546	● LTT EDTA whole blood (minimum 0.5 mL)		3-5 days	\$100.75
<i>Ureaplasma</i> spp (genus) PCR†	25001547	● RTT, Culture swab Reproductive tissue or swabs. Respiratory washes, tissue or swabs		3-5 days	\$100.75
Uveitis PCR Panel, Feline†	25001548	● RTT, Culture swab Aqueous fluid, vitreous fluid, corneal scraping, conjunctival swab	<i>Toxoplasma gondii</i> , <i>Feline herpesvirus</i> , Feline coronavirus, <i>Cryptococcus</i> , <i>Bartonella</i> , FeLV, FIV	3-5 days	\$160.29
<i>Yersinia</i> spp (genus) PCR†	25001549	● Feces Minimum 2 g feces in a clean leakproof container		3-5 days	\$100.75



See the complete Directory of Services
online at zoetisreflabs.com/DOS

Policies

Policies

Species Policy

In order to maintain our high quality standards, we accept samples for core testing (hematology, chemistry, endocrinology, urinalysis) from only the following species: dogs, cats, horses and other domestic hoofstock (cattle, sheep, goats, pigs, llamas, alpacas), rabbits, and ferrets. Exotic species samples are accepted, but are generally analyzed at a third-party specialty laboratory.

For microbiology and parasitology testing, we accept samples from all species, except non-human primates.

For cytopathology, air-dried slides are accepted from all species. Liquid samples are accepted for all species except non-human primates.

For histopathology, formalin-fixed samples can be accepted from any species. Non-formalin-fixed whole bodies or heads are not accepted from any species.

Fresh, non-fixed samples from non-human primates will not be accepted for any test. For these species, tissue for histopathology must be formalin-fixed and slides for cytopathology must be air-dried or alcohol-fixed.

Samples from species that fall outside of this policy will be returned to you or will be sent to a third party laboratory with associated charges applied.

Add-on Tests

We make it easy to add tests to samples you've already submitted. Simply call customer service at **1-888-965-9652**.

Or, if you order your tests electronically, you can add a test through ZoetisDx or via PIMS.

Sample Retention

After testing is completed, we retain biological specimens as follows:

Item	Retention time
Non-fixed biologic specimens (eg, blood, urine, feces)	1 week
EXCEPTION: citrated plasma in -20°C	2 weeks
Blood smears	3 months
Cytology slides	6 months
Histopathology specimens - wet tissues	1 month
Histopathology specimens - FFPE (formalin-fixed paraffin-embedded) tissues	Minimum 3 years
Histopathology slides	Minimum 3 years

You agree that Zoetis will own and may use for any purpose all data collected in connection with this transaction once such data has been either (i) rendered anonymous in such a way that Your identity is not or is no longer identifiable (“anonymized”) or (ii) changed, transformed, re-arranged or otherwise manipulated into a usable and desired form (“transformed”). Zoetis and its affiliates may use all collected data which has not been anonymized or transformed for internal business purposes, as described in Zoetis’ Privacy Policy available at zoetis.com/privacy-policy.

If applicable, you agree that you have the right to submit all specimens, animals, and/or biological materials, and photographs and/or documentation thereof (collectively “Specimens”) submitted to the applicable Zoetis Reference Lab (“ZRL”) location, and you transfer such right(s) to the Specimens to ZRL. Any test results, diagnoses, or other analyses produced by ZRL resulting from these submissions will become the property of Zoetis; provided, however, that ZRL shall not disclose such information in a manner that identifies you, unless specifically required by state or federal agencies. Specimens may be used for research, teaching and commercial purposes. Specimens will not be returned unless specific arrangements are made, at your expense, in advance of ZRL’s receipt.

Policies (cont'd)

Compromised/Suboptimal Samples

Samples that cannot be processed in a manner that upholds our quality standards may require additional attention or, in some cases, may be rejected. In the event that we receive such a sample, we will attempt to contact you. Examples of suboptimal samples include: lack of clear patient identifying information, damaged or leaking containers, inadequate volume for testing, incorrectly preserved samples, and contaminated samples.

Employee Pet Testing Program

We want to reward you for being a loyal customer of Zoetis Reference Laboratories. That's why we created our Employee Pet Testing Program.

Eligible customers can receive up to 10% of their monthly spend in free employee testing. For example, if a customer spends \$3,000 per month in Zoetis Reference Laboratories tests, they are eligible to receive \$300 in employee pet testing.

Any amount of Employee Pet testing sent in that exceeds 10% of their spend will be billed fully.

All existing customers who are eligible and meet the criteria will be automatically able to participate in this program. Any new Zoetis Reference Laboratories customers will be automatically enrolled in the program for 90 days.

Customers can indicate that a test is for a staff pet by checking a box on the printed PIMS form labeled "Employee Pet."

Who is eligible?

- Any new or existing customer who spends a minimum of \$2,000 per month with Zoetis Reference Laboratories. If a customer's monthly spend drops below \$2,000, they are dropped from the eligibility list the following month.
- OptimizeDx customers are not eligible.

How to Make a Payment

To view invoices and statements or pay online, go to: www.zoetisus.com/billpay

For Zoetis Reference Laboratories billing questions, call **1-888-965-9652**.

Processing & Handling Fee Policy

At Zoetis, it is a priority to provide our valued customers with pricing transparency. A monthly \$50.00 fee will be applied to customer accounts spending \$500 or less per month with Zoetis Reference Laboratories.

If you have questions about this fee, please reach out to your Zoetis representative or contact customer service. You can also ask them for more information about programs that can help you avoid the monthly fee.

Zoetis Terms & Privacy

Zoetis values your trust and respects your privacy. The resources on this link will help you understand how we use your personal information and how you can manage your privacy choices.

To view the full Zoetis Privacy Policy, visit <https://www.zoetis.com/privacy-center>.

Alphabetical Test Listing

Canine & Feline Test Listing

#

2M Antibody Masticatory Muscle Myositis [†] 6403	84
17-Hydroxyprogesterone, Pre & Post-(ACTH Stim) [†] 17HPROG	59

A

Acetylcholine Receptor Antibody/Myasthenia Gravis [†] 6402	84
Acid-Fast Stain, Non-Avian/Reptile ACIDFAST	68
ACTH Stimulation (2 samples) 6100	58
ACTH Stimulation (3 samples) 6101	58
Adrenal Panel (with ACTH Stimulation), Canine/Feline [†] ACTHADR	59
Aerobic and Anaerobic Culture 7135	67
Aerobic Culture 7156	67
Albumin 1001	51
Aldosterone Baseline [†] ALDOS	59
Aldosterone Pre & Post [†] ALDOSPOST	59
Allergy Panel, Canine/Feline [†] 25001385	84
ALP (Alk Phos) 1003	52
ALT 1002	51
Amylase (Mich. State Univ.) [†] 1004	62
Anaerobic Culture 7157	67
Anaerobic ID Add-on A7157	89
<i>Anaplasma phagocytophilum</i> (IFA) ANATITER	73
Anticoagulant (Rodenticide) [†] ANTICOAG	77

Anti-Müllerian Hormone [†] 5555	60
Anti-Nuclear Antibody (IFA) ANA	84
Anti-Platelet Flow Cytometry, Canine [†] PBA	43
Antithrombin (AT3) [†] ANTTHRO3	40
<i>Aspergillus</i> Galactomannan Antigen EIA [†] ASPGUR	74
<i>Aspergillus</i> Titer AGID [†] ASP	73
AST 1005	52
Australia Export Completion Panel [†] AUSECPN	73

B

<i>Babesia</i> Blood Smear [†] BABBLDSMR	46
<i>Babesia</i> IFA [†] BABESIAIFA	76
Baermann Test (Lungworm) BAER	65
<i>Bartonella</i> IFA [†] 6025	74
Basic CBC, Canine/Feline (Healthy, Non-geriatric) 2000	39
Bile Acid Pre/Resting Add-on A6180	89
Bile Acids Post 6181	50
Bile Acids Pre/Resting 6180	50
<i>Blastomyces</i> Quantitative Antigen EIA [†] BLASTANT	70
Blood Culture 7150	68
Blood Lead [†] 6240	78
Blood Typing BLTYP	40
<i>Borrelia</i> (Lyme) (Western Blot) LYMEWB	73

Canine & Feline Test Listing (cont'd)

Bromide 1025	77
<i>Brucella canis</i> Culture [†] BRUCLCULT	68
<i>Brucella canis</i> Slide AGGL/AGID II [†] 6315	73
<i>Brucella canis</i> Titer IFA BRUCIFA	73
Buffy Coat 6295	44
BUN 1006	52
C	
Calcium (Total Calcium) 1008	52
Canine Adenovirus Type 1 Antibody (VN) [†] K9ADNE1	75
Canine C-reactive Protein (CRP) [†] CRP	51
Canine Diarrhea PCR [†] with O&P, <i>Giardia</i> 25001623	37
Canine Distemper Virus Titer (Cornell) [†] DISTEMP	74
Canine Influenza Virus Titer [†] K9FLUAB	76
Canine Parvovirus Fecal Test 60002	76
Carnitine [†] CARNITINE	88
Centrifugal Fecal Flotation Test (O&P) & McMaster (Fecal Egg Counting) Test 25000487	65
Centronuclear Myopathy [†] CNM	86
Chloride 1007	52
Cholesterol 1009	52
C-Kit Mutation Assay [†] 6541	81
ClinChem 18 CHEM18	49
ClinChem 18 with SDMA 25001624	49

ClinChem 25 1000	49
ClinChem 25 with SDMA 25001625	49
Clotting Factor FIX [†] FACTORIX	45
Clotting Factor FVIII [†] FACTORVIII	46
Clotting Factor FXII [†] FACTORXII	46
Clotting Factor X [†] FACX	46
CO2 1020	53
Coagulation Panel and Fibrinogen, Canine/Feline [†] 25001390	47
Coagulation Panel, Canine/Feline (Cornell) [†] COAGFA	45
Cobalamin and Folate KFCOFO	63
Cobalamin (Vitamin B12) 6207	63
<i>Coccidioides</i> Antibody Screening Test 25001690	70
<i>Coccidioides</i> Titer AGID [†] 6129	70
<i>Coccidioides</i> Titer AGID Add-on A6129	91
Complete Amino Acid Profile [†] COMPAA	84
Complete Muscle Profile [†] CM	82
Comprehensive CBC, Canine/Feline CBCC	39
Comprehensive CBC Upgrade (add to Bundle) CBCCUP	39
Comprehensive (PCR and Serology) Vector Panel, Canine [†] K9FCOMP	74
Comprehensive (PCR and Serology) Vector Panel, Feline [†] FELVECTPNL	76
Congenital (CKCS) Macrothrombocytopenia [†] CKCSMACRO	86
Copper, Zinc, Iron Levels on Hepatic Tissue [†] COPZNCIRON	82
Core Vaccine Titer Panel (includes Rabies, non-export), Canine [†] 6160	72

Canine & Feline Test Listing (cont'd)

Corticosteroid-Induced Alk Phos (CIALP) [†] CIALP	51
Cortisol Add-on A5101	90
Cortisol (Baseline/Resting) 5101	58
Cortisol (Single Sample, Post-ACTH) CORTPOST	58
Creatinine 1010	52
Creatinine Kinase (CK) 1011	52
Critical Care ClinChem, Canine/Feline 25000488	48
Crossmatch, Additional Donor XMTCHX	41
Crossmatch, One Donor XMTCH1	40
Crossmatch, Three Donors XMTCH3	41
Crossmatch, Two Donors XMTCH2	40
<i>Cryptococcus</i> Antigen Titer by Latex Agglutination 7216	70
Cyclosporine [†] 6335	77
Cystinuria Nitroprusside [†] CYSTINURIA	86
Cytopathology Complex 7201	83
Cytopathology Simple/Single Site 7120	83
D	
D-Dimer 4104	44
Degenerative Myelopathy [†] DEGMYE	86
Dermatopathology Panel (includes Histopathology & Dermatologist Consult) ZN400	80
Dermatophyte Culture 7160	67
Dexamethasone Suppression (2 samples) 6103	58

Dexamethasone Suppression (3 samples) 6104	59
DGGR Lipase (ZRL Premium Lipase) 1030	62
DGGR Lipase (ZRL Premium Lipase) Add-on A1030	89
DGGR Lipase (ZRL Premium Lipase), TLI, Cobalamin, and Folate, Canine 7250	62
DIC Panel [†] 6565	41
Digoxin [†] 6320	77
Direct Coombs 2701	44
Direct Coombs Add-on A2701	90
Drugs of Abuse [†] DOA	78

E

<i>Ehrlichia canis</i> (IFA) EHRCH	75
Endogenous ACTH, Canine [†] 25001384	60
Erythropoietin (EPO) [†] 50227	46
Estradiol, Baseline [†] ESTRADIOL	60

F

Fanconi Syndrome Test [†] FANCONI	86
Fecal Culture 7155	67
Fecal Direct Smear Test 25001610	65
Fecal Occult Blood 7255	64
Fecal Sedimentation Test 25001612	65
Feline Coronavirus (IFA) CORONAIFA	73
Feline Diarrhea PCR [†] with O&P, <i>Giardia</i> 25001626	38

Canine & Feline Test Listing (cont'd)

Feline Heartworm Antibody (ELISA) [†] 6007	66
Feline Heartworm Antibody (ELISA) Add-on [†] A6007	90
Feline Hyperthyroid ClinChem 25000507	48
Feline Infectious Peritonitis Virus (FIP) IFA 25001570	72
Feline Leukemia Virus (FeLV) (ELISA) 25000531	76
Feline Vaccine Panel 2 [†] 7200	71
FeLV/FIV Rapid Test 6005	70
FELV/FIV Rapid Test Add-on A6005	89
Fibrinogen, Canine/Feline 4103	44
FIV by Western Blot & ELISA [†] FIVWB	75
Flow Cytometry [†] FLOW	42
Flow Cytometry with CBC [†] FLOWCBC	42
Folate 6208	63
Free (Ionized) Calcium IONISTAT	50
Free (Ionized) Calcium (Mich. State Univ.) [†] 6304	50
Free (Ionized) Calcium (Mich. State Univ.) Add-on [†] A6304	92
Free T4 5103	56
Free T4 Add-on A5103	90
Free T4 by ED [†] 6330	56
Fructosamine 1024	57
Fructosamine Add-on A1024	89
Fungal Culture 7145	67
Fungal Susceptibility 1 Drug [†] FUNGALSUS1	69

Fungal Susceptibility 2 Drugs [†] FUNGALSUS2	69
Fungal Susceptibility 3 Drugs [†] FUNGALSUS3	69
Fungal Titer Panel AGID [†] 6130	71

G

Gabapentin [†] GAB	78
Gastrin [†] GASTRIN	61
GGT 1013	53
<i>Giardia</i> Antigen Test 2321	64
<i>Giardia</i> Antigen Test, Canine/Feline Add-on A2321	90
<i>Giardia</i> & <i>Cryptosporidium</i> (IFA), Canine/Feline 2404	65
<i>Giardia</i> Focused O&P 25001611	64
Glomerular Filtration Rate (Iohexol Clearance) [†] GFRIC	54
Glucose 1014	53
Gross Parasite Identification 25001435	66

H

Heartworm Antigen Test 6001	65
Heartworm Antigen Test Add-on A6001	89
Heartworm Antigen Test (Alternate Methodology) 25000533	66
Heartworm Antigen with Heat Treatment HWHEAT	65
Hemophilia Panel, Canine/Feline [†] HEMO	43
Hemostasis Panel, Canine/Feline [†] 6560	45
Hepatopathology (Liver) Panel ZN500	80

Canine & Feline Test Listing (cont'd)

High-Sensitivity Cardiac Troponin I [†] TROPONIN	51
Histopathology Complex, Level 1 ZN102	79
Histopathology Complex, Level 2 ZN103	79
Histopathology Simple, Multiple Sites (Up to 5 sites) ZN101	79
Histopathology Simple, Single Site ZN100	79
<i>Histoplasma</i> Quantitative Antigen EIA [†] HISTOANT	75
Hold Urine Culture for 48 hours 25000534	67

I

IgA, Canine IGA	84
IgG, Canine [†] IGG	84
IgG, IgA, IgM, Canine [†] IGGAM	84
Immunofixation Free Lt Chain Add-On [†] 25001431	85
Insulin:Glucose Ratio [†] MSUIG	57
Insulin-like Growth Factor [†] ILGROFAC	61
Iron Profile [†] IRONPRO	43
Iron, Serum/Plasma, Canine [†] IRON	43

K

Keppra Levels Peak & Trough Samples [†] KEPPRA	78
---	----

L

LC/Mass Spectrometry Toxicology Screen [†] LCMS	78
LDH [†] LDH	50

Leflunomide Level [†] LEFLUNO	77
<i>Leishmania</i> Titer IFA [†] LEISH	75
<i>Leptospira</i> Antibodies (IgM), Canine 6008	72
<i>Leptospira</i> Antibodies (IgM), Canine Add-on A6008	91
<i>Leptospira</i> MAT, 5 Serovar [†] 6164	72
<i>Leptospira</i> MAT, 5 Serovar & PCR [†] 6163	72
Levetiracetam (Keppra) [†] 6321	77
Liver ClinChem 25000510	48
Low Molecular Weight Heparin [†] LMWHEPARIN	47
Luteinizing Hormone (LH), Canine/Feline LH	61
Lyme (<i>Borrelia</i>) Titer IFA [†] LYMEIFA	75
Lyme Multiplex Test, Canine [†] LYMEMPLX	73
Lymphoid Neoplasia PCR/PARR (T&B Cell) from Blood, Fluid, Cytopathology Sample [†] 6542	83
Lymphoid Neoplasia PCR/PARR (T&B cell) from Histopathology Sample [†] 6544	82

M

Magnesium 1016	53
Malignancy (Hypercalcemia) Panel [†] 6555	57
Mast Cell Tumor Prognostic Panel, Canine [†] 6540	82
Mature Cat ClinChem 25000482	48
McMaster (Fecal Egg Counting) Test 2328	65
Melanoma Diagnostic Panel, Canine [†] 6425	82

Canine & Feline Test Listing (cont'd)

Melanoma Prognostic Panel, Canine [†] 6426	81
Metabolic ClinChem 25000512	48
Metanephrines, Plasma-free [†] METAPLAS	61
Metanephrines, Urine Fractionated [†] METAUR	61
Methylmalonic Aciduria (Cobalamin Malabsorption) [†] COBMAL	86
Microbiota Dysbiosis Index, Canine [†] DYSINDEX	64
Mini Renal ClinChem 25000513	48
Modified Knott's Test for Microfilariae (Heartworm) 6225	66
Modified Knott's Test with Quantification MODKNOTTSQ	66
<i>Mycobacterium</i> Culture [†] MYCO	68
<i>Mycoplasma</i> Culture [†] MYCOPLASCU	68
N	
<i>Neospora caninum</i> IFA [†] NEO	74
Nephropathology (Renal) Panel ZN25000541	81
Neurology Infectious Disease Panel [†] NEUROID	73
New Zealand Export Panel I [†] NZECPN	76
New Zealand Export Panel II [†] NZECPN2	76
NSAID ClinChem 25000514	48
NT-proBNP [†] PROBNP	51
O	
OFA Thyroid Panel, Canine [†] K9THYROID	57
O&P (Centrifugal Fecal Flotation) Panel 25001613	64

O&P (Centrifugal Fecal Flotation) with Direct Smear 25001614	65
O&P with HW & Tick 25001627	66
O&P with HW & Tick 25001627	70
Osmolality [†] OSMOLAL	51
Ova and Parasite (Centrifugal Fecal Flotation) Test 2320A	64
Ova and Parasite (Centrifugal Fecal Flotation) Test and <i>Giardia</i> Antigen Test 2320	64
Ovarian Remnant Syndrome [†] OVRS	60
P	
Pancreatic Lipase Immunoreactivity (PLI), Canine (TX GI) [†] PLITX	62
Pancreatic Lipase Immunoreactivity (PLI), Feline (TX GI) [†] 7240	62
Parathyroid Hormone-Related Protein (PTHrP) [†] PTHRP	57
Parathyroid Panel [†] 6306	57
Pathologist Review PATH	39
Phenobarbital 1031	78
Phenobarbital Add-on A1031	91
Phenobarbital ClinChem 25000515	49
Phosphorus 1017	53
PLI and TLI, Feline (TX GI) [†] 7253	62
PLI, Cobalamin, and Folate, Canine/Feline (TX GI) [†] FEPLICOFOL	62
PLI, TLI, Cobalamin, and Folate, Canine (TX GI) [†] K9PTCFTX	62
PLI, TLI, Cobalamin, and Folate, Feline (TX GI) [†] 7251	63
Potassium 1018	53

Canine & Feline Test Listing (cont'd)

Pre-Anesthetic ClinChem 25000516	49
Pregnancy Test (Relaxin), Canine/Feline 6009	61
Premium Thyroid Panel + TSH, Feline† 6550	56
Premium Thyroid Panel, Canine† 6545	56
Progesterone 5105	60
Protein C† PROC	43
Protein Electrophoresis 6345	84
Protein Electrophoresis Immunofix† 25001430	85
PT & aPTT 4101	43
Pyruvate Kinase Deficiency DNA Test† PYRUDNA	86
<i>Pythium</i> Serology† OOPYTHIUM	75
Q	
Quantitative Copper (Blood or Tissue)† 6171	77
R	
Rabies Antibody Endpoint (Not for Export)† 7206	75
Rabies Screen Titer (not for export)† 7205	74
Rabies Titer for Export† 7205E	71
Renal ClinChem 25000479	48
Reproductive Hormone Panel† RHP	60
Reticulocyte RETIC	39
Reticulocyte Count Add-on 25000525	92
<i>Rickettsia rickettsii</i> (RMSF) IFA† RMSFIFA	75

S

SDMA 7037	50
SDMA Add-on AN7037	89
Selenium† SELENIUM	87
Skin Culture and Sensitivity (Dermatology) 25001434	69
Sodium 1019	53
Stone Analysis† 7032	55

T

Taurine:Plasma† TAURPLAS	87
Taurine (Plasma and Whole Blood)† TAUR	87
Taurine:Whole Blood† TAURWB	87
Technician Blood Smear Review SMREV	45
Testosterone, Baseline† TESTOS	60
Testosterone Post-HCG/GnRH† TESTPOST	60
Thiamine (Vitamin B1)† THIAMINE	87
Thyroglobulin Autoantibody† 6437	56
Thyroid Panel 25000484	56
Tick Panel Screening Test + Heartworm Antigen Test 25001389	76
Tick Panel Screening Test 6002	70
Tick Panel Screening Test Add-on A6002	90
Tick Serology Panel 7185	71
Tick with HW 25001628	71

Canine & Feline Test Listing (cont'd)

TLI, Canine 6206	63
TLI Canine Add-on A6206	90
TLI, Cobalamin, and Folate, Canine 7252	63
TLI, Cobalamin, and Folate, Feline (TX GI) [†] F7251	63
Total Bilirubin 1021	53
Total Protein 1022	54
Total T3 (TT3), Canine/Feline [†] 1396	56
Total T4 5102	56
Total T4 Add-on A5102	91
Toxic Heavy Metal Panel [†] 25001395	78
<i>Toxoplasma gondii</i> ELISA [†] TOXAB189	71
<i>Toxoplasma gondii</i> ELISA (Canine) & <i>Neospora caninum</i> IFA [†] K9TOXOPL	74
Trace Element Screen [†] TMINPAN	87
Triglycerides 1023	50
Triglycerides Add-on A1023	90
<i>Trypanosoma cruzi</i> (Chagas Disease) IFA Titer [†] TCRUIFA	74
Trypsin-Like Immunoreactivity (TLI), Feline (TX GI) [†] 6203	63
TSH 5106	56
TSH Add-on A5106	91
U	
Urinalysis 3000	54
Urinalysis Add-on 25001700	92

Urinalysis with Reflex UPC 25000521	55
Urinalysis with Reflex UPC Add-on 25001701	92
Urinalysis with Reflex UPC Upgrade 25001702	55
Urinalysis with Reflex UPC Upgrade 25001702	104
Urinalysis with Urine Culture & Sensitivity 25001629	37
Urinalysis with Urine Culture & Sensitivity 25001629	54
Urinalysis with Urine Culture & Sensitivity 25001629	67
Urine Cortisol:Creatinine Ratio 3310	59
Urine Cortisol:Creatinine Ratio Add-on [†] A3310	91
Urine Culture and Sensitivity 7140	67
Urine Culture and Sensitivity Add-on A7140	91
Urine Electrolyte Panel (for Fractional Excretion) [†] FREXCRE	55
Urine Protein:Creatinine Ratio 1555	54
Urine Protein:Creatinine Ratio Add-on A1555	91
V	
Vaccine Titer (Qualitative), Canine 6145	71
Vector Serology Panel, Canine [†] EXTENTICK	71
Visual Platelet Estimate 2211	45
Vitamin A [†] VITA	88
Vitamin D (25-hydroxyvitamin D3) [†] VITD	58
Vitamin D Profile [†] VITP	58
Vitamin E [†] 2331	88

Canine & Feline Test Listing (cont'd)

Vitamin E and Mineral Panel† VITEMINE	88
von Willebrand factor (vWF:Ag)† 4106	44
Z	
Zinc, Serum/Plasma† ZINC	77
Zonisamide† 6322	78
ZRL™ Adult 25000501	33
ZRL™ Adult Complete AWP	33
ZRL™ Adult Complete with O&P, HW, Tick 25001631	33
ZRL™ Adult with HW 25000502	33
ZRL™ Adult with O&P 25001630	33
ZRL™ Canine GI 25000491	36
ZRL™ Canine GI with O&P, <i>Giardia</i> 25001632	36
ZRL™ Coagulation Panel I COAG	47
ZRL™ Coagulation Panel II 25001391	47
ZRL™ Complete Coccidioides 25001711	36
ZRL™ Critical Care 25000508	36
ZRL™ Feline GI 25000492	36
ZRL™ Feline GI with O&P, <i>Giardia</i> 25001633	37
ZRL™ Feline Hyperthyroid 25000490	36
ZRL™ Junior 25001634	32
ZRL™ Junior Complete JWP	32
ZRL™ Junior Complete with HW 25001642	32

ZRL™ Junior with HW 25001638	32
ZRL™ Junior with HW, Tick 25001639	32
ZRL™ Junior with O&P 25001635	32
ZRL™ Junior with O&P, <i>Giardia</i> , HW, Tick 25001637	32
ZRL™ Junior with O&P, HW, Tick 25001636	32
ZRL™ Junior with PT/aPTT 25001640	33
ZRL™ Junior with T4 25001641	32
ZRL™ Mature Cat 2500050v0	36
ZRL™ Metabolic 25000493	37
ZRL™ Pre-Anesthetic 25000496	37
ZRL™ Premium Lipase (DGGR Lipase) 1030	62
ZRL™ Premium Lipase (DGGR Lipase) Add-on A1030	89
ZRL™ Renal 25000478	37
ZRL™ Senior SBASIC	34
ZRL™ Senior Complete SWP	33
ZRL™ Senior Complete O&P, HW, Tick, UPC 25001654	29
ZRL™ Senior Complete O&P, HW, Tick, UPC 25001654	35
ZRL™ Senior Complete O&P with FeLV/FIV 25001646	34
ZRL™ Senior Complete O&P with FeLV/FIV, HW 25001647	35
ZRL™ Senior Complete O&P with fPLI 25001649	35
ZRL™ Senior Complete O&P with <i>Giardia</i> 25001650	35
ZRL™ Senior Complete O&P with <i>Giardia</i> , HW, Tick 25001651	35
ZRL™ Senior Complete O&P with HW 25001648	35

Canine & Feline Test Listing (cont'd)

ZRL™ Senior Complete O&P with HW, Tick 25001652	35
ZRL™ Senior Complete O&P with HW, Tick, Zoetis Premium Lipase (DGGR), TG 25001653	36
ZRL™ Senior Complete with Zoetis Premium Lipase (DGGR), TG 25001655	34
ZRL™ Senior Complete with FT4 25001656	34
ZRL™ Senior Complete with HW, TG 25001657	34
ZRL™ Senior with FT4 25001645	34
ZRL™ Senior with O&P 25001643	34
ZRL™ Senior with O&P, HW 25001644	34

Large Animal Test Listing

A

Acid-Fast Stain, Non-Avian/Reptile ACIDFAST	112
Aerobic and Anaerobic Culture 7135	111
Aerobic Culture 7156	111
Albumin 1001	100
Aldosterone Baseline [†] ALDOS	108
Allergy Panel, Equine [†] 6307	122
ALP (Alk Phos) 1003	101
Amylase (Mich. State Univ.) [†] 1004	109
Anaerobic Culture 7157	111
Anaerobic ID Add-on A7157	125
Anticoagulant (Rodenticide) [†] ANTICOAG	117
Anti-Müllerian Hormone, Equine [†] ANTIMULL	107
Antithrombin (AT3) [†] ANTTHRO3	98
<i>Aspergillus</i> Galactomannan Antigen EIA [†] ASPGUR	115
<i>Aspergillus</i> Titer AGID [†] ASP	114
AST 1005	101

B

Baermann Test (Lungworm) BAER	110
Basic CBC, Large Animal (Healthy Animals) 7000	97
<i>Blastomyces</i> Quantitative Antigen EIA [†] BLASTANT	114
Blood Culture 7150	112

Blood Lead [†] 6240	118
Buffy Coat 6295	97
BUN 1006	101

C

Calcium (Total Calcium) 1008	101
Caseous Lymphadenitis (<i>Corynebacterium pseudotuberculosis</i>) Titer [†] CALYMTTR	116
Centrifugal Fecal Flotation Test (O&P) & McMaster (Fecal Egg Counting) Test 25000487	110
Chloride 1007	101
Cholesterol 1009	101
ClinChem 18, Large Animal CHEM18EQ	99
CO ₂ 1020	103
<i>Coccidioides</i> Titer AGID [†] 6129	114
Coggins AGID, Electronic Submission E1560	113
Coggins AGID, Paper Submission COGAGID	113
Coggins ELISA, Electronic Submission E1561	113
Coggins ELISA, Paper Submission COGELISA	113
Complete Amino Acid Profile [†] COMPAA	122
Comprehensive CBC, Large Animal 8000	97
Comprehensive CBC Upgrade, Large Animal 25001590	97
Cortisol Add-on A5101	125
Cortisol AM/PM, Equine AMPMEC	105

Large Animal Test Listing (cont'd)

Cortisol (Baseline/Resting) 5101	105
Creatinine 1010	101
Creatinine Kinase (CK) 1011	102
<i>Cryptococcus</i> Antigen Titer by Latex Agglutination 7216	115
Cryptorchid Panel, Equine [†] EQCRYPTO	107
Cytopathology Complex 7201	121
Cytopathology Simple/Single Site 7120	121

D

Dermatopathology Panel (includes Histopathology & Dermatologist Consult) ZN400	120
Dermatophyte Culture 7160	112
Drug Screen, Firocoxib [†] DGSCRFIRO	117
Drug Screen, NSAIDs, Equine [†] DSNSAIDS	118
Drug Screen with Tranquilizers, NSAIDs [†] DSNSAIDTQ	117
Drugs of Abuse [†] DOA	118

E

Endogenous ACTH ENDOACTH	106
Equine Influenza (HAI) [†] EQINFLUENZA	116
Equine Metabolic Panel (ACTH, Insulin, Glucose) EQCUSHING	106
Equine Protozoal Myelitis IFA Titer (NeoFluor™) [†] NEOIFAT	116
Equine Protozoal Myelitis IFA Titer (SarcoFluor™) [†] SARCOIFAT	116
Equine Protozoal Myelitis IFA Titers (SarcoFluor™ & NeoFluor™) [†] EPMIFAT	115

Equine Uterine Culture 7154	111
Equine Viral Arteritis [†] EVA	114
Erythropoietin (EPO) [†] 50227	98
Estrogens, Total [†] EQEST	108
Estrone sulfate, Equine 25001386	106
Exercise ClinChem, Large Animal 25000509	99

F

Fecal Culture 7155	112
Fecal Occult Blood 7255	110
Free (Ionized) Calcium IONISTAT	100
Free (Ionized) Calcium (Mich. State Univ.) [†] 6304	100
Free (Ionized) Calcium (Mich. State Univ.) Add-on [†] A6304	125
Free T4 5103	105
Free T4 Add-on A5103	125
Free T4 by ED [†] 6330	105
Fungal Culture 7145	112
Fungal Susceptibility 1 Drug [†] FUNGALSUS1	112
Fungal Susceptibility 2 Drugs [†] FUNGALSUS2	112
Fungal Susceptibility 3 Drugs [†] FUNGALSUS3	113

G

Gabapentin [†] GAB	118
Gastrin [†] GASTRIN	107

Large Animal Test Listing (cont'd)

GCT Panel II, Equine [†] EQGCTPII	107
GGT 1013	102
GLDH 1040	102
Glucose 1014	102
Gross Parasite Identification 25001435	110

H

High-Sensitivity Cardiac Troponin I [†] TROPONIN	101
Histopathology Complex, Level 1 ZN102	119
Histopathology Complex, Level 2 ZN103	119
Histopathology Simple, Multiple Sites (Up to 5 Sites) ZN101	119
Histopathology Simple, Single Site ZN100	119
<i>Histoplasma</i> Quantitative Antigen EIA [†] HISTOANT	115
Hold Urine Culture for 48 hours 25000534	111

I

IgG, Equine FOALIGG	122
IgG, Equine: Rapid, Semi-Quantitative FOALIGGRAP	122
IgM, Equine 25000540	122
Insulin, Equine 5104	108
Insulin:Glucose Ratio, Equine EQINSGLU	108
Insulin-like Growth Factor [†] ILGROFAC	106
Iron Profile [†] IRONPRO	98

L

Large Animal ClinChem 24 CHEM24EQ	99
Large Animal Exercise Bundle 25000489	95
Large Animal Liver ClinChem 25000511	99
Large Animal Metabolic Bundle 25000494	95
Large Animal Pre-Anesthetic Bundle 25000495	95
Large Animal Pre-Anesthetic ClinChem 25000486	99
Large Animal Renal Bundle I 25000497	95
Large Animal Renal Bundle II 25000498	95
Large Animal Renal ClinChem 25000485	99
Large Animal Wellness Bundle EQWP	95
LC/Mass Spectrometry Toxicology Screen [†] LCMS	118
LDH [†] LDH	100
<i>Leptospira</i> MAT, 5 Serovar [†] 6164	114
<i>Leptospira</i> MAT, 5 Serovar & PCR [†] 6163	114
Lipid Profile [†] LIPIDPROF	103
Low Molecular Weight Heparin [†] LMWHEPARIN	98
Lyme Multiplex Test, Equine [†] LYMEMPLXEQ	115

M

Magnesium 1016	102
McMaster (Fecal Egg Counting) Test 2328	110
Metabolic Syndrome Panel I, Equine [†] CUSHING	106

Large Animal Test Listing (cont'd)

Metanephrines, Plasma-free† METAPLAS	108
Metanephrines, Urine Fractionated† METAUR	109
<i>Mycobacterium</i> Culture† MYCO	112
<i>Mycoplasma</i> Culture† MYCOPLASCU	111
N	
Nephropathology (Renal) Panel ZN25000541	120
O	
Oral Sugar Test (Insulin Testing), Equine† EQGLUCOSE	108
Osmolality† OSMOLAL	100
Ova and Parasite (Centrifugal Fecal Flotation) Test 2320A	110
P	
Parathyroid Hormone-Related Protein (PTHrP)† PTHRP	109
Pathologist Review PATH	97
Phosphorus 1017	102
Potassium 1018	102
Potomac Horse Fever (<i>Neorickettsia risticii</i>) IFA† POTOMAC	116
Pregnancy Panel, Equine† 25001387	107
Pregnant Mare Serum Gonadotropin (PMSG)† 25001388	107
Progesterone 5105	105
Protein C† PROC	98
Protein Electrophoresis 6345	122

PT & aPTT 4101	98
<i>Pythium</i> Serology† OOPYTHIUM	115
Q	
Quantitative Copper (Blood or Tissue)† 6171	117
R	
Reflex Aerobic Culture (Contingent) 25000536	111
Reproductive Hormone Panel† RHP	107
Reticulocyte RETIC	97
Reticulocyte Count Add-on 25000525	125
<i>Rhodococcus equi</i> ELISA† RHODOCOCCUS	116
S	
<i>Sarcocystis neurona</i> by Western Blot† SNEURONA	115
Selenium† SELENIUM	123
Selenium and Trace Elements Screen† SELETREL	117
Serum Amyloid A (SAA), Equine† SAA	100
Skin Culture and Sensitivity (Dermatology) 25001434	113
Sodium 1019	102
<i>Streptococcus equi</i> ELISA† STREPEQUI	116
T	
Technician Blood Smear Review SMREV	97
Testosterone, Baseline† TESTOS	108

Large Animal Test Listing (cont'd)

Testosterone Post-HCG/GnRH [†] TESTPOST	108
Thiamine (Vitamin B1) [†] THIAMINE	123
Thyroid Panel EQTHYROID	105
Total Bilirubin 1021	103
Total Protein 1022	103
Total T3 (TT3), Equine TT3	105
Total T4 5102	105
Total T4 Add-on A5102	125
Toxic Heavy Metal Panel [†] 25001395	118
Trace Element Screen [†] TMINPAN	123
TRH Stimulation Test for Equine Cushing's/PPID 2222	106
Triglycerides 1023	103
Triglycerides Add-on A1023	125

U

Urinalysis 3000	104
Urinalysis Add-on 25001700	126
Urinalysis with Reflex UPC 25000521	104
Urinalysis with Reflex UPC upgrade 25001702	104
Urine Culture and Sensitivity 7140	111
Urine Electrolyte Panel (for Fractional Excretion) [†] FREXCRE	104
Urine Protein:Creatinine Ratio 1555	104
Urine Protein:Creatinine Ratio Add-on A1555	126

V

Vitamin A [†] VITA	124
Vitamin D (25-hydroxyvitamin D3) [†] VITD	107
Vitamin D Profile [†] VITP	109
Vitamin E [†] 2331	123
Vitamin E and Mineral Panel [†] VITEMINE	123

W

West Nile Virus IgM ELISA [†] WESTNILE	116
---	-----

Z

Zinc, Serum/Plasma [†] ZINC	117
--	-----

Avian, Exotic & Small Mammal Test Listing

A

Advanced Reptile Panel† 25001448	128
Adv Avian Chem Panel† 25001437	130
Adv Mammal Chem Panel† 25001438	131
Adv Reptile Chem Panel† 25001439	131
Aspergillus Panel – PenguinPLUS † 25001449	132
Aspergillus Panel PLUS † 25001440	132

B

Basic Avian Chemistry Panel† 25001450	131
Basic Mammal Chemistry Panel† 25001451	131
Basic Mammal Panel† 25001452	128
Basic Reptile Chemistry Panel† 25001453	131
Basic Reptile Panel† 25001454	128
Bile Acids † 25001442	131

C

CBC Automated Count Manual Diff Mammal† 25001443	130
CBC Hemocytometer Cnt Man Diff Reptile† 25001444	129

E

Elasmobranch CBC and Diff† 25001455	129
Elasmobranch Chemistry Panel† 25001456	131

F

Fibrinogen Heat Precipitation Method† 25001457	129
--	-----

L

Lipoprotein Chemistry Panel† 25001458	130
---	-----

M

Mini Avian Panel † 25001459	129
Mini Elasmobranch Panel† 25001460	129
Mini Mammal Panel† 25001461	130
Mini Reptile Panel† 25001462	129
Mini Teleost Panel† 25001463	130

P

Protein Electrophoresis (EPH) † 25001446	133
--	-----

R

Reptile CBC and Diff† 25001447	129
--	-----

S

Sex Determination (Psittacine)† 25001539	132
--	-----

T

T4, Total (Mammals Only)† 25001464	132
Teleost CBC and Diff† 25001465	130
Teleost Chemistry Panel† 25001466	131

Avian, Exotic & Small Mammal Test Listing (cont'd)

V

Vitamin D (25-Hydroxy) Mammals Only† 25001467	133
---	-----

W

Well Bird Panel Package† 25001468	128
-------------------------------------	-----

PCR Test Listing

A

Adult Canine Diarrhea PCR Panel† 25001618	137
Adult Equine Diarrhea PCR Panel† 25001469	135
<i>Anaplasma/Ehrlichia</i> PCR† 25001470	135
Anemia PCR Panel, Canine† 25001471	135
Anemia PCR Panel, Feline† 25001472	135

B

<i>Babesia</i> spp (genus) PCR† 25001473	135
<i>Bartonella</i> spp (genus) PCR† 25001474	135
Blood Bank Disease PCR Panel, Canine† 25001475	135
<i>Bordetella bronchiseptica</i> PCR† 25001476	135
<i>Borrelia burgdorferi</i> (Lyme) PCR† 25001477	136
<i>Brucella</i> spp (genus) PCR† 25001478	136

C

<i>Campylobacteri coli</i> PCR† 25001480	136
<i>Campylobacteri jejuni</i> PCR† 25001479	136
Canine Adenovirus PCR† 25001481	136
Canine Puppy Diarrhea PCR w O&P, Giardia 25001623	136
Canine Distemper Virus (CDV) PCR† 25001482	137
Canine Herpesvirus 1 (CHV-1) PCR† 25001483	137
Canine Influenza PCR† 25001484	137
Canine Mycoplasma PCR† 25001485	137

Canine Neurological Disease, Central Nervous System PCR Panel† 25001486	137
Canine Parvovirus (CPV) PCR† 25001487	138
Canine Puppy Diarrhea PCR Panel† 25001488	138
Canine Reproductive/Abortion PCR Panel† 25001489	138
Canine Respiratory PCR Panel† 25001490	138
<i>Chlamydia</i> spp (genus) PCR† 25001491	138
<i>Chlamydophila felis</i> PCR† 25001492	138
Comprehensive (PCR and Serology) Vector Panel, Canine† K9FCOMP	139
Comprehensive (PCR and Serology) Vector Panel, Feline† FELVECTPNL	139
Coronavirus - Pan/Consensus PCR† 25001493	139
<i>Cryptococcus neoformans</i> PCR† 25001494	139
<i>Cryptosporidium</i> spp (genus) PCR† 25001495	139
<i>Cytauxzoon felis</i> PCR † 25001496	139

E

<i>Ehrlichia</i> spp (genus) PCR† 25001498	139
<i>Encephalitozoon cuniculi</i> PCR† 25001499	139
Equine Abortion/Reproduction PCR Panel† 25001500	140
Equine Foal Diarrhea PCR Panel† 25001501	140
Equine Herpes Virus type 1 PCR† 25001502	140
Equine Herpes Virus type 4 PCR† 25001503	140
Equine Influenza PCR† 25001504	140
Equine Respiratory PCR Panel I† 25001505	140

PCR Test Listing (cont'd)

Equine Respiratory PCR Panel II [†] 25001506	141
Equine Tickborne Disease PCR Panel [†] 25001507	141
Eubacteria Universal DNA - 16S rRNA PCR [†] 25001508	141
F	
Feline Calicivirus PCR [†] 25001509	141
Feline Coronavirus PCR [†] 25001510	141
Feline Diarrhea PCR Panel [†] 25001511	141
Feline Diarrhea PCR [†] with O&P, <i>Giardia</i> 25001626	142
Feline Hemotropic (blood) <i>Mycoplasma</i> PCR [†] 25001512	142
Feline Herpes and <i>Chlamydomphila felis</i> PCR Combo [†] 25001513	142
Feline Herpesvirus PCR [†] 25001514	142
Feline Immunodeficiency Virus FIV PCR [†] 25001515	142
Feline Leukemia Virus (FeLV) PCR [†] 25001516	142
Feline Neurological Disease PCR Panel [†] 25001517	142
Feline Ocular PCR Panel [†] 25001518	143
Feline Panleukopenia (Parvovirus) PCR [†] 25001519	143
Feline Respiratory PCR Panel 1 [†] 25001520	143
Fever of Unknown Origin PCR Panel, Canine [†] 25001521	143
Fever of Unknown Origin PCR Panel, Feline [†] 25001522	143
<i>Francisella tularensis</i> (Tularemia) PCR [†] 25001523	143
Fungal Sequencing and PCR on Tissue Block FUNGSEQPCR	143

G

<i>Giardia intestinalis</i> PCR [†] 25001524	144
---	-----

H

<i>Helicobacter</i> spp (genus) PCR [†] 25001525	144
---	-----

L

<i>Leishmania</i> spp (genus) PCR [†] 25001526	144
---	-----

<i>Leptospira</i> MAT, 5 Serovar & PCR [†] 6163	144
--	-----

<i>Leptospira</i> spp (genus) PCR [†] 25001527	144
---	-----

<i>Listeria</i> spp (genus) PCR [†] 25001528	144
---	-----

M

<i>Mycobacterium</i> spp (genus) PCR [†] 25001529	145
--	-----

<i>Mycoplasma</i> spp (genus) PCR [†] 25001530	145
---	-----

N

<i>Neospora caninum</i> PCR [†] 25001531	145
---	-----

<i>Nocardia</i> spp (genus) PCR [†] 25001532	145
---	-----

P

Pan-Fungal PCR [†] 25001533	145
--	-----

PCR Formalin-Fixed Paraffin Embed Tissue [†] 25001534	145
--	-----

R

<i>Rhodococcus equi</i> PCR [†] 25001535	145
---	-----

PCR Test Listing (cont'd)

<i>Rickettsia</i> spp (genus) PCR [†] 25001536	146
Ringworm, PCR, Canine/Feline 25001616	146
Ringworm, PCR, Equine 25001617	146

S

<i>Salmonella enterica</i> PCR [†] 25001537	146
SARS-COV-2 PCR [†] 25001538	146
<i>Strep equi</i> subsp. <i>zooepidemicus</i> PCR [†] 25001541	146
<i>Streptococcus equi</i> subsp. <i>equi</i> PCR [†] 25001540	146

T

Tickborne Disease PCR Panel Non-equine [†] 25001542	146
<i>Toxoplasma gondii</i> PCR [†] 25001543	147
<i>Tritrichomonas foetus</i> PCR [†] 25001544	147
<i>Trueperella pyogenes</i> PCR [†] 25001545	147
<i>Trypanosoma cruzi</i> PCR [†] 25001546	147

U

<i>Ureaplasma</i> spp (genus) PCR [†] 25001547	147
Uveitis PCR Panel, Feline [†] 25001548	147

Y

<i>Yersinia</i> spp (genus) PCR [†] 25001549	147
---	-----



LOOK DEEPER